AMEDA FINESSE™

BREAST PUMP | TIRE-LAIT | SACALECHE

EN | Instructions for use
FR | Mode d'emploi
ES | Instrucciones de uso

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Congratulations on joining the Ameda family!

We are excited to help you meet your breastfeeding goals along your breastfeeding journey with your baby. With our Finesse™ breast pump and HygieniKit® Milk Collection System, you can express and collect your breast milk safely and effectively.

The Ameda HygieniKit is a closed system that contains a diaphragm that keeps the air from the pump and tubing separate from the pump parts that are exposed to your breast milk, so your tubing stays dry. This diaphragm helps prevent milk backup into tubing.

If you have any questions regarding your Finesse™ breast pump or HygieniKit, we have a certified customer service ParentCare team available to help you. Have your breast pump and parts nearby when you call. ParentCare is available by phone during business hours on weekdays at **1-866-992-6332(1-866-99-AMEDA)**. Please check www.ameda.com for up-to-date available times.

For more information about breastfeeding, breast pumping or Ameda products visit www.ameda.com. For more product information or feedback in other countries, visit www.ameda.com for a list of distributors outside of the United States.

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TABLE OF CONTENTS

1. WARNINGS AND CAUTIONS.....	7-9
2. INTENDED USE.....	9
3. CONTRAINDICATIONS FOR USE.....	10
4. PRODUCT DESCRIPTION.....	10
5. SET-UP INSTRUCTIONS.....	11-13
5.1 Before First Use: Cleaning The HygieniKit.....	11
5.2 Assembling The HygieniKit.....	12-13
6. INSPECTING BEFORE PUMPING.....	14
7. OPERATING INSTRUCTIONS.....	14-17
7.1 Plug into Power Source.....	14-15
7.2 Important to Know Before Breast Pumping.....	15
7.3 To Begin Pumping.....	16
7.4 How to Fit Your Flange.....	17
7.5 Storing Pumped Milk in Bottles.....	17
8. CLEANING THE HYGIENIKIT.....	18-19
9. MAINTENANCE.....	19-20
10. TROUBLESHOOTING.....	21-22
11. WARRANTY.....	23
12. TECHNICAL DATA.....	24-27
13. STORAGE AND OPERATIONAL CONDITIONS.....	27
14. TRAVEL OR INTERNATIONAL USE.....	27
15. AMEDA ACCESSORIES AND SPARE PARTS.....	28

Icon Glossary

	Warning		Temperature Limitation For Transport/Storage		Serial Number
	Caution		Humidity Limitation For Transport/Storage		Catalogue Number
	Important		Temperature Limitation For Operation		Unit of Measure
	Washing Instruction		Humidity Limitation For Operation		Class II/Double Insulated
	Yes		Pressure Limitation For Operation/Transport/Storage		Type BF Applied Part
	No		Expiration		Solid Particle/Liquid Ingress Classification
	Consult Instructions for Use		Manufacturer		UL Certified
	Fragile		This product meets the applicable European regulations and is under the supervision of a Notified Body.		Waste Electrical and Electronic Equipment
	Protect from sunlight		Authorized Representative in the European Community		
	Keep Dry				

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1. WARNINGS AND CAUTIONS

IMPORTANT SAFEGUARDS

When using electrical products, especially if children are present, the following basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT

DANGER: To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- **DO NOT** use while bathing, showering or swimming.
- **DO NOT** place or store where product can fall or be pulled into a bathtub, shower, sink, or pool.
- **DO NOT** place or drop product into water or other liquid.
- **DO NOT** reach for electrical product that has fallen into water. Immediately unplug from electrical outlet.

The **WARNING** symbol identifies important safety instructions that if not followed can lead to injury or damage to the breast pump. The warning symbols when used with the following words stand for:

	WARNING Can lead to serious injury or death.
	CAUTION Can lead to minor injury.
	IMPORTANT Can lead to material damage.

WARNING: to avoid fire, electrocution, heart arrhythmia, serious burns, or injury:

- **NEVER** leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when product is used near children or persons with disabilities.
- Use product for its intended use **ONLY**, as described in this manual.
- **DO NOT** use parts that are not recommended by the manufacturer.

- Prior to use, always inspect power adapter for damage or exposed wires. If damage and/or exposed wires are found, do not use power adapter and call Ameda Customer Service ParentCare at **1.866.992.6332 (1.866.99. AMEDA)**.
- **NEVER** operate electrical product: if it has a damaged cord or plug; if it is not working properly; if it has been dropped, damaged, or becomes wet.
- Keep cord away from heated surface.
- **DO NOT** install batteries backwards, charge, put in fire, or mix with other battery types. May explode or leak causing injury. Replace all batteries at the same time.
- **NEVER** use while sleeping or when drowsy.
- **NEVER** insert or drop any object into any opening or into tubing.

Environmental Conditions That Affect Use

- **DO NOT** use outdoors with electrical cord.
- **DO NOT** operate where oxygen or aerosol spray products are being used.

WARNING: To avoid health risk and reduce risk of injury including engorgement, mastitis, and minor lacerations:

This product is a personal care item and intended for use by only one user. This product should not be shared between users or resold or given from one user to another. Pumping does not prevent transmission of substances or infections carried via breast milk to baby.

Note: Do not use product if damaged. Use product only for intended use as described in this manual.

- **DO NOT** use pump while driving.
- Pumping can induce labor, **DO NOT** use while pregnant unless approved by your breastfeeding specialist or healthcare professional.
- Prior to first use and after every use, clean all parts (flange, diaphragm, valve and bottle) that come in contact with your breast and breastmilk.
- **DO NOT** use attachments other than those recommended by manufacturer
- **DO NOT** make mechanical or electrical modifications to breast pump base.
- Inspect all appropriate pump components before each use. Check the parts for cracks, chips, tears or wear. **DO NOT** use, contact Ameda Customer Service ParentCare at **1.866.992.6332 (1.866.99. AMEDA)** immediately.
- **DO NOT** continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- **DO NOT** continue to pump if the device is causing pain.
- **DO NOT** heat or thaw frozen breastmilk in a microwave or boiling water. Microwaving can heat liquids unevenly. Both microwaving and boiling water can result in possible burns.
- Undesirable side effects: Skin irritation, allergic reaction, discomfort, pain, and soreness. Follow all instructions to minimize undesirable side effects.

STRANGULATION HAZARD: Keep tubing and cords out of reach of children.

CHOKING HAZARD: Small parts, **NOT** for children under 3 years.

CAUTION: Can lead to minor injury:

- If batteries are in pump base, remove batteries before using the power adapter.
- **DO NOT** wrap cord around the body of the adapter.
- Use **ONLY** the AC power adapter that is recommended for the Finesse™ Breast Pump.
- **DO NOT** shorten tubing or modify the breast pump kit.
- Make sure the voltage of the power adapter is compatible with the power source.
- Plug the AC power adapter into the pump base prior to plugging into the wall socket.
- Remove AC power adapter from power source and remove batteries before wiping outside of pump base clean with a damp cloth.
- Keep away from dust, lint and direct sunlight, as damage may occur to breast pump where speed and suction may be reduced.
- **DO NOT** attempt to remove the breast flange from your breast while pumping. Turn the breast pump off and with your finger, break the seal between your breast and breast flange, then remove breast flange from your breast.
- **ALWAYS** place pump on a horizontal flat surface when using product.
- **DO NOT** change batteries or open battery cover when pump is on.
- Remove batteries before storing product for long periods of time.

IMPORTANT: Can lead to material damage.

- **NEVER** put electrical product in water or a sterilizer, as you can cause permanent damage to the breast pump.
- **DO NOT** use antibacterial or abrasive cleaners/detergents when cleaning breast pump or breast pump parts.
- When plastic bottles and plastic breast pump kit parts are frozen, they can become brittle and may break when dropped.
- **DO NOT** over-fill bottles. Replace bottles before they overflow.
- Bottles and breast pump kit parts may become damaged if mishandled, e.g. dropped, overtightened, or knocked over.
- Always check bottles or breast pump kit parts prior to pumping. Discard breastmilk if bottles or breast pump kit parts are damaged.

SAVE THESE INSTRUCTIONS

2. INTENDED USE


The Ameda Finesse™ Breast Pump is intended to express and collect the mother's milk from the breasts of a nursing woman for the purpose of feeding the collected milk to a baby.

3. CONTRAINDICATIONS FOR USE


There are no known contraindications for use with this product.

4. PRODUCT DESCRIPTION

Finesse Breast Pump & HygieneKit Milk Collection System

Pump parts may vary

1. Pump Base
2. Power Light
3. Speed Dial
4. Suction Dial
5. Tubing Adapter Port
6. Power Port
7. AC Power Adapter (included)
8. AA Battery Compartment (batteries not included)
9. Tubing Adapter
10. Tubing
11. Adapter Cap
12. Diaphragm
13. Breast Flange (25.0mm)
14. Valve (Extra Included)
15. Locking Ring
16. Locking Disc
17. Bottle (120ml)

5. SET-UP INSTRUCTIONS

5.1 Before First Use: Cleaning The HygieniKit

Before first use, unless otherwise directed by your healthcare provider, prepare the product for use as follows:

1. **POP OFF** the adapter cap(s) from the flange(s).

SET ASIDE adapter cap(s), tubing and tubing adapter—you do not need to take these apart from each other. These 3 parts **DO NOT** need to be washed before first use or after each use—milk does not come into direct contact with these parts.

Take apart the remaining parts, the diaphragm(s), flange(s), valve(s) and milk storage bottles, including locking ring(s) and locking disc(s), since these are the parts that can touch the breast milk once you begin pumping.

DO NOT wash adapter cap, tubing and tubing adapter. Parts shown above.

ALWAYS wash hands well with soap and water.

CLEANING: Cleaning by boiling

IMPORTANT: Set aside tubing, adapter cap and tubing adapter. These parts **DO NOT** need to be washed.

DO NOT allow boiling parts to contact bottom of pot. It may cause damage to parts.

WARNING: DO NOT leave boiling parts unattended.

Remove parts carefully from boiling water to avoid burning. Boiled water and steam from boiled water can cause severe skin burns, so use caution when boiling parts in water and removing parts from boiling water.

1.

- **BOIL** a pot of clean water.

2.

- **WASH** your hands.

3.

- **REMOVE** the parts carefully with clean tongs.

4.


- **Air DRY** on a clean surface.

- **PUT** flange(s), diaphragm(s), valve(s) and bottle(s), locking ring(s) and locking disc(s) into pot for **20 MINUTES**.

5.2 Assembling the HygieniKit

1.

Once your pump parts are dry, wash your hands and assemble the HygieniKit (single or double).

ALWAYS wash hands well with soap and water.

2.

HOLD by the rim or thicker sides of the valve.

3.

PUSH the valve(s) onto the bottom of the flange(s).

Pump parts may vary.

4.

ATTACH and **TIGHTEN** the bottle(s) or attach milk storage bag(s).

5.

INSERT the diaphragm(s) into the top of the flange.
*It is important that the inside of the diaphragm is completely **DRY** to avoid moisture from forming in the tubing.*

6.

ATTACH the adapter cap(s) by pressing it onto the top of the flange(s) until you hear a pop.

7.

The tubing should still be **CONNECTED** to the adapter cap(s) and tubing adapter.
If you have taken these apart, attach tubing to adapter cap(s) and the tubing to the tubing adapter.

8.a.

ASSEMBLED HYGIENIKIT
 for single pumping.

8.b.

ASSEMBLED HYGIENIKIT
 for double pumping.

Pump parts may vary.

6. INSPECTING BEFORE PUMPING


IMPORTANT: Before using the Finesse™ Breast Pump and HygieneKit system, check the parts for cracks, chips, tears, or wear. They can affect the function of the kit. If parts need replacing due to defects, please call your local distributor or the location where you purchased the pump. Do not use damaged parts. For a listing of distributors in your country, please visit www.ameda.com.

7. OPERATING INSTRUCTIONS

7.1 Plug into Power Source

Choosing and checking your electrical source

AC POWER ADAPTER

WARNING: Use **ONLY** the AC power adapter recommended for the Finesse™ Breast Pump. AC power adapters with Australia, United Kingdom or Europe or other international plug configurations are not intended for use in the USA or Canada.

WARNING: Always inspect the AC power adapter for damage such as kinking or twisting of cords or frayed or exposed wires. If damaged, discontinue use immediately and call customer service at **1-866-992-6332 (1-866-99-AMEDA)**.

1. If batteries are in pump base, **REMOVE** batteries before using AC power adapter plug.
2. **PUSH** small end of AC power adapter into left side of the pump base as shown.
3. **PLUG** other end into AC outlet.

TIP - USE new premium batteries (such as Energizer, Duracell), insert correctly as shown.

BATTERIES

WARNING: DO NOT

install backwards, charge, put in fire, or mix with other battery types. May explode or leak causing injury. Replace all batteries at the same time.

1. **OPEN** battery compartment on the bottom of the pump base.
2. **INSERT** six new AA batteries as shown.
3. **CLOSE** the battery compartment.

Low Batteries

When batteries are low, the Pump Light will blink slowly and the pump will stop working. Turn the Suction dial to the off position. You may then replace all batteries or remove batteries and use the AC Power Adapter. Make sure batteries are not installed in pump while the AC Power Adapter is plugged into the side of the pump or while the Suction dial is in the on position.

7.2 Important to Know Before Breast Pumping

Learn how your breasts release milk and why Ameda's Custom Control™ mimics the way your baby breastfeeds.

Let-down & Milk Flow

When babies breastfeed, they latch on and start sucking quickly to help your body release hormones that signal your breasts to release your milk. This is called let-down. Once let-down occurs, your baby sucks longer and slower to help manage the milk flow that is spraying into your baby's mouth. With the Finesse™ breast pump, the separate dials allow you to mimic your baby's comfortable sucking pattern.

7.3 To Begin Pumping

Note: Place pump base on a horizontal flat surface before use.

1.

INSERT TUBING

adapter into right side of pump base, then turn slightly to the right while pushing in firmly.

a.

For **SINGLE PUMPING:** Remove one tube from the tubing adapter and close one side.

b.

For **DOUBLE PUMPING:** Both tubes should be attached to the tubing adapter.

2.

CENTER flange on breast to create seal.

NOTE: Make sure upper body is upright.

3.

TURN SUCTION

dial to right to turn on. Turn up to your highest comfortable suction level. If discomfort, turn suction down.

4.

TURN SPEED

dial all the way to the right to highest speed to start milk flow.

5.

Once milk flows, **TURN SPEED** dial to left to slow pump speed to help drain milk. Try half way, adjust as needed.

6.

Once milk flow slows, try for another let-down. **TURN SPEED** all the way up, then slow speed when milk is flowing. Repeat as needed.

7.

TURN SUCTION

dial to left to turn off.


8.

After pump is off, if necessary, **INSERT** a finger between breast(s) and breast flange(s) to break suction.


7.4 How to Fit Your Flange

Check breast flange fit, if pumping causes discomfort a larger or smaller breast flange size may be needed.


GOOD FIT -

During pumping, your nipple moves freely in the breast flange tunnel. You see space around the nipple. Not much areola is drawn into the tunnel with the nipple.

TOO SMALL -

During pumping, some or all of your nipple rubs against the sides of the breast flange tunnel.

TOO LARGE -

During pumping, more areola is drawn into the breast flange with your nipple. Your areola may rub against the side of the breast flange tunnel.

Additional flange sizes sold separately.

7.5 Storing Pumped Milk in Bottles

Make sure to not overfill bottles. If you plan on freezing bottles, fill only up to $\frac{3}{4}$ full, as milk expands with freezing. After milk has been pumped into Ameda storage bottles, close by placing locking disc on top of bottle (make sure circular raised middle part of disc is facing up). Tighten locking ring onto bottle. Label bottle with date and time that milk was pumped. Use storage chart provided as a guideline.

Guidelines for Storing and Thawing Human Milk

Storage Time for Human Milk*	Deep Freeze	Refrigerator Freezer	Refrigerator	Cooler with Frozen Ice Packs	Room Temperature	
	(0°F/-18°C)	(variable 0°F/-18°C)	(39°F/4°C)	(59°F/15°C)	(66°F-72°F) (19°C-22°C)	(72°F-79°F) (22°C-26°C)
Fresh	12 Months	3-4 Months	8 Days	24 Hours	6-10 Hours	4 Hours
Frozen, Thawed in Fridge	Do Not Refreeze	Do Not Refreeze	24 Hours	Do Not Store	4 Hours	4 Hours
Thawed, Warmed, Not Fed	Do Not Refreeze	Do Not Refreeze	4 Hours	Do Not Store	Until Feeding Ends	Until Feeding Ends
Warmed, Fed	Discard	Discard	Discard	Discard	Until Feeding Ends	Until Feeding Ends

*Storage times may vary for premature or sick babies. Sources: Jones, F. *Best Practice for Expressing, Storing and Handling Human Milk in Hospitals, Homes and Child Care Settings*. Raleigh, NC: Human Milk Banking Association of North America, 2011.

8. CLEANING THE HYGIENIKIT

After Every Use: Cleaning The HygieneKit

1. **POP OFF** the adapter cap from the flange.
2. **SET ASIDE** adapter cap, tubing and tubing adapter— you do not need to take these apart from each other. These 3 parts **DO NOT** need to be washed before first use or after each use— milk does not come into direct contact with these parts.

DO NOT wash adapter cap, tubing and tubing adapter. Parts shown above.

Take apart the remaining parts, the diaphragm(s), the flange(s), valve(s) and milk storage bottles, including locking ring(s) and locking disc(s), since these are the parts that can touch the breast milk once you begin pumping.

ALWAYS wash hands well with soap and water.

CLEANING: *Cleaning with warm, soapy water*

After every use, unless otherwise directed by your healthcare provider, prepare the product for use as follows:

IMPORTANT: Set aside tubing, adapter cap and tubing adapter. These parts **DO NOT** need to be washed.

1.

RINSE any milk from pump parts with cool water.

2.

Hand **WASH** with warm, soapy water, the flange(s), diaphragm(s), valve(s) and (if applies) bottle(s), locking ring(s) and locking disc(s). Wash valve(s) and diaphragm(s) gently. Do not insert anything into valve(s) while cleaning.

3.

RINSE with warm, clean water the flange(s), diaphragm(s), valve(s) and (if applies) bottle(s), locking ring(s) and locking disc(s).

4.

Air **DRY** each part on a clean surface.

IMPORTANT: the flanges and milk storage containers can also be washed on upper rack of a standard, household dishwasher. **DO NOT** use solvents or abrasives. Allow to air dry on a clean surface. **DO NOT** put parts in Steam Bottle Sterilizers or Microwave Sterilizer Bags.

Cleaning Finesse Pump: Make sure it is unplugged prior to wiping the Finesse Breast Pump and power adapter with clean, damp cloth.

9. MAINTENANCE

The breast pump base does not require service or maintenance except for changing batteries, if applies.

9.1 Breast Pump Valve Check

Inspect the valves for tears or holes. Replace any damaged valves prior to pumping. Wash valves per cleaning section. Do not insert bottle brush or something similar into valves. We recommend replacing valves approximately every 8-12 weeks but this can vary depending on cleaning and use.

9.2 Storing and Transporting Finesse/HygieniKit

Disconnect HygieniKit and power adapter from pump before transporting or storing the product. When storing product, store in conditions as described in these instructions (see Section 13). Remove batteries before storing product for long periods of time. Handle with care when transporting product.

9.3 Power Adapter Care


To care for your adapter be sure to follow these storage directions.


DO NOT wrap the power cord around the power adapter body.

CORRECT way to bundle cord.


DO NOT pull the cord to unplug the power adapter.


CORRECT way to plug in the power adapter.

WARNING: Prior to use, always inspect power adapter for damage or exposed wires. If damage and/or exposed wires are found, **DO NOT** use power adapter and call Ameda Customer Service ParentCare Team at **1-866-992-6332 (1-866-99-AMEDA)**.

Twists, kinks or bends

Has bulges

Exposed wires

10. TROUBLESHOOTING

SLOW MILK FLOW, UNEVEN OR NO SUCTION	<p>VALVE CHECK Inspect the valves for tears or holes. Replace any damaged valves prior to pumping. Wash valves per cleaning section. Do not insert bottle brush or something similar into valves. We recommend replacing valves approximately every 8-12 weeks but this can vary depending on cleaning and use.</p> <p>CONNECTION CHECK Check assembly is correct and connection points are securely attached.</p> <ul style="list-style-type: none">• Tubing adapter to tubing port• Tubing adapter to tubing <p>Check adapter cap is attached to breast flange. It may be useful to remove adapter cap and reattach to breast flange.</p> <p>DIAPHRAGM AND FLANGE/BOTTLE CHECK Inspect the diaphragms for tears or holes. Replace any damaged diaphragms prior to pumping. Check flanges, adapter cap or bottles for chips, cracks, breaks, or other damage. Replace any damaged flanges/adapter caps/bottles prior to pumping.</p> <p>PUMP BASE MOTOR CHECK Slowly and steadily increase the suction dial from the OFF position all the way to the Maximum position. The motor sound should increase as you turn the dial up. Slowly and steadily increase the speed dial from Slow to Maximum. The motor should go faster as you turn the dial up.</p> <p>BREAST FLANGE FIT CHECK</p> <ul style="list-style-type: none">• Good Fit: Nipple is moving freely in tunnel of flange.• Too Small: Nipple is rubbing against tunnel of flange.• Too Large: Areola is being pulled into tunnel of flange.
AUDIBLE, HISSING NOISE OR AIR LEAKAGE	<p>VALVE CHECK Inspect the valves for tears or holes. Replace any damaged valves prior to pumping. Wash valves per cleaning section. Do not insert bottle brush or something similar into valves. We recommend replacing valves approximately every 8-12 weeks but this can vary depending on cleaning and use.</p> <p>CONNECTION CHECK Check assembly is correct and connection points are securely attached.</p> <ul style="list-style-type: none">• Tubing adapter to tubing port• Tubing adapter to tubing <p>Check adapter cap is attached to breast flange. It may be useful to remove adapter cap and reattach to breast flange.</p>

<p>POWER SOURCE</p>	<p>AC POWER ADAPTER Finesse™ is optimized for use with the AC Power Adapter. Make sure the AC power adapter is plugged all the way into the side of your pump base. Make sure there are no batteries installed in the pump base. Make sure you rotate the suction dial up to turn on the pump. Try your pump in a different electrical outlet.</p> <p>AA BATTERIES Be sure the batteries have been inserted correctly as shown inside the breast pump base's battery compartment. Use new premium batteries (such as Energizer, Duracell.) Do not use batteries with AC power adapter plugged into side of pump base. If problem persists, you can test that the batteries work in another battery powered appliance.</p>
<p>PUMP LIGHT INDICATOR</p>	<p>LOW POWER: The light blinks slowly and the pump stops working. Turn the Suction dial to the off position. If using batteries, replace batteries (see section 7.1) or remove all batteries and use AC power adapter. If using AC power adapter and a slow blink occurs, check it according to Troubleshooting Power Source instructions.</p> <p>HIGH POWER: The light blinks quickly and the pump stops working. Turn the Suction dial to the off position. Check AC power adapter according to Troubleshooting Power Source instructions, make sure to unplug and re-connect the AC power adapter into the side of the pump and check all connections are correct. If using batteries, check them according to Troubleshooting Power Source instructions.</p> <p>ERROR: The light changes to a double fast blink and pump stops working. Turn the Suction dial to the off position. Unplug and re-connect the AC power adapter into the side of the pump and check all connections are correct, then return to operation.</p>

If there are still questions or issues please call Ameda's certified customer service ParentCare team at **1-866-992-6332 (1-866-99-AMEDA)**.

11. WARRANTY

THE AMEDA WARRANTY IS NONTRANSFERABLE BETWEEN INDIVIDUALS. THE AMEDA FINESSE ELECTRIC BREAST PUMP (“PRODUCT”) IS WARRANTED TO THE ORIGINAL USER – WHETHER INSTITUTIONAL OR INDIVIDUAL (“FIRST USER”) – ONLY.

From the date of the Individual User's purchase of this product, Ameda Company, Inc. (“Ameda”) warrants the Product to the Individual User against defects in material or workmanship for a period of two years for the pump mechanism and 90 days on other pump components except for breast pump valves. Ameda's sole obligation under this limited warranty shall be to repair or replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Repair or replacement under this limited warranty is the sole and exclusive remedy of the First User. **PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE.** This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty service contact Ameda at 1.866.99.AMEDA (26332).

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR ANY OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

NON-APPLICABILITY OF WARRANTY

In addition to any other limitations on the warranty discussed above, the warranty shall expire and be of no force or effect upon the occurrence of any of the following:

- (i) the cause of damage to the Product is due to improper handling or use.
- (ii) opening of the pump motor, modifications, self-repair or use outside the instructions in this manual voids warranty.
- (iii) the Product is not properly maintained by the Individual User as determined by Ameda.
- (iiii) non-Ameda approved accessories or spare parts are used with the Product as determined by Ameda.

If any of the above occurs, Ameda shall not be responsible for damage to the Product, parts thereof or injuries arising therefrom, either directly or indirectly. The warranty terms set forth above replaces all prior warranties with respect to the purchased Product, whether in writing or otherwise.

Outside the United States: Call your local distributor or location where you purchased the Product. For a listing of distributors in your country, please visit www.ameda.com.

12. TECHNICAL DATA

AC POWER ADAPTER (included)

North America (4100003), Australia (4100007), United Kingdom (4100005), Europe (4100004) Plug Configurations available, Direct plug in AC power adapter. Output: 9.0 VDC/1.3A, Input: 100-240VAC 50/60HZ, 0.4A MAX

BATTERIES (NOT INCLUDED): Six AA batteries, 1.5 V alkaline

PRODUCT TECHNICAL INFORMATION

PROTECTION CLASS:

Class II/Double Insulated

APPLIED PART:

Type BF Applied Part

INGRESS CLASSIFICATION:

IP22

Solid Particle/Liquid Ingress Classification/
Protection against fingers and dripping water.

COMPLIANCE:

UL CERTIFIED E493585 MEDICAL – GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN
ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD1 (2012) AND CAN/
CSA-C22.2 NO. 60601-1 (2008) + (2014) AND IEC 60601-1-6 (2010) + AMD1
(2013) AND IEC 62366 (2007) + AMD1 (2014) AND IEC 60601-1-11:2015

CE 2797

This product meets the applicable European regulations
and is under the supervision of a Notified Body.

Authorized Representative in the European Community

SUCTION:

0 to 214.5 mmHg–Dual HygieniKit Configuration

0 to 297 mmHg–Single HygieniKit Configuration

SPEED:

30 to 60 cycles per minute

OPERATING LIFE:

560 hours

EMC Information

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC provided in these Instructions for use.

Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

Warning: The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of

the Ameda Finesse Breast Pump as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the Ameda Finesse Breast Pump.

Warning: The Ameda Finesse Breast Pump should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Ameda Finesse Breast Pump should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Ameda Finesse Breast Pump is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Finesse Breast Pump should assure that it is used in such an environment.
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ameda Finesse Breast Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The Ameda Finesse Breast Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Ameda Finesse Breast Pump is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Finesse Breast Pump should assure that it is used in such an environment.	
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ±1 kV line(s) to line(s)	± 0.5 kV, ±1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0,5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 s	<5% U _T (>95% dip in U _T) for 0,5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ameda Finesse Breast Pump requires continued operation during power mains interruptions, it is recommended that the Ameda Finesse Breast Pump be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The Ameda Finesse Breast Pump is intended for use in the electromagnetic environment specified below. The customer

or the user of the Ameda Finesse Breast Pump should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Ameda Finesse Breast Pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an

electromagnetic site survey should be considered. If the measured field strength in the location in which the Ameda Finesse Breast Pump is used exceeds the applicable RF compliance level above, the Ameda Finesse Breast Pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ameda Finesse Breast Pump.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Ameda Finesse Breast Pump

The Ameda Finesse Breast Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ameda Finesse Breast Pump can help prevent electromagnetic interference by maintaining a minimum

distance between portable and mobile RF communications equipment (transmitters) and the Ameda Finesse Breast Pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


PRODUCT DISPOSAL

This product contains electrical and electronic equipment.


1. The internal components of this product may contain hazardous materials. Dispose of this product according to local or regional waste administration systems and regulations.
2. Do not dispose of as unsorted municipal waste.
3. For further information regarding the environmental performance of this product please visit our website www.ameda.com.

13. STORAGE AND OPERATIONAL CONDITIONS

Operational conditions

Transport and Storage conditions

Operation, Transport and Storage conditions

Fragile

Protect from sunlight

Keep Dry

14. TRAVEL OR INTERNATIONAL USE

WARNING: Use **ONLY** the AC power adapter recommended for the Finesse™ Breast Pump.

AC power adapters with Australia, United Kingdom, Europe or other international plug configurations are not intended for use in North America.

If you will be traveling outside of your country where you received your Finesse™ breast pump you will need to use batteries or a Finesse™ AC power adapter specific to that area of travel in order to use your breast pump while traveling.

WARNING: You could cause serious harm to yourself or damage the breast pump.

15. AMEDA ACCESSORIES AND SPARE PARTS

SOME PUMPS MAY ALSO INCLUDE SOME OF THESE ACCESSORIES:

You may also be able to order these additional accessories or spare parts online, at a local Ameda retailer or Ameda insurance covered provider, visit www.ameda.com for more information.

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Tapa del adaptador

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To register your pump.

Please visit AMEDA.com and input the serial number (SN) from the bottom of your pump.

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Mothers Choice Products Ltd.
2133-11871 Horseshoe Way
Richmond, BC V7A5H5
1.800.604.6225

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Midmed Pty Ltd.
Unit 4, 62 Borthwick Avenue
Murarrie, QLD 4172
61 (7) 3348 9155

Distribution in other countries: For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit www.ameda.com

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The Ameda Logo, Ameda, Cool'N Carry, CustomControl, Flexishield, Finesse™, HygieniKit, Store'N Pour, and/or registered trademarks of Ameda, Inc. in the United States and other countries.

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CE 2797

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Manufactured for:
Ameda, Inc.
485 Half Day Road Ste 320
Buffalo Grove, IL 60089
USA