

*Marina Home* automatic recliner



The *Marina Home* automatic recliner is an advancement on DECAM's original product which has demonstrated outstanding reliability and durability through use by hundreds of customers in both private and public hospital settings.

The Marina Home automatic recliner embodies a novel concept, namely an automatic reclining chair for indoor use which is also particularly useful for patient mobility.

○ ● ● ● ● positions

**reclined position**

thanks to two independent motors, the Marina Home automatic recliner features multiple reclining positions controlled with an easy to use, hand held controller



**sleeper**

with a simple touch of a button and without any effort, the Marina Home automatic recliner can be converted into a bed which provides an exceptional level of comfort



**upright position**

the armrests can be raised up to transfer the patient to or from a bed without any effort. The armrests are movable in any chair position



**Lifting position**

Patients often suffer from reduced mobility. Helping patients get out of the chair will not only make life easier for the patients but will also reduce caregiver injuries



## ○ ● ● ● ● standard features

### push bar

an innovative design that facilitates the movement and rotation of the chair

### urethane arm caps

provides users with a comfortable grip and protection for their arms

### moisture barrier

offers protection from fluids such as blood, sweat and urine, with minimum maintenance..



### moveable armrests

the armrests can be lifted to transfer the patient to or from a bed without any effort



### controller

easy to use, hand held controller

### pull out food rest

for additional support for legs when seated

### zero turning radius

eases mobility in small spaces

### casters

silent and resistant double casters with individual locks

## ○ ● ● ● ● optional features

### battery

a rechargeable battery provides the user with a much greater level of autonomy.

### oxygen tank holder

attaches to the back of the chair to carry an oxygen cylinder.

### swing away table - side table - front table

3 different types of tables to satisfy the different needs of the user.

### head rest

the adjustable head rest is removable, and provides greater comfort to the user.

### side head rest

supports patient's head by stabilizing it in an upright position preventing side to side lateral leaning.

### upholstered armrests

for an elegant style with the extreme comfort of the upholstery.

### memory foam

provides a better comfort relieving pressure points, preventing pressure sores, etc...

### IV Pole holder

enable patients to take IV fluids with them.

○ ● ● ● ● technical specifications

*Marina Home automatic*

dimensions	overall		seat		arm height		full recline	
width	62,5 cm	24,6"	50 cm	19,7"	70 cm	27,5"		
depth	75 cm	29,9"	50 cm	19,7"			170 cm	66,9"
height	127 cm	50"	52 cm	20,5"				
weight	70Kg	154 lb						
max load	150 Kg	330 lb						

*Marina Home automatic (wider seat)*

dimensions	overall		seat		arm height		full recline	
width	68 cm	26,7"	56 cm	22"	70 cm	27,5"		
depth	75 cm	29,5"	50 cm	19,7"			170 cm	66,9"
height	127 cm	50"	52 cm	20,5"				
weight	75 Kg	165 lb						
max load	150 Kg	330 lb						

\* Bariatric and custom sizes available upon request.

○ ● ● ● ● fabrics

The Decam textile collection, based on Spradling® Diamante and Silvertex patterns is designed specifically to meet the strictest requirements of the healthcare environment by creating effective barriers against germs, abrasion and stains.



\* visit our website for a full list of available colors

Our standard collection meets the following features:

Abrasion resistance	Urine protection	Sweat Protection	100% Polyester
Antibacterial Protection	Blood Protection	Mildew Resistant	Sulphide Stain Resistant
Antifungal Protection	Stain Protection		

This brochure provides an overview of the main and benefits of the DECAM product. For a full list of products available, including specifications, fabrics, and finishes, please visit [www.decamlux.com](http://www.decamlux.com). For pricing enquiries contact your local DECAM sales representative.

Decam Llits Abatibles S.L. - Bujonis Industrial, Parcel·la 13 - 17220 Sant Feliu de Guíxols (Girona) - Spain - Phone: (+34) 972 32 04 00  
[www.decamlux.com](http://www.decamlux.com) - [info@decamlux.com](mailto:info@decamlux.com)