

# RESPONDER BANDAGE

The Responder Bandage is designed for quick and easy application by caregivers to provide effective multi-functional treatment when time is of the importance.

The Responder Bandage's sterile, non-adherent pad applies pressure to any site, is easily wrapped and secured. The Responder Bandage will restrict blood flow to help in the treatment of Traumatic Bleeds.

The Responder Bandage is available in three sizes:  
4in", 6in", 8in" width

4" bandage

Pad wt: three layer, 280g+280g+110g/m<sup>2</sup>

Pad size: 10x18cm

Bandage wt: 80-85g/m<sup>2</sup>

Width of bandage: 10cm

Natural length: 1.5m

Elastic len th: 4.5m

Double bag rolling packing, inner bag vaccum packing.

6" bandage

Pad wt: three layer, 280g+280g+110g/m<sup>2</sup>

Pad size: 15x18cm

Bandage wt: 80-85g/m<sup>2</sup>

Width of bandage: 15cm

Natural length: 1.5m

Elastic len th: 4.5m

Double bag rolling packing, inner bag vaccum packing.

