

526 Series

CLASSIC OVERNIGHTER CHAIR

RATED UP TO 136.08 KG

Part #526



Three positions available for optimal function and comfort.

Shown with optional oak armrests.

STANDARD FEATURES

Three positions

Steel and wood construction

Caster wheels (rear locking)

Open design to allow for body mobility

The open-concept design of the Classic Overnighter provides stylish daytime seating and spacious sleeping for overnight guests. One of our very first products, the Classic Overnighter features time-tested engineering with superior durability for years of use. Medical-grade vinyls clean easily and come in a wide range of colors to match any interior color scheme. Keep loved ones close and comfortable with the Classic Overnighter.

Caster wheels for easy conversion



Z-spring construction



Open-concept design




STANDARD VINYL SELECTION



Natural
T09



Aquamarine
T19



Bonnie Blue
T21



Iced Mint
T26



Tearose
T43



Colonial Blue
T45




Gray
T56



Cranberry
T96

Printed color swatches may not be an exact match to actual fabrics and are enlarged to show detail.

SPACE REQUIREMENTS



SPECIFICATIONS

Weight capacity	136.08 KG
Seat height/width	48.26cm/63.5cm
Seat depth	53.34cm
Arm width	6.35cm
Overall height	106.68cm
Overall width	78.74cm
Overall length (seated position)	86.36cm
Overall length (sleep position)	187.96cm
Optional wood arm cap width	7.62cm



Steel frame components are powdered-coated to resist corrosion.

OPTIONS

CAL 133 fire barrier (3)*

Vinyl upgrades

Champion, graded vinyl selections**

Customer's own material (COM)

PVC-free material

*Contact Midmed for a list of fire-tested vinyls.

**For available patterns and colors for graded vinyl selections, please contact your Midmed sales representative.



VINYL UPGRADES