

PORTABLE PROCEDURE TABLE (PPT) PATIENT TRANSPORT AND OPERATION PLATFORM

Field tested and outfitted with the most advanced materials available, this unit can serve a multitude of purposes. The PPT quickly assembles into a stable patient aid platform for medics and doctors. The device is versatile and easy to transport for fast changing scenarios.

INCLUDES :

- NATO Wheeled Litter Carrier
- Heavy duty transport container
- IV Arm Board (with removable pad)
- LED Exam Light
- IV Pole
- Mayo Tray (removable tray)
- Accessory Pole (electronics)
- Utility storage catcher
- Oxygen Rack, D Cylinder (optional)

*This item can be customized to fit your needs.

FEATURES :

- Heavy duty wheels that include solid inserts to prevent flat tires.
- Quick release pins to attach and detach tires in seconds.
- Fully adjustable rail clamps for perfect position, height and angle of accessories.

PPT PACKAGE FULL (includes everything shown)

-Part # 2221044

PPT PACKAGE (Minus Litter)

-Part # 2221040

PPT CONVERSION KIT

-Part #2221043

PPT HEAVY DUTY TRANSPORT CASE

-Part #2221110



BENEFITS :

- Easy to clean.
- Assembles/Disassembles in minutes.
- Adjusts to fit multiple sizes/shapes of blackboards and litters.
- Entire unit becomes mobile with press of pedal.

SPECIFICATIONS :

- Cold rolled steel frame
- Anti-corrosion base layer & powder coating paint process
- Stainless steel and anodized aluminum fittings.
- Case Size: 41.8" x 37.6" x 24.0"
- Total Weight: 160 lbs
- Weight limit: Current max. test level, 500 lbs.