


CO.ES480A
ELECTRIC
STRETCHER



CO.ES480A ELECTRIC STRETCHER

Cogent offers this Electric Stretcher which provides excellent safety, comfort, and quality both to the patient & care provider. This stretcher can do six positions electrically to enhance the comfort of patients. Patients' can adjust positions by themselves without asking for help from nursing staff. The stretcher provides durable quality and a 318kg weight capacity to accommodate patients of different sizes. It is also very suitable to be used as a temporary hospital bed.





STANDARD FEATURES

- Integrated storage tray with oxygen tank trough which is made with ABS plastic
- Electric ascending and descending
- Electric Fowler position
- Retractable steering fifth-wheel
- Integrated with conductive chain
- Structure with four-section bed board
- Foot controlled brake / steer pedal for four wheels in two ends of stretcher
- Drainage Bag clasp x 2 & Urinary Bag clasp x 2
- IV pole inserted pole x 5
- Electric Backboard articulation
- Collapsible side rail (PP plastic)
- Two stage movable IV pole x 1
- Foot end push bar
- Jack / Hook upright Oxygen tank holder
- Protective corner bumpers x4
- Electric knee gatch
- CPR device
- Upright oxygen tank holder x 1

SPECIFICATIONS

- Overall Length: 213cm
 - Overall Width: 81cm
 - Height Range (from ground to litter top):
(Electric) High: 93cm ± 2
(Electric) Low: 53cm ± 2
 - Weight Capacity: 318kg
 - Litter Positioning:
(Electric) Backrest: 0-80° ± 2°
(Electric) Trend. / Reverse Trend: 0-23° ± 2°
Electric knee gatch: 0-45° ± 2°
 - Mattress Base: 193 x 64cm
 - Mattress Dimensions: 195 x 66 x 8cm
 - Siderail: 141 x 48cm
 - Caster Diameter: 20cm
 - Power supply 110V-240V
- Two years warranty for stretcher under normal use after date of shipping

Ordering Code: CO.ES480A



www.midmed.com.au



1300 643 633