



TECHLEM
STRETCHERS

PRODUCT CATALOGUE



Techlem Stretchers delivers tailored patient handling solutions to hospital sites and clinics around the world.



TABLE OF CONTENTS

ABOUT US

1

ANATOMY OF THE TECHLEM STRETCHER

2 - 3

TRANSPORT & GENERAL USE 5000

4

BARIATRIC 5000BR

5

EXTENDED STAY 5000W

6

FIXED HEIGHT 1100 / 1000N

7

PATIENT LIFT 7600

8

TRAUMA & X-RAY 5000XR

9

FOLD RAIL 4500

10

COMPACT 4500N / 4500NS

11

GYNAECOLOGY 4500GYN

12

CADAVER & BARIATRIC CADAVER

13

COLOR GUIDE

14

I.V. STAND 8012A

15

OEM REPLACEMENT MATTRESSES

16 - 17

CUVERRO SHIELD BY AEREUS

18 - 19

ACCESSORIES GUIDE

20 - 21


ABOUT US

For 60 years, Techlem has provided patient handling solutions to hospitals and clinics around the world. From our first stretcher prototype to our most recent full hospital outfit, we'vestrived for value, efficiency, and care. Our products are proudly North American made with the highest quality components, delivering industry best costs of ownership.

Our passionatestaff work tirelessly to produce durable, high quality products with customer service to match. We create unique, user-friendly catalogs for our supply chain partners, bridging the gap between purchasers and end-users. Why choose to partner with Techlem? Because of our commitment to finding the most useful, cost effective, and personalized catalog program for every customer. Why choose Techlem stretchers? Because in quality control, durability, cost of ownership, and value over time, they are the envy of the industry.


We value our 60 years of partnerships. We have a dedicated staff with more than 100 combined years' experience in the industry, and nearly 100 years' combined service to Techlem.

Our unshakably caring corporate culture is our pride. We extend this culture in well established, long lasting relationships with customers in the Americas, Europe,Asia, and Australia. We maintain a commitment to a same or next-day response time, and to a highly competitive turnaround time at every step of the purchase process, from quotation to invoicing. We offer personal, conscientious customer service unmatched in the industry.



TRANSPORT & GENERAL USE 5000

The premier Techlem stretcher for transport and general use features an integrated oxygen tank holder and patient storage, 5th Wheel technology for steering assist, and more than 20 customizable options for design and functionality.



DIMENSIONS
TOTAL DIMENSIONS: 81" x 32 ¹ / ₂ "
LITTER DIMENSIONS: 81" x 29 ¹ / ₂ "
FULL LENGTH SIDE-RAILS: 75" x 12 ¹ / ₂ "
HIGH DENSITY FOAM MATTRESS: 76" x 28 ¹ / ₂ " x 3"
*WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT
HEIGHT ADJUSTMENT: 22" TO 35 ¹ / ₂ "
HANDS-FREE TRENDLENBERG: +18° TO -18°
PNEUMATIC ASSIST FOWLER: 0° TO 90°
8" ALL-SURFACE, LOCKING CASTERS
ACCU-STEER 5 TH WHEEL
POSI-LOCK BRAKING SYSTEM
DROP UNDER SIDE-RAILS

MATTRESS
STANDARD: HIGH DENSITY FOAM MATTRESS WITH CONDUCTIVE COVER
OPTIONAL: VISCOLASTIC, PRESSURE SENSITIVE MATTRESS, WITH FOUR-WAY, ELASTATEX STRETCH COVER



*See page 14 for color accent options

BARIATRIC 5000BR

Reliable, nimble, and equipped with Techlem's 5th Wheel technology and high tensile components.



DIMENSIONS
TOTAL DIMENSIONS: 81" X 35 ¹ / ₂ "
LITTER DIMENSIONS: 81" X 32"
FULL-LENGTH SIDE RAILS: 75" X 16"
HIGH-DENSITY FOAM MATTRESS: 76" X 32"
*WEIGHT CAPACITY: 750LBS / 341KG

MOVEMENT
HEIGHT ADJUSTMENT: 22" TO 32"
HANDS-FREE TRENDLENBERG: +12° TO -12°
ACCU-STEER 5 TH WHEEL
POSI-LOCK BRAKING SYSTEM
PNEUMATIC ASSIST FOWLER: 0° TO 90°
8" ALL-SURFACE, LOCKING CASTERS
DROP UNDER SIDE-RAILS
OPTIONAL FOUR-SIDED ACCESS

MATTRESS
STANDARD: HIGH DENSITY FOAM MATTRESS WITH CONDUCTIVE COVER
OPTIONAL: VISCOLASTIC, PRESSURE SENSITIVE MATTRESS, WITH FOUR-WAY, ELASTATEX STRETCH COVER

EXTENDED STAY 5000W

Versatile, functional, and comfortable. Our extended stay stretcher is equipped with extra-wide frame, DROP UNDERside-rails, and Techlem’s distinguished Viscolastic mattress technology.



DIMENSIONS
TOTAL DIMENSIONS:81" x 35 1/2"
LITTER DIMENSIONS:81" x 32"
MATTRESS DIMENSIONS:76" x 32"
WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT
HANDS-FREE TRENDLENBERG: +15° TO -15°
ACCU-STEER 5 TH WHEEL
POSH-LOCK BRAKING SYSTEM
PNEUMATIC ASSIST FOWLER: 0° TO 90°
8" ALL-SURFACE, LOCKING CASTERS
DROP UNDER SIDE-RAILS
OPTIONAL FOUR-SIDED ACCESS

MATTRESS
PRESSURESENSITIVE FOAM MATTRESS WITH VISCOLASTIC, FOUR-WAY STRETCH COVER



*See page 14 for color accent options

FIXED HEIGHT 1100 / 1000N

Techlem’s value stretcher option, our fixed height stretcher is equipped with many of the same features as our transport stretcher, without the adjustable height.



1100 DIMENSIONS
TOTAL DIMENSIONS:81" X 32 1/2"
MATTRESS DIMENSIONS:75" X 28 1/2"
WEIGHT CAPACITY: 700LBS / 318KG

1100N (NARROW)DIMENSIONS
TOTAL DIMENSIONS:81" X 29"
MATTRESS DIMENSIONS:75" X 25"
WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT
HANDS-FREE TRENDLENBERG: +18° TO -18°
ACCU-STEER 5 TH WHEEL
POSH-LOCK BRAKING SYSTEM
PNEUMATIC ASSIST FOWLER: 0° TO 90°
8" ALL-SURFACE, LOCKING CASTERS
DROP UNDER SIDE-RAILS

MATTRESS
FULL WIDTH, HIGH DENSITY FOAM MATTRESS WITH A CONDUCTIVE COVER



PATIENT LIFT 7600

Simple and efficient in design,our patient lift and transfer stretcher is adaptable to your clinical environment. Handcrafted from high grade materials, the 7600 is highly durable, maintaining its form and function over decades of use.



FRAME DIMENSIONS
OUTSIDE DIMENSIONS: 77" x 31"
INSIDE DIMENSIONS: 74" x 24 1/2"
WEIGHT CAPACITY: 300LBS / 138KG

MOVEMENT
HEIGHT ADJUSTMENT: 30" TO 42"
5" ALL-SURFACE, LOCKING CASTERS
5 TH WHEEL FOR STEER ASSIST
CRANK HANDLE FOR EASE OF USE

SPECIAL FEATURES
HIGH DENSITY PATIENT TRANSFER SHEET
HIGH GRADE STEEL SHEET BUCKLES
TELESCOPING I.V. POLE INCLUDED

TRAUMA & X-RAY 5000XR

For use in trauma and X-Ray environments. Equipped with a fiber resin film surface and high density radiolucent mattress. A Techlem total-view design, with full surface accessibility for cassette placement.



DIMENSIONS
TOTAL DIMENSIONS: 81" X 32 1/2"
LITTER DIMENSIONS: 81" X 28 1/2"
MATTRESS DIMENSIONS: 75" X 29" X 3"
DROP-RAIL SIDE-RAILS: 60 1/2" X 12 1/2"
*WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT
HEIGHT ADJUSTMENT: 22" TO 35 1/2"
HANDS-FREE TRENDLENBERG: +18° TO -18°
PNEUMATIC ASSIST FOWLER: 0° TO 90°
DROP UNDER SIDE-RAILS
8" ALL-SURFACE, LOCKING CASTERS
ACCU-STEER 5TH WHEEL
POSI-LOCK BRAKING SYSTEM

MATTRESS
STANDARD: HIGH DENSITY FOAM
MATTRESS WITH CONDUCTIVE

*See page 14 for color accent options

FOLD RAIL 4500

Reliable, mobile, and accessible. Our fold-rail stretcher is built on the same, durable chassis as the 5000, and is equipped with collapsing side-rails for unobstructed, four-sided access for medical professionals.

4500 DIMENSIONS
TOTAL DIMENSIONS: 81" x 32 1/2"
LITTER DIMENSIONS: 81" x 28 1/2"
MATTRESS DIMENSIONS: 75" x 29" x 3"
DROP-RAIL SIDE-RAILS: 60 1/2" x 12 1/2"
WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT
HEIGHT ADJUSTMENT: 22" TO 35 1/2"
HANDS-FREE TRENDLENBERG: +18° TO -18°
PNEUMATIC ASSIST FOWLER: 0° TO 90°
FOLD DOWN SIDE-RAILS
8" ALL-SURFACE, LOCKING CASTERS
ACCU-STEER 5TH WHEEL
POSH-LOCK BRAKING SYSTEM

MATTRESS
STANDARD: HIGH DENSITY FOAM MATTRESS WITH CONDUCTIVE COVER
OPTIONAL: VISCOLASTIC, PRESSURE SENSITIVE MATTRESS, WITH FOUR-WAY, ELASTATEX STRETCH COVER



*See page 14 for color accent options

COMPACT 4500N / 4500NS

Highly adaptable, our compact stretcher is suited to your small-space, clinical environment. Provides many of the same features as the 4500, including collapsing side-rails. Available in narrow and short & narrow.

4500N (NARROW) DIMENSIONS
NARROW
TOTAL DIMENSIONS: 81" X 29"
LITTER DIMENSIONS: 81" X 25"
MATTRESS DIMENSIONS: 78" X 25" X 3"
WEIGHT CAPACITY: 500LBS / 227KG

4500NS (SHORT & NARROW) DIMENSIONS
SHORT & NARROW
TOTAL DIMENSIONS: 74" X 29"
LITTER DIMENSIONS: 74" X 25"
MATTRESS DIMENSIONS: 74" X 25" X 3"
WEIGHT CAPACITY: 500LBS / 227KG

MOVEMENT
HEIGHT: 22" TO 35 1/2"
HANDS-FREE TRENDLENBERG: +18° TO -18°
ACCU-STEER 5TH WHEEL
POSH-LOCK BRAKING SYSTEM
PNEUMATIC ASSIST FOWLER: 0° TO 90°
8" ALL-SURFACE, LOCKING CASTERS
DROP UNDER SIDE-RAILS

MATTRESS
STANDARD: HIGH DENSITY FOAM MATTRESS WITH CONDUCTIVE COVER
OPTIONAL: VISCOLASTIC, PRESSURE SENSITIVE MATTRESS, WITH FOUR-WAY, ELASTATEX STRETCH COVER



GYNAECOLOGY 4500GYN

For use in obstetrics and gynaecology, the 4500GYN stretcher is equipped with a two-piece sectional, Viscolastic mattress. Versatile, it boasts supine to full-chair positioning with pneumatic assist for ease of use.



DIMENSIONS

TOTAL DIMENSIONS: 78" x 32"
COLLAPSIBLE SIDE RAILS: 48" x 12" SECTIONAL
MATTRESS DIMENSIONS: 74" x 26" x 4"
WEIGHT CAPACITY: 700LBS / 318KG

MOVEMENT

SLIDING TOP DECK: STRETCHER TO CHAIR
SUPINE TO FULL CHAIR POSITIONING
INTEGRATED, ADJUSTABLE HEEL STIRRUPS
HEIGHT ADJUSTMENT: 24" TO 36 1/2"
HANDS-FREE TRENDLENBERG: +18° TO -18°
PNEUMATIC ASSIST FOWLER: 0° TO 90°
PNEUMATIC ASSIST LEG SECTION: 0° TO 105°
8" ALL-SURFACE, LOCKING CASTERS
ACCU-STEER 5TH WHEEL
POSI-LOCK BRAKING SYSTEM

MATTRESS

SECTIONAL, PRESSURE SENSITIVE FOAM
MATTRESS WITH VISCOLASTIC, FOUR-WAY
STRETCH COVER

CADAVER & BARIATRIC CADAVER 5000C & 5000BRC

Featuring flap access for personal belongings, and an easy-clean base with integrated storage, our cadaver stretcher is a simple, functional design for use in cadaver handling up to 600lbs (C) or 750lbs (BRC).



DIMENSIONS(5000C)

TOTAL DIMENSIONS: 78" X 29"
RECESSED LITTER: 77" X 29"
CANOPY: 76" X 28" X 14"
WEIGHT CAPACITY: 700LBS / 318KG

DIMENSIONS(5000BRC)

TOTAL DIMENSIONS: 81" X 32"
RECESSED LITTER: 80" X 32"
CANOPY: 80" X 32" X 24"
WEIGHT CAPACITY: 750LBS / 341KG

MOVEMENT

HEIGHT: 22" TO 35 1/2" (C) OR 22" TO 32" (BRC)
ACCU-STEER 5TH WHEEL
POSI-LOCK BRAKING SYSTEM
8" ALL-SURFACE, LOCKING CASTERS
OPTIONAL FOUR-SIDED ACCESS

SPECIAL FEATURES

- VINYL CANOPY COVER WITH FLAP ACCESS
- LIGHTWEIGHT, REMOVABLE CANOPY FRAME
- STAINLESS STEEL RECESSED LITTER SHEET

*See page 14 for color accent options

COLOR GUIDE

12 colors available




* Color swatches are approximate. Techlem Medical Corporation reserves the right to change specifications without notice

I.V. STAND

8012A

Stainless Brushed and optional CuVerro Shield.



OEM REPLACEMENT MATTRESSES

Transforming the standard of stretcher support systems

Techlem patient surfaces and OEM replacement mattresses are hand-crafted in North America with patients and medical professionals in mind. We use high grade materials sourced domestically to produce top quality products that are the envy of the industry.

Techlem is proud to offer a catalogue of customizable sizing and material options to suit your clinical requirements.

- + WARRANTY PROTECTED
- + INDUSTRY BEST VALUE OVER TIME
- + HIGHLY DURABLE
- + CHEMICAL-RESISTANT
- + CUSTOM + STANDARD DESIGNS
- + EASY TO CLEAN
- + TEAR-RESISTANT
- + AVAILABLE IN 3", 4", AND 5"

TECHLEM QUALITY:



LATEX-FREE

Safe for use in hospital and clinical settings.



ANTI-BACTERIAL

Foam and cover material.



DOUBLE INSULATED

Keeping the core material intact and fully functional.



FLAME RETARDENT

Meets or exceeds CTB 117.

These features
are available on
all Techlem OEM
replacement
mattresses!

CUVERROSHIELD BY AEREUS ANTIMICROBIAL COPPERCOATING


CuVerro Shield™ by Aereus Technologies is an antimicrobial permanent coating that is applied to create surfaces with 24/7 bacteria protection. It's one-time application bonds to almost any hard surface mechanically, without altering its design. This Aereus solution is cost-effective, medical grade and easy to clean.

- Kills bacterial matter virtually on contact. Remains effective after repeated contaminations.
- Inhibits the buildup and growth of bacteria between routine cleanings.
- Permanent, durable, excellent wear resistance. No need to reapply for the life of the product.



CUVERRO COPPER
CONTINUOUSLY KILLS
BACTERIA EVEN AFTER
RECONTAMINATION


Why CuVerro Shield?



ACCESSORIES	5000	5000BR	5000W	1100/N	5000XR	7600	4500-N/NS	4500GYN	5000C/BRC
Aereus Shield SideRail Upgrade (B500A)	O	O	O	O	O	N/A	N/A	N/A	N/A
KneeLift (5000KG)	O	O	O	O	O	N/A	O-N/A	N/A	N/A
Vertical O2Tank Holder (4012)	O	O	O	O	O	S	O	O	N/A
Removable,Adjustable I.V.Pole (3113)	O	O	O	O	O	S	O	O	N/A
Push-Pull Handles Set (5009)	O	O	O	O	O	N/A	O	O	O
Push-Pull Handles,Aereus Set (5009A)	O	O	O	O	O	N/A	O	O	O
Permanent I.V.Pole (5010)	O	O	O	O	O	N/A	O	O	N/A
Permanent I.V.Pole,Aereus Shield (5010A)	O	O	O	O	O	N/A	O	O	N/A
Head/Footboard/Chart Holder (3101)	O	O	O	O	O	N/A	O	O	N/A
Litter Extension(3102)	O	O	O	O	O	N/A	O	N/A	N/A
I.V.Armboard (5003)	O	O	O	O	O	N/A	O	O	N/A
Articulating Armboard (5005)	O	O	O	O	O	N/A	O	O	N/A
I.V.Tow Clamp (5004)	O	O	O	O	O	N/A	O	O	N/A
Restraint Straps (3104)	O	O	O	O	O	N/A	O	O	N/A
Safety Side Padding (5110)	O	O	O	O	O	N/A	N/A	N/A	N/A
PaperHolder (2109)	O	O	O	O	O	N/A	O	N/A	N/A
Food/Instrument Tray (5007)	O	O	O	O	O	N/A	O	O	N/A
Department I.D.Labels (5008)	S	S	S	S	S	O	S	S	N/A
Monitor Stand(5011)	O	O	O	O	O	N/A	O	O	N/A
Stirrups With Adapter (4108)	O	O	O	O	O	N/A	O	S	N/A
IntegratedUtility Tray (5019)	S	N/A	N/A	S	S	N/A	S	S	N/A
DoubleArticulating Headrest (4015)	O	O	O	O	N/A	N/A	O	O	N/A
Utility Hitch For I.V.Pole (5004)	O	O	O	O	N/A	N/A	O	O	N/A
Lateral Wrist Rest(Optionalfor 4015) (401502)	O	O	O	O	N/A	N/A	O	O	N/A
SuperiorWrist Rest (Optional For 4015) (401501)	O	O	O	O	N/A	N/A	O	O	N/A
Crank Fowler To ReplacePneumatic (5000C)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
X-RayCassette Holder - Back Section (5014)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
X-RayCassette Holder - Lateral C-Spine (5015)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
Crank Fowler To ReplacePneumatic (5000XLC)	N/A	O	O	N/A	N/A	N/A	N/A	N/A	N/A
Extra Height Side-Rails Set (5111)	O	O	O	O	O	N/A	N/A	N/A	N/A
Head/Footboard/Chart Holder- Extra Wide (3101XL)	N/A	O	O	N/A	O	N/A	N/A	N/A	N/A
Litter Extension- Extra Wide (3102XL)	N/A	O	O	N/A	O	N/A	N/A	N/A	N/A
IntegratedUtility Tray - Extra Wide (5019XL)	N/A	S	S	N/A	O	N/A	N/A	N/A	N/A
PoweredActuator - Fowler (5020)	O	O	O	O	O	N/A	O	O	N/A
PoweredActuator - KneeLift (5021)	O	O	O	O	O	N/A	O	N/A	N/A

ACCESSORIES	5000	5000BR	5000W	1100/N	5000XR	7600	4500-N/NS	4500GYN	5000C/BRC
X-RayFowler With Cassette Holder (5018)	O	N/A	N/A	O	O	N/A	O	O	N/A
X-RayFowler With Cassette Holder (5018XL)	N/A	O	O	N/A	N/A	N/A	N/A	N/A	N/A
Lift Sheet - Standard (7601C)	N/A	N/A	N/A	N/A	N/A	S	N/A	N/A	N/A
Lift Sheet - Cysto-Obstetrical (7602C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Lift Sheet - Cataract Sheet (7603C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Lift Sheet - Bath (7608C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
RodFor Cataract Sheet (7621)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Integrated Scale (5022)	O	O	O	O	O	N/A	O	O	O

*TECHLEM MEDICAL CORPORATIONRESERVES THE RIGHT TO MODIFY PARTS AND ACCESSORIESLIST, INCLUDING ADDITIONS, DISCONTINUATIONS.



Distributed in Australia by Midmed



1300 643 633
midmed.com.au
sales@midmed.com.au
Unit 4, 62 Borthwick Avenue
Murarrie QLD 4172