

## DUCANTO CATHETER™



The SSCOR DuCanto Catheter is a new catheter that is unlike traditional Yankauer suction tips that are designed for controlled surgical settings. The SSCOR DuCanto Catheter was developed with both routine and emergency airway management in mind. It has a larger internal diameter making it significantly less likely to clog and is shaped to fit the oropharyngeal anatomy. The hyper curved catheter is shaped more like a patient's airway making it easy to position the catheter.

Additionally, the SSCOR DuCanto Catheter does not have a thumb port. By removing the vacuum control, the caregiver does not need to locate and continuously occlude the thumb port to clear a patient's airway. This allows continuous suction and keeps hands free.

**SKU:** SCR.200-00002C

# INDEX

## D

DuCanto Catheter™ 1