



SUCTION EASY

Suctioning a patient is not a glamorous task, but in the presence of a blocked airway...

Undoubtedly, the ability to effectively clear the pharynx is one of the most critical of basic rescue skills. **Suction-Easy is a simple, inexpensive solution to facilitate training and use for this aspect of airway management.** It also provides you with the quickest, safest and most effective form of immediate oral suction. Suction-Easy is a must have in your first responder airway kit.

Suction-Easy has a unique double-valved design that allows continuous operation with one hand, while the other hand directs the suction tip throughout your patient's oropharynx. A consistent vacuum force is applied with each squeeze of the bulb and the material drawn into the bulb is expelled into an attached collection bag. You maintain excellent visibility throughout the procedure.

SKU: EMI610

INDEX

S

Suction Easy 1