



UP TREATMENT CHAIR

Thanks to an advanced geometric solution, the “leg” section disappears under the seat. This allows the patient to sit on the chair with the lower limbs perfectly perpendicular to the seat itself. Comfort and postural benefits will derive from this position. Therefore, all patients, even those with motor difficulties, can easily sit down and get up from the chair, even from the front, always with a correct posture. UP chair facilitates the patient’s exit with a correct postural approach. Backrest, seat and leg sections move simultaneously and automatically, raising the patients a semi-erect position according to their height. The device is equipped with three memory functions indeed, each of these positions facilitates the exit of patients according to their size. The following images represent an example of the Patient’s “Lift” or Verticalization function.

SKU: GB.UP

INDEX

U

UP Treatment Chair 1