

# TABLE OF CONTENTS

**Furniture** ..... 2

**Equipment** ..... 2

    Trolleys & Carts ..... 2



## PROCEDURE CART

### Procedure Cart for Healthcare with 6 Locking Drawers

- Large Procedure Cart for flexible storage
- New option to traditional metal carts
- Non-metal, rustproof polymer construction
- One-piece body for strength
- 6 full extension drawers (4 small, 2 large)
- Key lock secures all drawers
- Optional electronic keypad lock
- 121.92cm tall cart with 6 large drawers
- Easy to push or pull built in handle
- Durable chemical resistance
- Four 12.7cm Total Lock Casters (locks rotation and swivel)
- White polyethylene & ABS construction
- 6 full-extension drawers
- Four drawers are 7.62 inches deep
- Two drawers are 17.78 inches deep
- Built-in 3 inch deep tray top

**SKU:** TR.XXXX

# INDEX

## P

Procedure Cart 1