

TABLE OF CONTENTS

- Techlem Stretchers** 2
- Consumables** 2
 - Airway Management 2
 - Curtains & Wall Displays 4
 - Enteral Feeding 5
 - Enteral Feeding Syringes 5
 - Fluid Warming 5
 - Hygiene 6
 - Patient Warming 6
 - Securement Devices 7
 - Suction 10
- Furniture** 11
 - Stretchers 11
- Pre-Hospital** 12
 - Airway Management 12
 - Bags & Kits 13
 - Patient Handling 14
 - Patient Warming 15
 - Triage 15
- Equipment** 16
 - Patient Transfer 16
 - Lighting 17
 - Trolleys & Carts 20
- Cair & Micsafe Specials** 21
- Medical Furniture & Equipment Specials** 21



TRANSPORT/GENERAL STRETCHER

The premier Techlem stretcher for Transport and General Use features an integrated oxygen tank holder and patient storage, 5th wheel technology for steering assist, and more than 20 customisable options for design and functionality. The Premier Series 5000 is the preferred stretcher in emergency and trauma centres worldwide, due to its durability and ease of use. It also features a high density foam mattress with conductive cover.

SKU: T.5000SAHF



SUCTION EASY

Suctioning a patient is not a glamorous task, but in the presence of a blocked airway...

Undoubtedly, the ability to effectively clear the pharynx is one of the most critical of basic rescue skills. **Suction-Easy is a simple, inexpensive solution to facilitate training and use for this aspect of airway management.** It also provides you with the quickest, safest and most effective form of immediate oral suction. Suction-Easy is a must have in your first responder airway kit.

Suction-Easy has a unique double-valved design that allows continuous operation with one hand, while the other hand directs the suction tip throughout your patient's oropharynx. A consistent vacuum force is applied with each squeeze of the bulb and the material drawn into the bulb is expelled into an attached collection bag. You maintain excellent visibility throughout the procedure.

SKU: EMI610




QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230



PYLANT MONITOR - ENDOTRACHEAL TUBE CUFF PRESSURE MANOMETER

SKU: INO.IM.100



TRACHEOSTOMY TUBE HOLDER

SKU: DA.240



ENDOTRACHEAL TUBE HOLDER

SKU: DA.273



DISPOSABLE CURTAINS

Each EcoMed disposable curtain is treated with the 3rd Generation Ultra-Fresh built-in antimicrobial protection.

SKU: EN.



SINGLE USE ENFIT SYRINGES


Syringes from 1ml to 50/60 ml for enteral medication, flushing, gastric residuals and breast milk or enteral feeding administration by gravity, bolus or syringe-pump.

SKU: N/A



QINFLOW WARRIOR - BLOOD/FLUID WARMER

SKU: QF.Q211050000SUC



WRASFER™ TRANSFER SHEET

The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

SKU: WRASFER



QINFLOW WARRIOR - BLOOD/FLUID WARMER


SKU: QF.Q211050000SUC



NASOGASTRIC TUBE HOLDER

The NasoGastric Tube Holder is ideal for holding nasogastric feeding and aspiration tubes, such as Salem Sump, Levin Tubes, and Feeding Catheter Tubes (Dobbhoff). The unique dual interlocking tabs spiral around the nasogastric tube ensuring superior tube security. The non-adherent plastic tips on both tabs make tube removal safer and easier without scissors.

SKU: DA.161




CATHETER HOLDER ADHESIVE PATCH

The Hold-n-Place Adhesive Patch is a skin friendly, breathable adhesive device, versatile enough for all of your indwelling Foley securement needs. The adhesive patch, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.

All Hold-n-Place Foley Catheter Holders feature a locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma. The locking devices speed up application and simplify adjustment. The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters.

SKU: DA.150



CATHETER HOLDER LEG BAND

The Leg Band, made entirely of stretch material, has a soft, supple backing for added skin comfort. The elasticity and 2" width of the leg band distribute compression evenly around the high to minimize circumferential pressure.

SKU: DA.316



BENDABLE IV ARMBOARD

The Bendable ArmBoard is a simple, cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

It can be used in the Operating Room, Recovery Room, ICU, Patient Care Units, and more.

SKU: DA.650



IV ARMOR

Midmed's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering. It has several benefits including reducing the introduction of infection associated with IV reinsertion, and reducing downstream occlusion alarms. The IV-ARMOR also reduces risk of catheter kinking in a bend such as a wrist or elbow, and reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm. IV-ARMOR uses skin-friendly adhesive, and comes in a one size fits most box of 10.

SKU: DA.675



TRACHEOSTOMY TUBE HOLDER

SKU: DA.240



ENDOTRACHEAL TUBE HOLDER

SKU: DA.273



SUCTION EASY

Suctioning a patient is not a glamorous task, but in the presence of a blocked airway...

Undoubtedly, the ability to effectively clear the pharynx is one of the most critical of basic rescue skills. **Suction-Easy is a simple, inexpensive solution to facilitate training and use for this aspect of airway management.** It also provides you with the quickest, safest and most effective form of immediate oral suction. Suction-Easy is a must have in your first responder airway kit.

Suction-Easy has a unique double-valved design that allows continuous operation with one hand, while the other hand directs the suction tip throughout your patient's oropharynx. A consistent vacuum force is applied with each squeeze of the bulb and the material drawn into the bulb is expelled into an attached collection bag. You maintain excellent visibility throughout the procedure.

SKU: EMI610



QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230



BARIATRIC STRETCHER

The Premier Series 5000BR is a safe, reliable and nimble stretcher equipped with Techlem's 5th wheel technology and high tensile components. The Bariatric stretcher boasts a 750 pound weight capacity. It also features a high density foam mattress with conductive cover.

SKU: T.5000BRSAHF



TRANSPORT/GENERAL STRETCHER

The premier Techlem stretcher for Transport and General Use features an integrated oxygen tank holder and patient storage, 5th wheel technology for steering assist, and more than 20 customisable options for design and functionality. The Premier Series 5000 is the preferred stretcher in emergency and trauma centres worldwide, due to its durability and ease of use. It also features a high density foam mattress with conductive cover.

SKU: T.5000SAHF



QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230



STOR-IT BAG

The SMART Stor-It Reflective Triage Bag is an extremely well organised solution for bulk storage and rapid deployment of medical equipment. Made to the highest specification the Stor-It is designed specifically for disaster and emergency response but has been shown to be highly effective at aid stations during mass gathering events.

The rugged storage system has triple reinforced stitching and exterior padding to protect fragile medical equipment during storage, transport and shipment. It is easy to stack multiple units and the inside divides into two easy to hang sections with highly visible contents.

NOTE: Bag Only

SKU: TS.STOR-IT.

SMARTSAFE VALUABLE STORAGE BAGS



SMARTSAFE small property bags are an effective way of increasing patient comfort and security when they need it most. These bags can be used to store many items, including valuables such as jewellery or medication.

An evident tamper seal and barcode are used to track the bag's location, and for the patient's peace of mind, the barcode is attachable to patient notes. The bags also include a wristband that is to be placed around the patient's wrist. Clearly labelled forms further help identify the belongings of each patient. The bags also include a wristband that is to be placed around the patient's wrist.


SKU: N/A

BANANA SLIDE SINGLE PATIENT TUBULAR



The GBUK Banana Single Patient Tubular Slides are strong and durable with a high slip performance. It has a strong, red, hemmed slide sheet and a silicone coat on both sides. With a low profile seam, they are long-lasting and have high-slip characteristics, enabling easy and safe patient transfers and re-positioning. Available in two sizes: 72 cm by 72 cm and 100 cm by 200 cm.

SKU: N/A



WRASFER™ TRANSFER SHEET

The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

SKU: WRASFER



QINFLOW WARRIOR - BLOOD/FLUID WARMER

SKU: QF.Q211050000SUC




STOR-IT BAG

The SMART Stor-It Reflective Triage Bag is an extremely well organised solution for bulk storage and rapid deployment of medical equipment. Made to the highest specification the Stor-It is designed specifically for disaster and emergency response but has been shown to be highly effective at aide stations during mass gathering events.

The rugged storage system has triple reinforced stitching and exterior padding to protect fragile medical equipment during storage, transport and shipment. It is easy to stack multiple units and the inside divides into two easy to hang sections with highly visible contents.

NOTE: Bag Only

SKU: TS.STOR-IT.



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48

PATIENT TRANSFER SCALE



Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS



SL430 LED MINOR SURGICAL LIGHT

The Daray SL430 LED minor surgical light is specifically designed for minor procedures under local anaesthetic in health centres, GP surgeries and outpatient departments. It is ideal for minor surgery, maxillofacial, obstetrics, diagnostics and gynaecology as well as for use within casualty departments and intensive care units. This light is consistent, accurate and cost effective. With great shadow reduction and accurate colour rendering, it is an invaluable tool.

SKU: D.SL430LC



X340 LED EXAMINATION LIGHT

The Daray X340 examination light is specifically designed for use during medical consultation and for diagnosis. Featuring a wide range of mount options, it ensures excellent placement options. The smooth spring balanced arm and telescopic ceiling tube makes this light extremely versatile. It has an intense beam and excellent colour rendering properties. This light has a low cost of ownership and around 50,000 hours of life, equating to approximately 25 years of clinical use. Requires minimal maintenance.

SKU: D.X340



X200 LED EXAMINATION LIGHT

The Daray X200 examination light range is perfectly suited for use during medical consultations and for diagnosis. With LED technology, this light range provides an intense beam with excellent colour rendering properties. The range has a low total ownership cost and around 50, 000 hours of life, equating to approximately 25 years of clinical use. It features a fully enclosed BioProtect head to assist with the management of infection control. Minimal maintenance is required.

SKU: D.X200



X700 LED EXAMINATION LIGHT

The Daray X700 light range are premium quality, powerful and incredibly flexible examination lights, featuring variable colour adjustment. The range is available with a wide range of mount options to ensure your light is exactly where you need it. Additionally, there are two arm options – flexible gooseneck or spring-balanced arm – and a large range of hardwired and plug-top mount options to suit all clinical needs. They have a low cost ownership and around 50,000 hours of life, equating to approximately 25 years of clinical use. These lights are ideal for GP surgeries and clinics where flexibility and performance are crucial.

SKU: D.X740



X400 LED EXAMINATION LIGHT

The Daray X400 examination light is a premium quality, extremely efficient general examination light. With a large diameter light-patch, this light provides excellent illumination. Single-handed lighting control, on/off and variable intensity, can all be controlled from a single switch located next to the light head. A fully enclosed light head and arm system assists with the management of infection control. Available as ceiling, wall and mobile models, this product is a perfect complement to any examination room.

SKU: D.X400



SL180 LED MINOR SURGICAL LIGHT

The Daray SL180 LED minor surgical light features a strong beam of light that has an even intensity across the viewing area. With various mounting options, this light is ideal for gynaecology, A&E departments, day surgery, GP treatment rooms, podiatry, dermatology and other clinical practices. It has variable light intensity and around 50,000 hours of life, equating to approximately 25 years of clinical use. BioProtect anti-microbial solution coats all painted finishes, assisting in infection control. This light is a great complement to any treatment room.

SKU: D.SL180



SL720 LED MINOR SURGICAL LIGHT

The Daray SL720 LED minor surgical light is the newest addition to the range of minor procedure lighting. It is suitable for multiple minor surgical environments including health centres, GP surgeries and maxillofacial departments, as well as within casualty or intensive care. Designed for use in minor procedures, the SL720 provides clean, accurate light where it's needed. Featuring variable colour temperature, it is incredibly versatile and easily adapted to suit your specific needs. The wall and ceiling editions of this light come on spring balanced arms with a large range of movement as standard. This, combined with shadow reduction technology and low running costs, make this light a smart choice for your clinic.

SKU: D.SL720



AURION CARTS

Aurion Medical carts are innovative, high quality and built to last. These medical trolleys have a tough design, numerous compartments and are safe for use in any medical field. The Aurion range includes carts for use in emergency, anaesthetics, therapy, paediatrics, operating rooms and more.

SKU: AU.CP.



SHORT-TERM ENFIT FEEDING TUBES

Our unique PVC enteral tubes provide compatibility across the brand's range of enteral feeding products. They meet ISO 80369-3 standards and come with a number of distinct features.

- Readymade male EnFit™ connectors
- Numerical markings (in centimetres)
- Closed rounded tip
- 2 lateral distal openings
- 2 radiopaque stripes over transparent tube
- Available in a range of lengths and sizes from 4 Ch/Fr to 10 Ch/Fr and in a variety of widths

SKU: CA.NCE110P



PATIENT TRANSFER SCALE

Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS

INDEX

A

Aurion Carts 20

B

Banana Slide Single Patient Tubular 13

Bariatric Stretcher 11

Bendable IV ArmBoard 8

C

Catheter Holder Adhesive Patch 7

Catheter Holder Leg Band 7

D

Disposable Curtains 4

E

Endotracheal Tube Holder 3, 9

I

IV Armor 8

N

Nasogastric Tube Holder 6

P

Patient Transfer Scale 16, 21

Pylant Monitor - Endotracheal Tube Cuff Pressure Manometer 2

Q

QinFlow Warrior - Blood/Fluid Warmer 5, 6, 14

Quickdraw® Portable Suction Unit 2, 10, 12

S

Short-Term ENFIT Feeding Tubes 20

Single Use ENFIT Syringes 4

SL180 LED Minor Surgical Light 19

SL430 LED Minor Surgical Light 16

SL720 LED Minor Surgical Light 19

SmartSafe Valuable Storage Bags 13

Standard Evacuation Sled 15

Stor-It Bag 12, 15

Suction Easy 1, 10

T

Tracheostomy Tube Holder 3, 9

Transport/General Stretcher 1, 11

W

WRASFER™ Transfer Sheet 5, 14

X

X200 LED Examination Light 17

X340 LED Examination Light 17

X400 LED Examination Light 18

X700 LED Examination Light 18