

TABLE OF CONTENTS

- Consumables** 2
 - Airway Management 2
 - Curtains & Wall Displays 4
 - Enteral Feeding 4
 - Fluid Warming 6
 - Hygiene 7
 - Patient Warming 7
 - Securement Devices 8
 - Suction 11
- Furniture** 12
 - Bariatric Furniture 12
 - Beds 13
 - Mattresses 14
 - Stretchers 15
- Pre-Hospital** 16
 - Airway Management 16
 - Bags & Kits 17
 - Patient Handling 17
 - Patient Transfer 20
 - Patient Warming 21
- Equipment** 22
 - Lighting 22
 - Trolleys & Carts 27
- Specials** 28
 - Cair & Micsafe Specials** 29
 - Medical Furniture & Equipment Specials** 29



SUCTION EASY

Suctioning a patient is not a glamorous task, but in the presence of a blocked airway...

Undoubtedly, the ability to effectively clear the pharynx is one of the most critical of basic rescue skills. **Suction-Easy is a simple, inexpensive solution to facilitate training and use for this aspect of airway management.** It also provides you with the quickest, safest and most effective form of immediate oral suction. Suction-Easy is a must have in your first responder airway kit.

Suction-Easy has a unique double-valved design that allows continuous operation with one hand, while the other hand directs the suction tip throughout your patient's oropharynx. A consistent vacuum force is applied with each squeeze of the bulb and the material drawn into the bulb is expelled into an attached collection bag. You maintain excellent visibility throughout the procedure.

SKU: EMI610




QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230



MANOMETER

SKU: INO.IM.100



TRACHEOSTOMY TUBE HOLDER

SKU: DA.240



ENDOTRACHEAL TUBE HOLDER

SKU: DA.273



DISPOSABLE CURTAINS

Each EcoMed disposable curtain is treated with the 3rd Generation Ultra-Fresh built-in antimicrobial protection.

SKU: EN.



NASOGASTRIC TUBE HOLDER

The NasoGastric Tube Holder is ideal for holding nasogastric feeding and aspiration tubes, such as Salem Sump, Levin Tubes, and Feeding Catheter Tubes (Dobbhoff). The unique dual interlocking tabs spiral around the nasogastric tube ensuring superior tube security. The non-adherent plastic tips on both tabs make tube removal safer and easier without scissors.

SKU: DA.161



PH TESTING STRIPS

This pH Test strip pack is CE marked and meets the special demands of health care professionals to ensure rapid and safe results for medical pH-testing. These Test strips range from 2.0 to 9.0 with 0.5 increments that make for accurate testing.

SKU: MN.92118



SINGLE USE ENFIT SYRINGES

Syringes from 1ml to 50/60 ml for enteral medication, flushing, gastric residuals and breast milk or enteral feeding administration by gravity, bolus or syringe-pump.

SKU: MS.SE60



SHORT-TERM ENFIT FEEDING TUBES

Our unique PVC enteral tubes provide compatibility across the brand's range of enteral feeding products. They meet ISO 80369-3 standards and come with a number of distinct features.


- Readymade male EnFit™ connectors
- Numerical markings (in centimetres)
- Closed rounded tip
- 2 lateral distal openings
- 2 radiopaque stripes over transparent tube
- Available in a range of lengths and sizes from 4 Ch/Fr to 10 Ch/Fr and in a variety of widths

SKU: CA.NCE110P



QINFLOW WARRIOR - BLOOD/FLUID WARMER

SKU: QF.Q211050000SUC



WRASFER™ TRANSFER SHEET

The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

SKU: WRASFER



QINFLOW WARRIOR - BLOOD/FLUID WARMER


SKU: QF.Q211050000SUC



NASOGASTRIC TUBE HOLDER

The NasoGastric Tube Holder is ideal for holding nasogastric feeding and aspiration tubes, such as Salem Sump, Levin Tubes, and Feeding Catheter Tubes (Dobbhoff). The unique dual interlocking tabs spiral around the nasogastric tube ensuring superior tube security. The non-adherent plastic tips on both tabs make tube removal safer and easier without scissors.

SKU: DA.161



CATHETER HOLDER ADHESIVE PATCH

The Hold-n-Place Adhesive Patch is a skin friendly, breathable adhesive device, versatile enough for all of your indwelling Foley securement needs. The adhesive patch, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.

All Hold-n-Place Foley Catheter Holders feature a locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma. The locking devices speed up application and simplify adjustment. The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters.

SKU: DA.150



CATHETER HOLDER LEG BAND

The Leg Band, made entirely of stretch material, has a soft, supple backing for added skin comfort. The elasticity and 2" width of the leg band distribute compression evenly around the high to minimize circumferential pressure.

SKU: DA.316



BENDABLE IV ARMBOARD

The Bendable ArmBoard is a simple, cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

It can be used in the Operating Room, Recovery Room, ICU, Patient Care Units, and more.

SKU: DA.650



IV ARMOR

Midmed's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering. It has several benefits including reducing the introduction of infection associated with IV reinsertion, and reducing downstream occlusion alarms. The IV-ARMOR also reduces risk of catheter kinking in a bend such as a wrist or elbow, and reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm. IV-ARMOR uses skin-friendly adhesive, and comes in a one size fits most box of 10.

SKU: DA.675



TRACHEOSTOMY TUBE HOLDER

SKU: DA.240



ENDOTRACHEAL TUBE HOLDER

SKU: DA.273



SUCTION EASY

Suctioning a patient is not a glamorous task, but in the presence of a blocked airway...

Undoubtedly, the ability to effectively clear the pharynx is one of the most critical of basic rescue skills. **Suction-Easy is a simple, inexpensive solution to facilitate training and use for this aspect of airway management.** It also provides you with the quickest, safest and most effective form of immediate oral suction. Suction-Easy is a must have in your first responder airway kit.

Suction-Easy has a unique double-valved design that allows continuous operation with one hand, while the other hand directs the suction tip throughout your patient's oropharynx. A consistent vacuum force is applied with each squeeze of the bulb and the material drawn into the bulb is expelled into an attached collection bag. You maintain excellent visibility throughout the procedure.

SKU: EMI610



QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230




BARIATRIC SAPPHIRE AIR MATTRESS

SKU: SW.33060420



BARIATRIC BED TITAN

SKU: PA.300S-020-00-105



SOLO ULTRA LOW HOSPITAL BED

SKU: SOLO.SPECIAL



BARIATRIC BED TITAN

SKU: PA.300S-020-00-105



BARIATRIC SAPPHIRE AIR MATTRESS

SKU: SW.33060420



COMPACT STRETCHER

Techlem's Premier Series 4500CFS is a highly adaptable, compact stretcher designed to suit your small-space, clinical environment. It provides many features as the 4500, including collapsing side-rails. It's available in narrow and short & narrow, and features a high density foam mattress with conductive cover.

SKU: T.4500CFS



BARIATRIC STRETCHER

The Premier Series 5000BR is a safe, reliable and nimble stretcher equipped with Techlem's 5th wheel technology and high tensile components. The Bariatric stretcher boasts a 750 pound weight capacity. It also features a high density foam mattress with conductive cover.

SKU: T.5000BRSAHF



TRANSPORT/GENERAL STRETCHER

The premier Techlem stretcher for Transport and General Use features an integrated oxygen tank holder and patient storage, 5th wheel technology for steering assist, and more than 20 customisable options for design and functionality. The Premier Series 5000 is the preferred stretcher in emergency and trauma centres worldwide, due to its durability and ease of use. It also features a high density foam mattress with conductive cover.

SKU: T.5000SAHF



QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

Available in both rechargeable and battery operated variations

SKU: SCR.2403-230



SMARTSAFE VALUABLE STORAGE BAGS

SKU: TS.7325(0/1)



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48



EXTRA ERGOLIFT STAIR CHAIR

Extra Ergolift is the powered version of the Extra track assisted stair chair with the same main structural features. The real innovation of Extra Ergolift is given by the electric motor that moves the tracks up and down the stairs without any effort by the operator. The motor is activated by a wireless console. The battery ensures a high autonomy and, in case of sudden discharging, this device will allow you to safely complete the service.

Extra Ergolift loads up to 200 kg in powered transport configuration and 230 kg in manual transportation. It can be used under rain or adverse weather conditions.

SKU: MR.671



STAIR CHAIR

SKU: B.STAIRCHAIR.NEW



BANANA SLIDES

The GBUK Banana Single Patient Tubular Slides are strong and durable with a high slip performance. It has a strong, red, hemmed slide sheet and a silicone coat on both sides. With a low profile seam, they are long-lasting and have high-slip characteristics, enabling easy and safe patient transfers and re-positioning. They are available in two sizes: 72 cm by 72 cm and 100 cm by 200 cm.

SKU: EGT



WRASFER™ TRANSFER SHEET

The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

SKU: WRASFER



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48

PATIENT TRANSFER SCALE



Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS



EXTRA ERGOLIFT STAIR CHAIR

Extra Ergolift is the powered version of the Extra track assisted stair chair with the same main structural features. The real innovation of Extra Ergolift is given by the electric motor that moves the tracks up and down the stairs without any effort by the operator. The motor is activated by a wireless console. The battery ensures a high autonomy and, in case of sudden discharging, this device will allow you to safely complete the service.

Extra Ergolift loads up to 200 kg in powered transport configuration and 230 kg in manual transportation. It can be used under rain or adverse weather conditions.

SKU: MR.671



QINFLOW WARRIOR - BLOOD/FLUID WARMER

SKU: QF.Q211050000SUC



LEO MINOR LED LIGHT

Design innovation meets the latest advancements in LED lighting technology in the Leo Minor® multipurpose LED surgical light. The Leo Minor features a unique patented optical solutions that produces a uniform and homogeneous light pattern with excellent shadow dilution and control. When you combine these features with daylight-like colour rendering and colour temperature, the end result is a true surgical light performance meeting or exceeding IEC requirements.

SKU: SU.LMI-



SL430 LED MINOR SURGICAL LIGHT

The Daray SL430 LED minor surgical light is specifically designed for minor procedures under local anaesthetic in health centres, GP surgeries and outpatient departments. It is ideal for minor surgery, maxillofacial, obstetrics, diagnostics and gynaecology as well as for use within casualty departments and intensive care units. This light is consistent, accurate and cost effective. With great shadow reduction reduction and accurate colour rendering, it is an invaluable tool.

SKU: D.SL430LC



TRI-STAR LED

The Sunnex Tri-Star LED light is excellent for minor surgical and examination applications. The cool, bright light can be easily directed where needed and the caster base allows for easy portability. This light is ideal for medical procedures in environments where space is limited. Additionally, with the telescoping pole, height adjustment in the mobile configuration is possible. Wall configuration offers extensive reach with the innovative arm design. Ideal for outpatient surgery centers, examination and emergency rooms and veterinarian surgical suites.

SKU: SU.TSLED

CELESTIAL STAR LED

SKU: SU.CS2050



X340 LED EXAMINATION LIGHT

The Daray X340 examination light is specifically designed for use during medical consultation and for diagnosis. Featuring a wide range of mount options, it ensures excellent placement options. The smooth spring balanced arm and telescopic ceiling tube makes this light extremely versatile. It has an intense beam and excellent colour rendering properties. This light has a low cost of ownership and around 50,000 hours of life, equating to approximately 25 years of clinical use. Requires minimal maintenance.

SKU: D.X340





X200 LED EXAMINATION LIGHT

The Daray X200 examination light range is perfectly suited for use during medical consultations and for diagnosis. With LED technology, this light range provides an intense beam with excellent colour rendering properties. The range has a low total ownership cost and around 50,000 hours of life, equating to approximately 25 years of clinical use. It features a fully enclosed BioProtect head to assist with the management of infection control. Minimal maintenance is required.

SKU: D.X200



X700 LED EXAMINATION LIGHT

The Daray X700 light range are premium quality, powerful and incredibly flexible examination lights, featuring variable colour adjustment. The range is available with a wide range of mount options to ensure your light is exactly where you need it. Additionally, there are two arm options – flexible gooseneck or spring-balanced arm – and a large range of hardwired and plug-top mount options to suit all clinical needs. They have a low cost ownership and around 50,000 hours of life, equating to approximately 25 years of clinical use. These lights are ideal for GP surgeries and clinics where flexibility and performance are crucial.

SKU: D.X740



X400 LED EXAMINATION LIGHT

The Daray X400 examination light is a premium quality, extremely efficient general examination light. With a large diameter light-patch, this light provides excellent illumination. Single-handed lighting control, on/off and variable intensity, can all be controlled from a single switch located next to the light head. A fully enclosed light head and arm system assists with the management of infection control. Available as ceiling, wall and mobile models, this product is a perfect complement to any examination room.

SKU: D.X400



SL180 LED MINOR SURGICAL LIGHT

The Daray SL180 LED minor surgical light features a strong beam of light that has an even intensity across the viewing area. With various mounting options, this light is ideal for gynaecology, A&E departments, day surgery, GP treatment rooms, podiatry, dermatology and other clinical practices. It has variable light intensity and around 50,000 hours of life, equating to approximately 25 years of clinical use. BioProtect anti-microbial solution coats all painted finishes, assisting in infection control. This light is a great complement to any treatment room.

SKU: D.SL180



SL720 LED MINOR SURGICAL LIGHT

The Daray SL720 LED minor surgical light is the newest addition to the range of minor procedure lighting. It is suitable for multiple minor surgical environments including health centres, GP surgeries and maxillofacial departments, as well as within casualty or intensive care. Designed for use in minor procedures, the SL720 provides clean, accurate light where it's needed. Featuring variable colour temperature, it is incredibly versatile and easily adapted to suit your specific needs. The wall and ceiling editions of this light come on spring balanced arms with a large range of movement as standard. This, combined with shadow reduction technology and low running costs, make this light a smart choice for your clinic.

SKU: D.SL720



PROCEDURE CART

Procedure Cart for Healthcare with 6 Locking Drawers

- Large Procedure Cart for flexible storage
- New option to traditional metal carts
- Non-metal, rustproof polymer construction
- One-piece body for strength
- 6 full extension drawers (4 small, 2 large)
- Key lock secures all drawers
- Optional electronic keypad lock
- 121.92cm tall cart with 6 large drawers
- Easy to push or pull built in handle
- Durable chemical resistance
- Four 12.7cm Total Lock Casters (locks rotation and swivel)
- White polyethylene & ABS construction
- 6 full-extension drawers
- Four drawers are 7.62 inches deep
- Two drawers are 17.78 inches deep
- Built-in 3 inch deep tray top

SKU: TR.XXXX



AURION CARTS

Aurion Medical carts are innovative, high quality and built to last. These medical trolleys have a tough design, numerous compartments and are safe for use in any medical field. The Aurion range includes carts for use in emergency, anaesthetics, therapy, paediatrics, operating rooms and more.

SKU: AU.CP.



BARIATRIC SAPPHIRE AIR MATTRESS

SKU: SW.33060420

BARIATRIC BED TITAN

SKU: PA.300S-020-00-105



SHORT-TERM ENFIT FEEDING TUBES

Our unique PVC enteral tubes provide compatibility across the brand's range of enteral feeding products. They meet ISO 80369-3 standards and come with a number of distinct features.

- Readymade male EnFit™ connectors
- Numerical markings (in centimetres)
- Closed rounded tip
- 2 lateral distal openings
- 2 radiopaque stripes over transparent tube
- Available in a range of lengths and sizes from 4 Ch/Fr to 10 Ch/Fr and in a variety of widths

SKU: CA.NCE110P





PATIENT TRANSFER SCALE

Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS



BARIATRIC SAPPHIRE AIR MATTRESS

SKU: SW.33060420



BARIATRIC BED TITAN

SKU: PA.300S-020-00-105

INDEX

A

Aurion Carts 27

B

Banana Slides 18

Bariatric Bed TITAN 12, 13, 28, 30

Bariatric Sapphire Air Mattress 12, 14, 27, 29

Bariatric Stretcher 15

Bendable IV ArmBoard 9

C

Catheter Holder Adhesive Patch 8

Catheter Holder Leg Band 8

Celestial Star LED 23

Compact Stretcher 14

D

Disposable Curtains 3

E

Endotracheal Tube Holder 3, 10

EXTRA ERGOLIFT Stair Chair 17, 20

I

IV Armor 9

L

Leo Minor LED Light 21

M

Manometer 2

N

Nasogastric Tube Holder 4, 7

P

Patient Transfer Scale 20, 29

pH Testing Strips 4

Procedure Cart 26

Q

QinFlow Warrior - Blood/Fluid Warmer 6, 7, 21

Quickdraw® Portable Suction Unit 1, 11, 16

S

Short-Term ENFIT Feeding Tubes 5, 28

Single Use ENFIT Syringes 5

SL180 LED Minor Surgical Light 25

SL430 LED Minor Surgical Light 22

SL720 LED Minor Surgical Light 26

SmartSafe Valuable Storage Bags 16

Solo Ultra Low Hospital Bed 13

Stair Chair 18

Standard Evacuation Sled 17, 19

Suction Easy 1, 11

T

Tracheostomy Tube Holder 2, 10

Transport/General Stretcher 15

Tri-Star LED 22

W

WRASFER™ Transfer Sheet 6, 19

X

X200 LED Examination Light 24

X340 LED Examination Light 23

X400 LED Examination Light 25

X700 LED Examination Light 24