

TABLE OF CONTENTS

- Consumables** 2
 - Curtains & Wall Displays 2
 - Enteral Feeding 2
 - Maternity 4
 - Securement Devices 6
- Furniture** 9
 - Bariatric Furniture 9
 - Beds 10
 - Examination Tables 12
 - Recliners 12
- Pre-Hospital** 14
 - Patient Handling 14
 - Patient Transfer 15
- Equipment** 17
 - Breast Pumps & Milk Warmers 17
 - Trolleys & Carts 20
- Cair & Micsafe Specials** 21
- Medical Furniture & Equipment Specials** 21



DISPOSABLE CURTAINS

Each EcoMed disposable curtain is treated with the 3rd Generation Ultra-Fresh built-in antimicrobial protection.

SKU: EN.



NASOGASTRIC TUBE HOLDER

The NasoGastric Tube Holder is ideal for holding nasogastric feeding and aspiration tubes, such as Salem Sump, Levin Tubes, and Feeding Catheter Tubes (Dobbhoff). The unique dual interlocking tabs spiral around the nasogastric tube ensuring superior tube security. The non-adherent plastic tips on both tabs make tube removal safer and easier without scissors.

SKU: DA.161



PH TESTING STRIPS

This pH Test strip pack is CE marked and meets the special demands of health care professionals to ensure rapid and safe results for medical pH-testing. These Test strips range from 2.0 to 9.0 with 0.5 increments that make for accurate testing.

SKU: MN.92118



SINGLE USE ENFIT SYRINGES

Syringes from 1ml to 50/60 ml for enteral medication, flushing, gastric residuals and breast milk or enteral feeding administration by gravity, bolus or syringe-pump.

SKU: MS.SE60

SHORT-TERM ENFIT FEEDING TUBES

Our unique PVC enteral tubes provide compatibility across the brand's range of enteral feeding products. They meet ISO 80369-3 standards and come with a number of distinct features.

- Readymade male EnFit™ connectors
- Numerical markings (in centimetres)
- Closed rounded tip
- 2 lateral distal openings
- 2 radiopaque stripes over transparent tube
- Available in a range of lengths and sizes from 4 Ch/Fr to 10 Ch/Fr and in a variety of widths

SKU: CA.NCE110P



PENGUIN NUTRITIONAL WARMER

Going above and beyond just warming breast milk, The Ameda Penguin Nutritional warmer ensure that each feeding is always safely delivered within the ideal human physiological temperature range of 36-37°C, regardless of the starting temperature. A feature which is unique to the Ameda Penguin Warmer is our gentle vibration technology which prevents potential nutritional degradation by continuously stirring the milk throughout the warming cycle for a more even and consistent warming. Our technology then eliminates the need to shake or roll the milk to even out any dangerous hot spots.

SKU: AM.PNW00115AU



FLOWCUP

The Materni flowcup is a simple, convenient accessory for feeding infants who have difficulty sucking.

SKU: MI.MEP



FLOWSPOON

A simple and convenient accessory for safely giving the first breast milk (colostrum).

SKU: MI.MEPAO112



NON-STERILE BABY BOTTLES

Materni bottles are suitable for the hygienic collection and storage of breast milk and the preparation of formula. They are bottle warmer, microwave, dishwasher and freezer safe making them highly versatile. Bottles have a screw on cap to ensure micro-biological cleanliness and have clear volume markings to make milk preparation quick and easy. They are available in four sizes: 50 ml, 80 ml, 120 ml and 240 ml.

SKU: MI.MEPA0109



STERILE BABY BOTTLES

Our disposable baby bottles come in a variety of configurations.


SKU: CA.NCB1240A



NASOGASTRIC TUBE HOLDER

The NasoGastric Tube Holder is ideal for holding nasogastric feeding and aspiration tubes, such as Salem Sump, Levin Tubes, and Feeding Catheter Tubes (Dobbhoff). The unique dual interlocking tabs spiral around the nasogastric tube ensuring superior tube security. The non-adherent plastic tips on both tabs make tube removal safer and easier without scissors.

SKU: DA.161




CATHETER HOLDER ADHESIVE PATCH

The Hold-n-Place Adhesive Patch is a skin friendly, breathable adhesive device, versatile enough for all of your indwelling Foley securement needs. The adhesive patch, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.

All Hold-n-Place Foley Catheter Holders feature a locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma. The locking devices speed up application and simplify adjustment. The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters.

SKU: DA.150



CATHETER HOLDER LEG BAND

The Leg Band, made entirely of stretch material, has a soft, supple backing for added skin comfort. The elasticity and 2" width of the leg band distribute compression evenly around the high to minimize circumferential pressure.

SKU: DA.316



BENDABLE IV ARMBOARD

The Bendable ArmBoard is a simple, cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

It can be used in the Operating Room, Recovery Room, ICU, Patient Care Units, and more.

SKU: DA.650



IV ARMOR

Midmed's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering. It has several benefits including reducing the introduction of infection associated with IV reinsertion, and reducing downstream occlusion alarms. The IV-ARMOR also reduces risk of catheter kinking in a bend such as a wrist or elbow, and reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm. IV-ARMOR uses skin-friendly adhesive, and comes in a one size fits most box of 10.

SKU: DA.675



56 SERIES BARIATRIC RECLINER

SKU: CH.56/CH.56P/CH.566-2



65 SERIES BARIATRIC RECLINER

The Ascent XL Bariatric Recliner is ruggedly built like our classic recliners, with added style and elegance. Gently contoured lines enhance its profile while comforting patients up to 500 lbs. Standard swing-arms provide easy cleaning, maintenance and patient transfer. Four-inch, twin wheel casters are corrosion resistant and roll with ease.

Suitable for Oncology, Surgery, Infusion, Dialysis Therapies, Orthopedics, Paediatrics, and more.

Purchase the Ascent XL patient chair through Midmed today for an easily transportable, well-designed recliner that always puts comfort first. Available with manual and power options in a wide range of colours.

SKU: CH.65/CH.65P/CH.660-2



OVERNIGHTER RECLINER

Founded on a strong steel frame, the 526 Overnighter Series patient chair is a hard-working sleeper that provides the greatest possible level of comfort to overnight patients. Mobile and impressively flexible given its size, the 526 has high-quality casters that make transport easy and non-intrusive.

It's open space, all-in-one design helps reduce claustrophobic sensations often caused by constant sitting. The quality foam and fabrics of the 526 Series make for real comfort, as well as a better-looking and longer-lasting product. Maximize patients' comfort by helping them choose between three positions. Years of research and testing have gone into the perfection of the Overnighter Series. It's specially designed to ensure that each has every opportunity to get a good night's sleep, whatever their medical situation.

Available with manual and power options in a wide range of colours.

SKU: CH.526-2



NEWBORN BABY BED

SKU: PA.300S-302-03-035



DANIEL PAEDIATRIC BED

The Gardhen Balance Daniel Paediatric Bed has an accurate, integrated weighing system. Weight display automatically displays weight. It has a anti-blackout function with rechargeable long life batteries, an automatic and manual tare mode and easy storage and transfer of data. The bed has a quick bilateral release in the back section for CPR, an ergonomic patient surface and a bed section protected with HPL panels (easy to be cleaned).

SKU: GB.DANIEL



CHILDREN'S BED - TEDDY

SKU: PA.S-306-00-002



PAEDIATRIC EXAM TABLES

PediaPals Paediatric Exam Tables are the #1 selling pediatric exam room tables world-wide. ZooPals Exam Tables are custom shaped and designed creating an environment of healing through comfort and familiarity to children. This allows for physicians, nurses, and pediatricians to accurately assess and consult with children and parents to diagnose and prescribe the best solutions. Customers will return, confident that you understand their children's needs.

SKU: PP.4XXXX



56 SERIES BARIATRIC RECLINER

SKU: CH.56/CH.56P/CH.566-2



54 SERIES RECLINER

SKU: CH.54/CH.54P/CH.547-2



OVERNIGHTER RECLINER

Founded on a strong steel frame, the 526 Overnighter Series patient chair is a hard-working sleeper that provides the greatest possible level of comfort to overnight patients. Mobile and impressively flexible given its size, the 526 has high-quality casters that make transport easy and non-intrusive.

It's open space, all-in-one design helps reduce claustrophobic sensations often caused by constant sitting. The quality foam and fabrics of the 526 Series make for real comfort, as well as a better-looking and longer-lasting product. Maximize patients' comfort by helping them choose between three positions. Years of research and testing have gone into the perfection of the Overnighter Series. It's specially designed to ensure that each has every opportunity to get a good night's sleep, whatever their medical situation.

Available with manual and power options in a wide range of colours.

SKU: CH.526-2



65 SERIES BARIATRIC RECLINER

The Ascent XL Bariatric Recliner is ruggedly built like our classic recliners, with added style and elegance. Gently contoured lines enhance its profile while comforting patients up to 500 lbs. Standard swing-arms provide easy cleaning, maintenance and patient transfer. Four-inch, twin wheel casters are corrosion resistant and roll with ease.

Suitable for Oncology, Surgery, Infusion, Dialysis Therapies, Orthopedics, Paediatrics, and more.

Purchase the Ascent XL patient chair through Midmed today for an easily transportable, well-designed recliner that always puts comfort first. Available with manual and power options in a wide range of colours.

SKU: CH.65/CH.65P/CH.660-2



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48



EXTRA ERGOLIFT STAIR CHAIR

Extra Ergolift is the powered version of the Extra track assisted stair chair with the same main structural features. The real innovation of Extra Ergolift is given by the electric motor that moves the tracks up and down the stairs without any effort by the operator. The motor is activated by a wireless console. The battery ensures a high autonomy and, in case of sudden discharging, this device will allow you to safely complete the service.

Extra Ergolift loads up to 200 kg in powered transport configuration and 230 kg in manual transportation. It can be used under rain or adverse weather conditions.

SKU: MR.671



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48



PATIENT TRANSFER SCALE

Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS



PAEDIATRIC/NICU EVACUATION SET

SKU: MSTOD2-Set



EXTRA ERGOLIFT STAIR CHAIR

Extra Ergolift is the powered version of the Extra track assisted stair chair with the same main structural features. The real innovation of Extra Ergolift is given by the electric motor that moves the tracks up and down the stairs without any effort by the operator. The motor is activated by a wireless console. The battery ensures a high autonomy and, in case of sudden discharging, this device will allow you to safely complete the service.

Extra Ergolift loads up to 200 kg in powered transport configuration and 230 kg in manual transportation. It can be used under rain or adverse weather conditions.

SKU: MR.671



PENGUIN NUTRITIONAL WARMER

Going above and beyond just warming breast milk, The Ameda Penguin Nutritional warmer ensure that each feeding is always safely delivered within the ideal human physiological temperature range of 36-37°C, regardless of the starting temperature. A feature which is unique to the Ameda Penguin Warmer is our gentle vibration technology which prevents potential nutritional degradation by continuously stirring the milk throughout the warming cycle for a more even and consistent warming. Our technology then eliminates the need to shake or roll the milk to even out any dangerous hot spots.

SKU: AM.PNW00115AU



FLOWCUP

The Materni flowcup is a simple, convenient accessory for feeding infants who have difficulty sucking.

SKU: MI.MEP

FLOWSPOON

A simple and convenient accessory for safely giving the first breast milk (colostrum).

SKU: MI.MEPAO112



NON-STERILE BABY BOTTLES

Materni bottles are suitable for the hygienic collection and storage of breast milk and the preparation of formula. They are bottle warmer, microwave, dishwasher and freezer safe making them highly versatile. Bottles have a screw on cap to ensure micro-biological cleanliness and have clear volume markings to make milk preparation quick and easy. They are available in four sizes: 50 ml, 80 ml, 120 ml and 240 ml.

SKU: MI.MEPAO109





STERILE BABY BOTTLES

Our disposable baby bottles come in a variety of configurations.

SKU: CA.NCB1240A



PROCEDURE CART

Procedure Cart for Healthcare with 6 Locking Drawers

- Large Procedure Cart for flexible storage
- New option to traditional metal carts
- Non-metal, rustproof polymer construction
- One-piece body for strength
- 6 full extension drawers (4 small, 2 large)
- Key lock secures all drawers
- Optional electronic keypad lock
- 121.92cm tall cart with 6 large drawers
- Easy to push or pull built in handle
- Durable chemical resistance
- Four 12.7cm Total Lock Casters (locks rotation and swivel)
- White polyethylene & ABS construction
- 6 full-extension drawers
- Four drawers are 7.62 inches deep
- Two drawers are 17.78 inches deep
- Built-in 3 inch deep tray top

SKU: TR.XXXX



AURION CARTS

Aurion Medical carts are innovative, high quality and built to last. These medical trolleys have a tough design, numerous compartments and are safe for use in any medical field. The Aurion range includes carts for use in emergency, anaesthetics, therapy, paediatrics, operating rooms and more.

SKU: AU.CP.



SHORT-TERM ENFIT FEEDING TUBES

Our unique PVC enteral tubes provide compatibility across the brand's range of enteral feeding products. They meet ISO 80369-3 standards and come with a number of distinct features.

- Readymade male EnFit™ connectors
- Numerical markings (in centimetres)
- Closed rounded tip
- 2 lateral distal openings
- 2 radiopaque stripes over transparent tube
- Available in a range of lengths and sizes from 4 Ch/Fr to 10 Ch/Fr and in a variety of widths

SKU: CA.NCE110P



PATIENT TRANSFER SCALE

Weigh Immobile Patients In Seconds

A revolutionary solution for weighing immobile patients without wasting valuable minutes.

The Patient Transfer Scale is the new way to achieve an accurate reading without compromise.

SKU: WB.PTS

INDEX

5

54 Series Recliner 12

56 Series Bariatric Recliner 8, 12

6

65 Series Bariatric Recliner 9, 13

A

Aurion Carts 20

B

Bendable IV ArmBoard 7

C

Catheter Holder Adhesive Patch 6

Catheter Holder Leg Band 7

Children's Bed - Teddy 11

D

Daniel Paediatric Bed 10

Disposable Curtains 1

E

EXTRA ERGOLIFT Stair Chair 14, 16

F

Flowcup 4, 17

Flowspoon 4, 18

I

IV Armor 8

N

Nasogastric Tube Holder 1, 6

Newborn Baby Bed 10

Non-Sterile Baby Bottles 5, 18

O

Overnighter Recliner 9, 13

P

Paediatric Exam Tables 11

Paediatric/NICU Evacuation Set 16

Patient Transfer Scale 15, 21

Penguin Nutritional Warmer 3, 17

pH Testing Strips 2

Procedure Cart 19

S

Short-Term ENFIT Feeding Tubes 3, 20

Single Use ENFIT Syringes 2

Standard Evacuation Sled 14, 15

Sterile Baby Bottles 5, 19