

TABLE OF CONTENTS

- Consumables** 2
 - Curtains & Wall Displays 2
 - Securement Devices 2
- Furniture** 3
 - Bariatric Furniture 3
 - Treatment Chairs 3
- Pre-Hospital** 4
 - Bags & Kits 4
 - Patient Handling 4
 - Patient Transfer 5
- Equipment** 6
 - Trolleys & Carts 6



DISPOSABLE CURTAINS

Each EcoMed disposable curtain is treated with the 3rd Generation Ultra-Fresh built-in antimicrobial protection.

SKU: EN.



BENDABLE IV ARMBOARD

The Bendable ArmBoard is a simple, cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

It can be used in the Operating Room, Recovery Room, ICU, Patient Care Units, and more.

SKU: DA.650



IV ARMOR

Midmed's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering. It has several benefits including reducing the introduction of infection associated with IV reinsertion, and reducing downstream occlusion alarms. The IV-ARMOR also reduces risk of catheter kinking in a bend such as a wrist or elbow, and reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm. IV-ARMOR uses skin-friendly adhesive, and comes in a one size fits most box of 10.

SKU: DA.675



STEPHEN H BARIATRIC TREATMENT CHAIR

The design, providing a high level of comfort for the patient, has been developed to provide a customised solution for the bariatric (high weight) patient care. The chair is able to move patients up to 300 kg. It is a multipurpose chair designed for hospital use during assessment, diagnosis, care and treatment. It is suitable for many different applications and may also be used to transport the patient from one department to another.

The backrest is produced with sufficient stability to carry out CARDIO PULMONARY RESUSCITATION. This medical device, made with stainless materials, represents the "ultimate" solution regarding the elimination of rust, thus significantly reducing the likelihood of bacterial contamination.

SKU: GB.SHA.STANDARD



STEPHEN H BARIATRIC TREATMENT CHAIR

The design, providing a high level of comfort for the patient, has been developed to provide a customised solution for the bariatric (high weight) patient care. The chair is able to move patients up to 300 kg. It is a multipurpose chair designed for hospital use during assessment, diagnosis, care and treatment. It is suitable for many different applications and may also be used to transport the patient from one department to another.

The backrest is produced with sufficient stability to carry out CARDIO PULMONARY RESUSCITATION. This medical device, made with stainless materials, represents the “ultimate” solution regarding the elimination of rust, thus significantly reducing the likelihood of bacterial contamination.

SKU: GB.SHA.STANDARD

SMARTSAFE VALUABLE STORAGE BAGS

SKU: TS.7325(0/1)





STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48



STAIR CHAIR

SKU: B.STAIRCHAIR.NEW



STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

SKU: MS.28, MS.36, MS.48



PROCEDURE CART

Procedure Cart for Healthcare with 6 Locking Drawers

- Large Procedure Cart for flexible storage
- New option to traditional metal carts
- Non-metal, rustproof polymer construction
- One-piece body for strength
- 6 full extension drawers (4 small, 2 large)
- Key lock secures all drawers
- Optional electronic keypad lock
- 121.92cm tall cart with 6 large drawers
- Easy to push or pull built in handle
- Durable chemical resistance
- Four 12.7cm Total Lock Casters (locks rotation and swivel)
- White polyethylene & ABS construction
- 6 full-extension drawers
- Four drawers are 7.62 inches deep
- Two drawers are 17.78 inches deep
- Built-in 3 inch deep tray top

SKU: TR.XXXX



AURION CARTS

Aurion Medical carts are innovative, high quality and built to last. These medical trolleys have a tough design, numerous compartments and are safe for use in any medical field. The Aurion range includes carts for use in emergency, anaesthetics, therapy, paediatrics, operating rooms and more.

SKU: AU.CP.

INDEX

A

Aurion Carts 6

B

Bendable IV ArmBoard 1

D

Disposable Curtains 1

I

IV Armor 2

P

Procedure Cart 5

S

SmartSafe Valuable Storage Bags 3

Stair Chair 4

Standard Evacuation Sled 4, 5

Stephen H Bariatric Treatment Chair 2, 3