

# TABLE OF CONTENTS

- Consumables** ..... 2
  - Airway Management ..... 2
  - Fluid Warming ..... 5
  - Hygiene ..... 7
  - Patient Warming ..... 7
  - Securement Devices ..... 8
  - Suction ..... 10
- Furniture** ..... 13
- Pre-Hospital** ..... 13
  - Airway Management ..... 13
  - Bags & Kits ..... 15
  - Bleeding Control ..... 18
  - Patient Handling ..... 18
  - Patient Warming ..... 22
  - Splints ..... 23
  - Triage ..... 26
- Equipment** ..... 29
  - Patient Transfer ..... 29
  - Refrigeration ..... 31

## DUCANTO CATHETER™



The SSCOR DuCanto Catheter is a new catheter that is unlike traditional Yankauer suction tips that are designed for controlled surgical settings. The SSCOR DuCanto Catheter was developed with both routine and emergency airway management in mind. It has a larger internal diameter making it significantly less likely to clog and is shaped to fit the oropharyngeal anatomy. The hyper curved catheter is shaped more like a patient's airway making it easy to position the catheter.

Additionally, the SSCOR DuCanto Catheter does not have a thumb port. By removing the vacuum control, the caregiver does not need to locate and continuously occlude the thumb port to clear a patient's airway. This allows continuous suction and keeps hands free.

**SKU:** SCR.200-00002C

## HI-D BIG STICK SUCTION TIP



The SSCOR Hi-D Big Stick Suction Tip aims to combat the most common complaint for those using suction equipment for resuscitation procedures... clogged tips.

**SKU:** SCR.44241



## **VX-2® SUCTION UNIT**

**SKU:** SCR.2310



## **DUET SUCTION UNIT**

**SKU:** SCR.2314-230




## QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

**Available in both rechargeable and battery operated variations**

**SKU:** SCR.2403-230



## PYLANT MONITOR - ENDOTRACHEAL TUBE CUFF PRESSURE MANOMETER

**SKU:** INO.IM.100



## TRACHEOSTOMY TUBE HOLDER

**SKU:** DA.240



## ENDOTRACHEAL TUBE HOLDER

**SKU:** DA.273

## LIFEBOX MEDICAL BOX

The next generation of EMS transport containers is here. The key lies within the Carbon panels that produce insulation values of R50; 2x more powerful than anything else on the market. This allows for more payload space and greater efficiency in both hot and cold environments.

**SKU:** F.6001024



## QINFLOW WARRIOR - BLOOD/FLUID WARMER

**SKU:** QF.Q211050000SUC



## QINFLOW WARRIOR LITE - BLOOD/FLUID WARMER

High performance yet compact and light blood and IV fluid warmer for space and weight constrained rescue gears, first responders, and critical care transports. It has field proven technology and affordable consumables with a unique continuum of emergency care proposition.

**SKU:** QF.QIF03-BUA1000SUC

## WRASFER™ TRANSFER SHEET

The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

**SKU:** WRASFER



## **QINFLOW WARRIOR - BLOOD/FLUID WARMER**

**SKU:** QF.Q211050000SUC



## **QINFLOW WARRIOR LITE - BLOOD/FLUID WARMER**

High performance yet compact and light blood and IV fluid warmer for space and weight constrained rescue gears, first responders, and critical care transports. It has field proven technology and affordable consumables with a unique continuum of emergency care proposition.

**SKU:** QF.QIF03-BUA1000SUC





## BENDABLE IV ARMBOARD

The Bendable ArmBoard is a simple, cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

It can be used in the Operating Room, Recovery Room, ICU, Patient Care Units, and more.

**SKU:** DA.650



## IV ARMOR

Midmed's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering. It has several benefits including reducing the introduction of infection associated with IV reinsertion, and reducing downstream occlusion alarms. The IV-ARMOR also reduces risk of catheter kinking in a bend such as a wrist or elbow, and reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm. IV-ARMOR uses skin-friendly adhesive, and comes in a one size fits most box of 10.

**SKU:** DA.675



## TRACHEOSTOMY TUBE HOLDER

**SKU:** DA.240



## ENDOTRACHEAL TUBE HOLDER

**SKU:** DA.273

## DUCANTO CATHETER™



The SSCOR DuCanto Catheter is a new catheter that is unlike traditional Yankauer suction tips that are designed for controlled surgical settings. The SSCOR DuCanto Catheter was developed with both routine and emergency airway management in mind. It has a larger internal diameter making it significantly less likely to clog and is shaped to fit the oropharyngeal anatomy. The hyper curved catheter is shaped more like a patient's airway making it easy to position the catheter.

Additionally, the SSCOR DuCanto Catheter does not have a thumb port. By removing the vacuum control, the caregiver does not need to locate and continuously occlude the thumb port to clear a patient's airway. This allows continuous suction and keeps hands free.

**SKU:** SCR.200-00002C

## HI-D BIG STICK SUCTION TIP



The SSCOR Hi-D Big Stick Suction Tip aims to combat the most common complaint for those using suction equipment for resuscitation procedures... clogged tips.

**SKU:** SCR.44241



## **VX-2® SUCTION UNIT**

**SKU:** SCR.2310



## **DUET SUCTION UNIT**

**SKU:** SCR.2314-230



## QUICKDRAW® PORTABLE SUCTION UNIT

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

**Available in both rechargeable and battery operated variations**

**SKU:** SCR.2403-230



## VX-2® SUCTION UNIT

**SKU:** SCR.2310



## **DUET SUCTION UNIT**

**SKU:** SCR.2314-230



## **QUICKDRAW® PORTABLE SUCTION UNIT**

Combining the convenience of handheld suction with the performance of powered suction. The most efficient suction pump you have ever held in your hand. At only 1.2 kg, the SSCOR Quickdraw® is more efficient than the most powerful traditional suction units. It is compact and light enough to fit in an airway bag and be readily transported to the patient's side on every call.

**Available in both rechargeable and battery operated variations**

**SKU:** SCR.2403-230



## SMART XTRACTSR STRETCHER

The Smart XtractSR Stretcher provides safe and rapid casualty extraction from hazardous locations, reducing exposure to hostile locations and secondary injury.

This stretcher is extremely robust despite it's lightweight material: it weighs less than 1.5 kg in total and can safely carry up to 300 kg in weight. It can be carried in operational equipment, making it the most convenient and easy-to-use stretcher for rescuers and patients in temperamental weather or terrain conditions.

**SKU:** TS.72121



## SMART TRIAGE PACK

**SKU:** TS.73212-AUS



## SMART COMMANDER

The SMART Commander is a triage system designed for the structured response and management of multiple casualty incidences. It is sized to fit your arm, allowing you to move around the scene and providing you with a solid surface to write on. The waterproof control boards can be written on with a sharpie marker, providing a durable record of actions on scene that can be cleaned post incident with alcohol based cleaner or wipes. It is built for mass casualty environments, as the rugged pack is suitable for all weather and environmental conditions. The SMART Commander features all of the essential scene management and response tools.

**SKU:** TS.73211



## SMART TRIAGE BAG

The SMART T.R.I. Pack has been designed around the equipment most frequently selected to address immediately life threatening conditions typical at a multiple casualty incident. Whether you are looking to ensure MUCC Compliance or improve your ability to respond to multiple trauma patients the SMART T.R.I. Pack will help you ensure you have the right equipment in the right place at the right time.

Guided by customer feedback and data from recent incidents, we have designed the SMART T.R.I. Pack to be a valuable resource to the responder, providing them with the space for medical equipment such as trauma dressings, airways and chest seals, without the need to take a complete trauma bag into the scene. This allows the delivery of Rapid Life Saving interventions at point of wounding.

**SKU:** TS.TRI-PK-AUS





## STOR-IT BAG

The SMART Stor-It Reflective Triage Bag is an extremely well organised solution for bulk storage and rapid deployment of medical equipment. Made to the highest specification the Stor-It is designed specifically for disaster and emergency response but has been shown to be highly effective at aide stations during mass gathering events.

The rugged storage system has triple reinforced stitching and exterior padding to protect fragile medical equipment during storage, transport and shipment. It is easy to stack multiple units and the inside divides into two easy to hang sections with highly visible contents.

**NOTE:** Bag Only

**SKU:** TS.STOR-IT.



## SMART MCI BAG

Midmed's SMART Triage products offer the definitive collection of items that all triage professionals need. In high-pressure emergency situations that call for triage, you can feel secure in calling on their products to seamlessly aid your decision-making process and save you valuable time. TSG have developed their incredible resources with the express aim of giving you quick, easy access to key information.

As an all-inclusive emergency kit, the SMART MCI Bag™ represents this aim most comprehensively. It's a complete grab bag. Ideal for use at large planned events where emergency help is needed, the SMART MCI Bag™ will provide the perfect base for your medical command post and is capable of treating up to 80 casualties.

**SKU:** TS.73210-AUS

## SMARTSAFE VALUABLE STORAGE BAGS



SMARTSAFE small property bags are an effective way of increasing patient comfort and security when they need it most. These bags can be used to store many items, including valuables such as jewellery or medication.

An evident tamper seal and barcode are used to track the bag's location, and for the patient's peace of mind, the barcode is attachable to patient notes. The bags also include a wristband that is to be placed around the patient's wrist. Clearly labelled forms further help identify the belongings of each patient. The bags also include a wristband that is to be placed around the patient's wrist.

**SKU:** N/A

## RESPONDER COMPRESSOR BANDAGE



The FareTec Responder Compression Bandages are designed for quick and easy application by caregivers to provide effective, multi-functional treatment in emergencies. The bandages are sterile, non-adherent pads that apply pressure to any site and are easily wrapped and secured. They restrict blood flow to aid in the treatment of traumatic bleeds. Available in three sizes: 4 inch, 6 inch and 8 inch widths.

**SKU:** F.222109(6/8)



## TACTICAL RESCUE SLED

**SKU:** MED.MS30



## COMBICARRIER® II

The CombiCarrier® II is the simplest and most ergonomic backboard on the market.

Working as a fully functional spinal immobilisation device, as well as a scoop stretcher and backboard, it allows for safe, easy extrication of injured patients. Its updated version further improves the inimitable 2-in-1 design of the original, and brings it fully up to date with a thinner, lower profile, allowing the whole board to fit inside all standard storage compartments.

Design is everything with the CombiCarrier® II. The concave surface ensures patients' lateral support and comfort, and its clear centre section makes it completely X-ray compatible, providing a clear view of the patient's spine. It comes with four restraint straps, but no pesky bits and pieces like rivets or screws. Just a clean, minimal, brilliant design.

Its simple intuitive design minimises the training requirements of EMS and hospital personnel, and its lower profile fits in standard storage compartments. Start using a CombiCarrier® II now to make the extrication of patients easier. The revolutionary medical equipment saves on maintenance and replacement costs.

**SKU:** HM.CC2200PY



## SMART XTRACTSR STRETCHER

The Smart XtractSR Stretcher provides safe and rapid casualty extraction from hazardous locations, reducing exposure to hostile locations and secondary injury.

This stretcher is extremely robust despite its lightweight material: it weighs less than 1.5 kg in total and can safely carry up to 300 kg in weight. It can be carried in operational equipment, making it the most convenient and easy-to-use stretcher for rescuers and patients in temperamental weather or terrain conditions.

**SKU:** TS.72121



## BANANA SLIDE SINGLE PATIENT TUBULAR

The GBUK Banana Single Patient Tubular Slides are strong and durable with a high slip performance. It has a strong, red, hemmed slide sheet and a silicone coat on both sides. With a low profile seam, they are long-lasting and have high-slip characteristics, enabling easy and safe patient transfers and re-positioning. Available in two sizes: 72 cm by 72 cm and 100 cm by 200 cm.

**SKU:** N/A

## MEGAMOVER® TRANSPORT UNIT




The MegaMover offers a great solution for transporting patients from tight, confined areas that are inaccessible to stretchers. The unit has 14 handles for ergonomic lifting, and can hold up to 450 kg (1,000 pounds). It is 80" (203 cm) x 40" (101 cm), and made of durable, latex-free non-woven/poly material.

The MegaMover Transport Unit is **not** designed for use when there is need for spinal stabilisation or traction of the patient being transported. It is recommended to use at least four people to provide safe transport/transfer, thus reducing the possibility of accident or injury.

**SKU:** 65459

## WRASFER™ TRANSFER SHEET



The WRASFER™ Patient Transfer wrap sheet is designed for swiftly and hygienically transporting immobile or injured patients in high-pressure emergency situations. WRASFER™ disposable patient wrap and transfer sheet is designed for use in paramedic services, hospital emergency wards and operating theatres, in situations where fluids need to be easily absorbed.

With a fluid absorption capacity of 5 litres and a max. carrying weight of 200 KG, WRASFER™ is an extremely versatile product. It addresses four critical issues: contamination control, cleaning time/ costs, ambulance turnaround times, and patient transfer.

**SKU:** WRASFER



## PORTABLE SURGICAL TABLE

Midmed stocks the Portable Surgery Table (PST), a patient transport and operation platform. Field tested and outfitted with the most advanced materials available, this unit can serve a multitude of purposes. The PST quickly assembles into a stable patient aid platform for medics and doctors. The device is versatile and easy to transport for fast changing scenarios.

**SKU:** F.2221027



## LIFEBOX MEDICAL BOX

The next generation of EMS transport containers is here. The key lies within the Carbon panels that produce insulation values of R50; 2x more powerful than anything else on the market. This allows for more payload space and greater efficiency in both hot and cold environments.

**SKU:** F.6001024



## **QINFLOW WARRIOR - BLOOD/FLUID WARMER**

**SKU:** QF.Q211050000SUC



## **QINFLOW WARRIOR LITE - BLOOD/FLUID WARMER**

High performance yet compact and light blood and IV fluid warmer for space and weight constrained rescue gears, first responders, and critical care transports. It has field proven technology and affordable consumables with a unique continuum of emergency care proposition.

**SKU:** QF.QIF03-BUA1000SUC



## CT7 FEMORAL TRACTION SPLINT

After 2 years of development, the next-generation CT-7 is finally available. It has been tested under the harshest conditions around the world by an array of first responders and medics. The CT-7 femoral traction splint has been simplified by removing straps, velcro and all loose ends. It has a smaller footprint and no longer extends beyond the patient's foot. It has been made more rigid and indestructible with thicker carbon tubes. The traction mechanism is easy, fast and precise by utilizing an acme threaded polycarbonate rod system.

**SKU:** F.1126629



## ADULT PELVIC BINDER - YELLOW

Simple, Fast, Easy – PelvicBinder® for pelvic fracture stabilization. Pelvic Binders are an essential piece of a multidisciplinary approach which has been shown to save lives in a case of a pelvic fracture.

- Radiolucent
- MRI Compatible
- Latex Free
- Made in the U.S.A.

The Adult Standard PelvicBinder is double the width and length of the Pediatric PelvicBinder. The Pediatric Pelvic Binder is designed for Hip Circumferences ranging between 14" to 32".

**SKU:** PB.Y-3000



## FASPLINT EXTREMITY SPLINT



The FASPLINT is an innovative type of extremity vacuum splint that contours around the injured area, evacuates air and molds perfectly around the damaged extremity. From dislocated shoulders to angulated wrist fractures, simply grab your vacuum splint, place it around the injured area, evacuate the air, and you're done! And if equipment retrieval is a major factor facing your department, then the FASPLINT is ideal. It is easy-to-use and lightweight to carry, making it the simplest kind of splint to use in the event of a medical emergency.

**SKU:** HM.FS80

## EVAC-U-SPLINT



The Evac-U-Splint provides a solution for almost every conceivable splinting situation. From dislocated shoulders to angulated wrist fractures, simply grab your vacuum splint, form it around the injured area, evacuate the air and you're done! You can be transporting your patient to the hospital while other rescuers are still bending, padding, and struggling to apply traditional fixation splints. Versatile and adaptable, the number of uses for the Evac-U-Splint is limited only by the creativity of the rescuer.

It is a very effective and reliable medical extremities splint. This device extracts the air in order for the material to mold around the injured limb and to the exact shape of injury without applying pressure. It is available in a small, medium, and large size.

**SKU:** EV3000.SR

## CT6 FEMORAL TRACTION SPLINT



The CT6 Femoral Traction Leg Splint is a compact, powerful, reliable device that helps relieve pain and prevents further damage to surrounding muscle and tissue while reducing the risk of death caused from severed arteries. It has been designed for pre hospital and patient transport. It is primarily used on patients with mid-line femur fractures. It can be used on adult and pediatric patients and is extremely light, compact and strong.

**SKU:** F.1126628P



## SMART TAPE

SMART Tape is designed specifically for children at a multiple casualty scene. This evidence-based method of assessing patient priority is simple and easy to use in high stress emergency situations. Adult based assessment methods are not suitable for the triage of children. Emotional reaction to an incident with multiple child casualties often sees children being over prioritised.

SMART Tape is specifically designed to assess the severity of pediatric casualties. Triage accuracy using SMART Tape allows for the fast identification of those most critical, both at and away from the scene, and prevents the misdirection of vital triage resources. SUITED for 5 different Age ranges

**SKU:** TS.73214



## SMART TAGS

SMART Tags are highly durable, waterproof and tear-resistant and suited to triage treatment and tracking. With SMART Tags effective triage is quick and simple, allowing the responder to reprioritise the casualty quickly without losing information.

**SKU:** TS.73215-20



## SMART TRIAGE PACK

**SKU:** TS.73212-AUS



## SMART COMMANDER

The SMART Commander is a triage system designed for the structured response and management of multiple casualty incidences. It is sized to fit your arm, allowing you to move around the scene and providing you with a solid surface to write on. The waterproof control boards can be written on with a sharpie marker, providing a durable record of actions on scene that can be cleaned post incident with alcohol based cleaner or wipes. It is built for mass casualty environments, as the rugged pack is suitable for all weather and environmental conditions. The SMART Commander features all of the essential scene management and response tools.

**SKU:** TS.73211



## SMART TRIAGE BAG

The SMART T.R.I. Pack has been designed around the equipment most frequently selected to address immediately life threatening conditions typical at a multiple casualty incident. Whether you are looking to ensure MUCC Compliance or improve your ability to respond to multiple trauma patients the SMART T.R.I. Pack will help you ensure you have the right equipment in the right place at the right time.

Guided by customer feedback and data from recent incidents, we have designed the SMART T.R.I. Pack to be a valuable resource to the responder, providing them with the space for medical equipment such as trauma dressings, airways and chest seals, without the need to take a complete trauma bag into the scene. This allows the delivery of Rapid Life Saving interventions at point of wounding.

**SKU:** TS.TRI-PK-AUS



## STOR-IT BAG

The SMART Stor-It Reflective Triage Bag is an extremely well organised solution for bulk storage and rapid deployment of medical equipment. Made to the highest specification the Stor-It is designed specifically for disaster and emergency response but has been shown to be highly effective at aide stations during mass gathering events.

The rugged storage system has triple reinforced stitching and exterior padding to protect fragile medical equipment during storage, transport and shipment. It is easy to stack multiple units and the inside divides into two easy to hang sections with highly visible contents.

**NOTE:** Bag Only

**SKU:** TS.STOR-IT.



## SMART MCI BAG

Midmed's SMART Triage products offer the definitive collection of items that all triage professionals need. In high-pressure emergency situations that call for triage, you can feel secure in calling on their products to seamlessly aid your decision-making process and save you valuable time. TSG have developed their incredible resources with the express aim of giving you quick, easy access to key information.

As an all-inclusive emergency kit, the SMART MCI Bag™ represents this aim most comprehensively. It's a complete grab bag. Ideal for use at large planned events where emergency help is needed, the SMART MCI Bag™ will provide the perfect base for your medical command post and is capable of treating up to 80 casualties.

**SKU:** TS.73210-AUS

## DECON / HAZMAT SLED



The Med Sled® Decon/Hazmat Sled ensures easy movement of a non-ambulatory person through the Decon Shower process and safely transports to a cleared zone. (Not used for vertical evacuation.)

**SKU:** MS.36DC

## VERTICAL LIFT RESCUE SLED



The Med Sled® Vertical Lift Rescue Sled provides all of the functionality of conventional CASEVAC/Tactical Evacuation Care methods with a less complicated, more durable system, all for a great price.

**SKU:** MS.36VL



## STANDARD EVACUATION SLED

When timing is critical, the Med Sled® Evacuation Sled enables a person to evacuate another person out of a single-story building or down a stairwell of a multi-floor building. Stairwell braking system included on every sled. Evacuations - Made Easier.

**SKU:** MS.28, MS.36, MS.48



## LIFEBOX MEDICAL BOX

The next generation of EMS transport containers is here. The key lies within the Carbon panels that produce insulation values of R50; 2x more powerful than anything else on the market. This allows for more payload space and greater efficiency in both hot and cold environments.

**SKU:** F.6001024

# INDEX

## A

Adult Pelvic Binder - Yellow 23

## B

Banana Slide Single Patient Tubular 19

Bendable IV ArmBoard 8

## C

CombiCarrier® II 18

CT6 Femoral Traction Splint 25

CT7 Femoral Traction Splint 23

## D

Decon / Hazmat Sled 29

DuCanto Catheter™ 1, 10

Duet Suction Unit 2, 11, 13

## E

Endotracheal Tube Holder 4, 9

Evac-U-Splint 24

## F

Fasplint Extremity Splint 24

## H

Hi-D Big Stick Suction Tip 1, 10

## I

IV Armor 8

## L

LifeBox Medical Box 5, 21, 30

## M

MegaMover® Transport Unit 20

## P

Portable Surgical Table 21

Pylant Monitor - Endotracheal Tube Cuff Pressure Manometer 3

## Q

QinFlow Warrior - Blood/Fluid Warmer 5, 7, 22

QinFlow Warrior Lite - Blood/Fluid Warmer 6, 7, 22

Quickdraw® Portable Suction Unit 3, 12, 13

## R

Responder Compressor Bandage 17

## S

SMART Commander 15, 27

SMART MCI Bag 16, 28

SMART Tags 26

SMART Tape 25

SMART Triage Bag 15, 27



SMART Triage Pack 14, 26  
SMART XtractSR Stretcher 14, 19  
SmartSafe Valuable Storage Bags 17  
Standard Evacuation Sled 30  
Stor-It Bag 16, 28

## **T**

Tactical Rescue Sled 18  
Tracheostomy Tube Holder 4, 9

## **V**

Vertical Lift Rescue Sled 29  
VX-2® Suction Unit 2, 11, 12

## **W**

WRASFER™ Transfer Sheet 6, 20