

TABLE OF CONTENTS

- Gardhen Balance** 2
 - Gardhen Balance Dialysis Chairs 2
- Furniture** 3
 - Bariatric Furniture 3
 - Beds 4
 - Stools 7
 - Treatment Chairs 7
- Equipment** 13
 - Trolleys & Carts 13
- Specials** 14
- Cair & Micsafe Specials** 14
- Medical Furniture & Equipment Specials** 15



STEPHEN H COMFORT 9.0

SKU: N/A



UP TREATMENT CHAIR

Thanks to an advanced geometric solution, the “leg” section disappears under the seat. This allows the patient to sit on the chair with the lower limbs perfectly perpendicular to the seat itself. Comfort and postural benefits will derive from this position. Therefore, all patients, even those with motor difficulties, can easily sit down and get up from the chair, even from the front, always with a correct posture. UP chair facilitates the patient’s exit with a correct postural approach. Backrest, seat and leg sections move simultaneously and automatically, raising the patients a semi-erect position according to their height. The device is equipped with three memory functions indeed, each of these positions facilitates the exit of patients according to their size. The following images represent an example of the Patient’s “Lift” or Verticalization function.

SKU: GB.UP



STEPHEN H ANATOMICAL TREATMENT CHAIR

Beyond the business class

SKU: GB.SHA



STEPHEN H COMFORT 7.0 TREATMENT CHAIR

Multifunctional chair

SKU: GB.SH.COMFORT.7.0



STEPHEN H ANATOMICAL

SKU: GB.SHA-




STEPHEN H BARIATRIC TREATMENT CHAIR

The design, providing a high level of comfort for the patient, has been developed to provide a customised solution for the bariatric (high weight) patient care. The chair is able to move patients up to 300 kg. It is a multipurpose chair designed for hospital use during assessment, diagnosis, care and treatment. It is suitable for many different applications and may also be used to transport the patient from one department to another.

The backrest is produced with sufficient stability to carry out CARDIO PULMONARY RESUSCITATION. This medical device, made with stainless materials, represents the “ultimate” solution regarding the elimination of rust, thus significantly reducing the likelihood of bacterial contamination.

SKU: GB.SHA.STANDARD



DANIEL-MINI PAEDIATRIC BED

SKU: N/A



CLOUD BARELLA ELETTRICA

SKU: N/A



MOSE MINI PAEDIATRIC BED

Italian design and technology in the world

SKU: GB.MINIMOSE



ALEX HOSPITAL BED - TILT

Italian design and technology in the world.

SKU: GB.ALEX



DANIEL PAEDIATRIC BED

The Gardhen Balance Daniel Paediatric Bed has an accurate, integrated weighing system. Weight display automatically displays weight. It has an anti-blackout function with rechargeable long life batteries, an automatic and manual tare mode and easy storage and transfer of data. The bed has a quick bilateral release in the back section for CPR, an ergonomic patient surface and a bed section protected with HPL panels (easy to be cleaned).

SKU: GB.DANIEL



SURGIMOVE SURGEON CHAIR

SurgiMove is perfectly suited for a wide range of surgical applications. Its advanced design delivers multiple positions including 360 degree adjustable armrests that can be provided as optional. Pedals are user friendly, the backrest and the armrests are fully adjustable allowing to the SurgiMove to reach superior comfort and adjustability.

SKU: GB.SURGIMOVE



UP TREATMENT CHAIR

Thanks to an advanced geometric solution, the “leg” section disappears under the seat. This allows the patient to sit on the chair with the lower limbs perfectly perpendicular to the seat itself. Comfort and postural benefits will derive from this position. Therefore, all patients, even those with motor difficulties, can easily sit down and get up from the chair, even from the front, always with a correct posture. UP chair facilitates the patient’s exit with a correct postural approach. Backrest, seat and leg sections move simultaneously and automatically, raising the patients a semi-erect position according to their height. The device is equipped with three memory functions indeed, each of these positions facilitates the exit of patients according to their size. The following images represent an example of the Patient’s “Lift” or Verticalization function.

SKU: GB.UP



THEIA 2.0 SURGERY CHAIR

Italian design and technology in the world

SKU: GB.THEIA2.0



THEIA 1.0 SURGERY CHAIR

Italian design and technology in the world

SKU: GB.THEIA1.0



STEPHEN H FREEDOM TREATMENT CHAIR

Feel free like at home

SKU: GB.FREEDOM



STEPHEN H FLAT TREATMENT CHAIR

Beyond the business class

SKU: GB.SHFLAT



STEPHEN H ANATOMICAL TREATMENT CHAIR

Beyond the business class

SKU: GB.SHA



STEPHEN H ECO 2.0 TREATMENT CHAIR

Highly performing device

SKU: GB.SH.ECO.2.0



STEPHEN H COMFORT 7.0 TREATMENT CHAIR

Multifunctional chair

SKU: GB.SH.COMFORT.7.0



STEPHEN H ECO CHAIR

An economical solution suitable for multiple situations, it may also be used for the transportation of patients between departments. The armchairs are electrically adjustable and have three separate cushions. The backrest is designed to support CARDIO PULMONARY RESUSCITATION. Thick-gauge stainless steel profiles, carefully manufactured and assembled, make up the strong frame that supports the cushions in an ergonomic shape.

SKU: GB.ECO



STEPHEN H ANATOMICAL

SKU: GB.SHA-



STEPHEN H BARIATRIC TREATMENT CHAIR

The design, providing a high level of comfort for the patient, has been developed to provide a customised solution for the bariatric (high weight) patient care. The chair is able to move patients up to 300 kg. It is a multipurpose chair designed for hospital use during assessment, diagnosis, care and treatment. It is suitable for many different applications and may also be used to transport the patient from one department to another.

The backrest is produced with sufficient stability to carry out CARDIO PULMONARY RESUSCITATION. This medical device, made with stainless materials, represents the “ultimate” solution regarding the elimination of rust, thus significantly reducing the likelihood of bacterial contamination.

SKU: GB.SHA.STANDARD

FREDDY THERAPY CART

Compact size, able to adapt to all situations

SKU: GB.FREDDYTHERAPY



FREDDY CART

Small but capacious easy to handle and upgradeable

SKU: GB.FREDDY



STEPHEN H ANATOMICAL

SKU: GB.SHA-






STEPHEN THERAPY CHAIR

SKU: N/A



STEPHEN H COMFORT 9.0

SKU: N/A



DANIEL-MINI PAEDIATRIC BED

SKU: N/A



CLOUD BARELLA ELETTRICA

SKU: N/A

INDEX

A

Alex Hospital Bed - Tilt 5

C

Cloud Barella Elettrica 4, 15

D

Daniel Paediatric Bed 6

Daniel-Mini Paediatric Bed 4, 15

F

Freddy Cart 13

Freddy Therapy Cart 12

M

Mose Mini Paediatric Bed 5

S

Stephen H Anatomical 3, 11, 13

Stephen H Anatomical Treatment Chair 2, 9

Stephen H Bariatric Treatment Chair 3, 12

Stephen H Comfort 7.0 Treatment Chair 2, 10

Stephen H Comfort 9.0 1, 14

Stephen H ECO 2.0 Treatment Chair 10

Stephen H ECO Chair 11

Stephen H Flat Treatment Chair 9

Stephen H Freedom Treatment Chair 8

Stephen Therapy Chair 14

Surgimove Surgeon Chair 6

T

Theia 1.0 Surgery Chair 8

Theia 2.0 Surgery Chair 7

U

UP Treatment Chair 1, 7