

## Portering Chair - Bariatric, Rear Steer, Sliding Footrest

### Features

- Heavy duty bariatric portering chair
- Seamless & contoured upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & porters
- Sliding footrest aids access as well as supporting feet during transit
- X-ray / notes pocket on rear of seat
- Hook for patients' belongings
- Fittings to accept infusion poles
- Retractable diagonal seatbelt
- 200mm fixed wheels / 200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg • Maximum capacity 380kg

### Dimensions (w x d x h)

- Overall Dimensions - 930 x 1000 x 1090mm
- Seat Width - 740mm
- Seat Depth - 670mm
- Seat Height - 530mm
- Backrest Height - 480mm

### Standard Colours

- Paint - Grey White
- Upholstery - Bristol Blue

### Options / Accessories

- BM.G/CLS/500 - Calf / Leg Supports (pair)
- BM.G/DP - Infusion Pole
- BM.G/CARD/STATION - Head of Station
- BM.G/CARD/LOCK - Card / Coin Lock, Factory Fitted
- BM.2LOKE011 - Release Card

Alternative colours are available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer Portering Chairs - Bariatric, Sliding Foot Rest in Aquamarine - G/500/RS/VAQ

### Vinyl - Ambla New Richmond in Colour



BM.G/500/RS/VBB

Item	Description	Weight
BM.G/500/RS/VBB	Portering Chair - Bariatric, Rear Steer, Sliding Footrest	43.8kg