

# Double Electric Breast Pump

with Large Tote and Accessories

The hospital-strength Mya Joy has been designed with versatile features, safety and comfort, offering the best value for moms who are looking for a convenient pumping experience. Compact, lightweight and quiet, the Mya Joy also features the world's only FDA-cleared pump kit that helps protect breastmilk from contaminants. This pump delivers on the ultimate promise of comfort by offering a wide range of settings, mom-friendly features, and flange sizes. Whether pumping on-the-go, at home or at work, Mya Joy offers pumping without compromise.



### Safe

Ameda's HygieniKit true closed system not only helps protect breastmilk from backing up into the tubing but also prevents the system from potential contaminants like bacteria, viruses, and mold growth



### Comfort

Wide range of suction and speed settings, functions, and flange sizes



### Ultra-Quiet

Operates at less than 45dB and allows for discreet pumping anywhere



### Powerful

Hospital Strength Performance and two modes of pumping



### Portable

Battery operated, lightweight and clips right onto a waist band, purse pocket or favorite yoga pants



SKU 131M22  
UPC 810725022160  
SKU 131T22  
UPC 810725022177

### Mya Joy Features:

- Hospital strength performance
- Up to 250 mmHg suction strength
- Backlit digital LCD display
- Timer display easily tracks pumping time
- Smart auto power off shuts down after 45 minutes
- 2-phase adjustable suction modes:
  - Stimulation (6 Levels)
  - Expression Mode (12 Levels)
- Manual or Automatic 2-phase, 2-minute auto switch from stimulation to expression
- Double or single electric pumping options

### Contents Include:

- Ameda Mya Joy Breast Pump
- AC Power Adapter (US 131X07, AU 131X10)
- Dual HygieniKit® Milk Collection System
- CustomFit™ Flanges (25mm, 28.5mm)
- 4 Extra Valves (for HygieniKit)
- 2 Extra Diaphragms (for HygieniKit)
- 2 Bottle Holders
- Tote Bag
- Accessory Bag (for Parts)
- Manual Breast Pump Handle
- Breast Milk Storage Guidelines
- Quick Start Guide