MEDI©GUARDPREMIUM



Pharmacy Refrigerator

MG300PR &MG400PR



Frequency Conversion Energy Saving

Variable frequency compressor with optimized cooling system, more energy efficient.



Environmentally Friendly

Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine which is super friendly to the environment



Large Space, Superior Uniformity

The space is efficiently used, and the utilization rate is high and the temperature uniformity is +3 C.



Low Noise

Optimized system and noise reduction design to reduce the noise level to 43dB(A) (220V/50Hz).



Stable performance, multiple alarm modes

Three alarm modes (light sound remote); six types of alarm



Microprocessor Control System



Self Closing Door with 90° Stay Open Feature



Precise Temperature Control



Security lock



Optional USB and Printer





Factory 1

13-15 Kevlar Close Braeside, VIC 3195 1300 459 140

Scope of Application



MEDI©GUARD PREMIUM


This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at $2 \sim 8$ °C.

Specifications



MG300PR

Technical Data			ENTA	0
Cabinet Type	Upright, g	lass door	EATING PUALITY ACCRE	/
Climate Class	N			
Cooling Type	Forced ai	r cooling		
Defrost Mode	Au	to		
Refrigerant	Н	С		
Noise (dB (A))	43	43		



Control	
Controller	Microprocessor
Display	LED

Electrical Data			
Power Supply (V/Hz)	220~240/50&60	115/60	
Power (W)	350	200	
Electrical Current (A)	1.8	2.2	

Dimensions			
Capacity (L/Cu.Ft)		290/10.2	
Net/Gross	(kg)	105/127	
Weight (approx)	(lbs)	231.6/280.1	
Interior	(mm)	530*555*1080	
Dimension (W*D*H)	(in)	20.9*21.9*42.5	
Exterior Dimension (W*D*H)	(mm)	665*710*1665	
	(in)	26.2*28.0*65.6	
Packing Dimension (W*D*H)	(mm)	690*790*1810	
	(in)	27.1*31.0*71.1	
Container Load (20'/40'/40'H)	22/44/44		

	7	
	£ .	
		II.
*		N.
	1	
		ll
	parks.	

Accessories			
Caster	Y	′	
Foot	Y	′	
Porthole	Y		
Shelves/Drawers	5/-		
USB Interface	Optional	Υ	
Temperature Recorder	Optional		

Alarms	
High/Low Temp	Υ
Remote System	Υ
Power Failure	Υ
Sensor Error	Υ
Low Battery	Υ
Door Ajar	Υ

Other			
Certificate	CE	UL	
EURO CHILL			

Scope of Application



MEDI©GUARD PREMIUM

This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at $2 \sim 8$ °C.

Specifications



MG400PR

Technical Data		
Cabinet Type	Upright, o	glass door
Climate Class	1	V
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
Noise((dB(A))	43	43





Performance	
Temp Range(°C)	2~8

Control	
Controller	Microprocessor
Display	LED

Electrical Data		
Power Supply(V/Hz)	220~240/50&60	115/60
Power(W)	360	240
Electrical Current(A)	1.8	2.4

Dimensions				
390/13.8				
(kg)	116/139			
(lbs)	255.8/306.6			
(mm)	530*555*1380			
(in)	20.9*21.9*54.3			
(mm)	665*710*1965			
(in)	26.2*28.0*77.4			
(mm)	690*790*2110			
(in)	27.1*31.0*82.9			
21/42/42				
	(lbs) (mm) (in) (mm) (in) (mm)			

Accessories		
Caster	Υ	
Foot	Υ	
Porthole	Υ	
Shelves/Drawers	7/-	
USB Interface	Optional	Υ
Temperature Recorder	Optional	
	Ориона	

Alarms	
High/Low Temp	Υ
Remote System	Υ
Power Failure	Υ
Sensor Error	Υ
Low Battery	Y
Door Ajar	Y

Other				
Certificate	CE	UL		
ROPECIALISTS IN MEDICAL REFRIGERATION AND EQUIPMENT				