



**“Your HotBag saved her life.**

Our crew began warming a hypothermia victim with a blanket from your HotBag after she fell through the ice. By the time they reached the hospital her core temperature had been raised from the low ~26.6°C to 35°C.”

Brett Murphy  
Regional EMS

# BLANKET WARMERS

DESIGNED FOR PRE-HOSPITAL TRAUMA AND MEDICAL CARE



## 3 Sizes To Choose From

**2-3**  
blankets

### SMALL

Order #  
MED-BLANKET-1LS

**4-5**  
blankets

### MEDIUM

Order #  
MED-BLANKET-XWLS

**6-7**  
blankets

### LARGE

Order #  
MED-BLANKET-2LS

All of our Blanket Warmers are

**IEC 60601-1**

Medical Accessory Compliant

- Interior Air Temperature of 54-60°C
- **Low Voltage** (12vDC) <50W
- Removable Heating Elements
- Flat Velcro Closure

- Durable Nylon Construction (Bag)
- **Washable** (Bag Only)
- Lightweight/Portable
- 12 Month Warranty

## Always ready

When plugged into an available power source The Electric HotBag™ keeps blankets at a safe, constant temperature, always ready for use. (Power supply and adaptor sold separately.)

## AC Power Supply

Order #  
HBPS-1



## 12v Cord Adaptor

Order #  
HBPC-AMP-72



## Orders, Quotes & Questions

A real person is waiting to speak with you. Contact us today!  
8am-4:30pm AEST, Mon-Fri  
[www.midmed.com.au](http://www.midmed.com.au)  
[sales@midmed.com.au](mailto:sales@midmed.com.au)

**1300 643 633**

Office Number



 **midmed**  
providing positive solutions