

MEDI GUARD 1001 PLUS


VACCINE REFRIGERATORS



BONUS

24 MONTHS PARTS

WHEN YOU REGISTER WARRANTY ONLINE





MEDI GUARD 401 PLUS

MEDI GUARD 601 PLUS

MEDI GUARD 1001 PLUS

- QCPP Compliant
- Operating Range 2-8°C
- Gross Volume 1000 LT
- Exterior Dimensions W1130 x D720 x H2023 mm
- Audible & Visual Hi/Lo Temp. Alarms
- Door Open Alarm
- Self-Closing Lockable Glass Doors
- 6 Fully Adjustable Shelves per door with Front & Rear Upstands
- Easy to Use Customised Digital Controller with High / Low Memory & Easy Reset Function
- Dual Temperature Probes With Simulated Product Temperature & Internal Air Temperature
- White Powder Coated Steel Internal & External
- 4 Castors, with Two Front Locking
- External Access Port
- 24 Months Warranty
- Energy Consumption 9.903kWh
- Bonus Medi Log 2 Product Simulation Data Logger Included With Audible & Visual Alarms
- Available in 3 sizes



MODEL	MEDI GUARD 401 PLUS	MEDI GUARD 601 PLUS	MEDI GUARD 1001 PLUS
Capacity Litre	400	600	1000
Ext Dimensions Mm	606W * 620D * 1950H	700W * 790D * 2058H	1130W * 720D * 2023H
Int Dimensions Mm	513W * 495D * 1501H	603W * 665D * 1596H	1033W * 595D * 1561H
Shelves	7	6	12
Power Consumption	4.935kWh	6.209kWh	9.903kWh