

# MEDI GUARD 401 PLUS

## VACCINE REFRIGERATORS



**BONUS**

**24 MONTHS PARTS**

WHEN YOU REGISTER WARRANTY ONLINE





**MEDI GUARD 401 PLUS**



**MEDI GUARD 601 PLUS**



- QCPP Compliant
- Independently Tested By The Pharmacy Guild Of NSW
- Operating Range 2-8°C
- Gross Volume 400 LT
- Exterior Dimensions W606 x D620 x H1950 mm
- Audible & Visual Hi/Lo Temp. Alarms
- Door Open Alarm
- Self-Closing Lockable Glass Door
- 7 Fully Adjustable Shelves with Front & Rear Upstands
- Easy to Use Customised Digital Controller with High / Low Memory & Easy Reset Function
- Dual Temperature Probes With Simulated Product Temperature & Internal Air Temperature
- White Powder Coated Steel Internal & External
- 4 Castors, with Two Front Locking
- External Access Port
- 24 Months Warranty
- Energy Consumption 4.935kWh
- Bonus Medi Log 2 Product Simulation Data Logger Included With Audible & Visual Alarms



MODEL	MEDI GUARD 401 PLUS	MEDI GUARD 601 PLUS	MEDI GUARD 1001 PLUS
Capacity Litre	400	600	1000
Ext Dimensions Mm	606W * 620D * 1950H	700W * 790D * 2058H	1130W * 720D * 2023H
Int Dimensions Mm	513W * 495D * 1501H	603W * 665D * 1596H	1033W * 595D * 1561H
Shelves	7	6	12
Power Consumption	4.935kWh	6.209kWh	9.903kWh