

Technology with style

engineered to meet the
most demanding
medical environments



Bariatric Cadaver Stretcher Lifter

Rago

First choice for professional refrigeration

Bariatric Cadaver Stretcher Lifters

Features

- ❖ 100% Stainless Steel
- ❖ 24V. Battery operated silent actuator
- ❖ Dual Push button panel rocker switch & remote handset
- ❖ IP rated for hosing down
- ❖ 4 heavy duty 200 mm swivel lockable rubber castors
- ❖ High Torque DC reduction drive motor - forward & reverse
- ❖ Sturdy polypropelene rollers
- ❖ Lock brace & guide rails for cadaver tray
- ❖ Corner bump rails
- ❖ 10A male receptacle for simple charging
- ❖ Fully welded 38 mm continuous hand rail
- ❖ Battery Charge Level Indicator
- ❖ Battery Charging Indicator
- ❖ Floor hugging drive wheel ideal for uneven surfaces
- ❖ Design rated for up to 1:7 slope upward and downward
- ❖ Australian Design & Manufactured
- ❖ Auto Stop Function
- ❖ Emergency Stop Button as standard
- ❖ Dimensions 2150 x 800
- ❖ Travel : 425 mm min. height
1725 mm max. height
- ❖ Capacity: 350 kg SWL



Note: Designed to be used with Rago Mortuary Refrigerators & Freezers models SM,SMF,SMB, SMBF 1, 2, & 3 berth



Rago Refrigeration Pty Ltd, primarily refrigeration engineers, has concentrated on the manufacturing of the highest quality refrigeration equipment for 40 years, incorporating the latest designs and concepts which is backed by qualified engineers and technical staff. Manufacture and installation focuses upon the individual needs of our clients.

Manufactured by:



Rago Refrigeration Pty Ltd

PO Box 3297

Wareemba NSW 2046 Australia

p 0411 399 636

e info@rago.com.au

w www.rago.com.au