

# Surgimove

## Surgeon chairs



### Benefits

SurgiMove is perfectly suited for a wide range of surgical applications. Its advanced design delivers multiple positions including 360 degree adjustable armrests that can be provided as optional. Pedals are user friendly, the backrest and the armrests are fully adjustable allowing to the SurgiMove to reach superior comfort and adjustability.

Additional benefits include:



- Effective central braking system on two wheels
- Excellent maneuvering thanks to the double rim wheels
- Generous space for legs thanks to the frame design
- Comfortable height adjustment via foot pedal

**2D armrest** can be fitted only with seat fixed angle.

**3D, 4D and revolving armrests** can be fitted only with the adjustable seating angle.



### Available models



**SH** Seat height

**WH** Working height

### Armrests



### Standard configuration, alternative, accessories

	515.7	525.7	531.7	558.7
Ergonomic seat in Stamskin	✓			
Curve seat in Stamskin		✓		
Saddle seat in Stamskin			✓	
Callista seat in Stamskin				✓
Stamskin backrest adjustable also in height	✓	✓	✓	✓
Fixed seating angle	✓	✓	✓	
Adjustable seating angle (not compatible with swivable armrest)	✓	✓	✓	✓
Aluminium base with centralized brake	✓	✓	✓	✓
Height adjustable armrest	✓	✓	✓	✓
4D armrest	✓	✓	✓	✓
Swivel and height adjustable armrest	✓	✓		
3D adjustable padded armrest	✓	✓		
Ring mechanism for height adjustment	✓	✓	✓	✓
Wheels:				
• soft castors for hard floor	✓	✓	✓	✓

✓ Standard

✓ Accessory





**Swivel and height adjustable armrest**



**Seat**



**Height adjustment ring mechanism**



**Wheels**



**3D adjustable padded armrest**



**Backrest**

OFFICIAL DISTRIBUTOR IN AUSTRALIA:



Unit 4/62 Borthwick Ave,  
Murarrie QLD 4172 (Australia)

Tel: 1300.643.633  
sales@midmed.com.au  
www.midmed.com.au