

SERENITY NEXT

- GYNECOLOGICAL EXAMINATIONS & DIAGNOSTIC
- HYSTEROSCOPY
- LASER TREATMENTS
- MEDICALLY ASSISTED PROCREATION
- OBSTETRICS EXAMINATIONS & MONITORING
- PROCTOLOGY
- SCREENING & COLPOSCOPY MONITORING
- ULTRASOUND DIAGNOSTIC
- UROLOGY

Reassuring design, ergonomics and functionality. Lifting capacity up to 250 kg. Easy treating and advanced features. Maximum comfort for both operator and patients.

Standard features

- Electric movement up, down, Trendelenburg, reverse Trendelenburg and backrest recline
- Paper roll holder 60 cm
- Pushbutton control
- Security stops
- Soft start and stop
- Stainless steel 2L basin
- Wrapping seamless upholsteries with antimicrobial feature (different colors available)



Options

- 2 auxiliary electric sockets IPX4 in the base (300W)
- Back-up battery (for brief movement in the absence of power supply)
- Colposcope support arm, with anti-crushing safety feature
- Doctor footrest
- Electrically controlled footrests with arc supports
- Electrically heated seat upholstery
- External wheels Ø 50mm
- Fireproof and antistatic fabric available
- Foot pedal with 3 user profiles, 1 preset program + 3 programmable positions for each profile
- Footrests with arc supports and 3-lever adjustment
- Göpel leg holders with 3 levers
- Göpel leg holders with 2 levers and vertical lowering
- Infusion transfusion armrest
- Instrument tray
- Integrated extractable leg-rest
- Magnetic headrest cushion
- Memory Foam upholstery
- Movable patient side handles
- Pedal control with display
- Quick Release system for emergency backrest manual lowering
- Round basin for urological visits with valve
- Slow mode for precision, ultra-slow patient positioning / available with foot pedal only
- Tool holder lateral bars



Auxiliary sockets



Colposcope-support-arm



Doctor footrest



Electrical footrest



Extractable leg-rest



Foot pedal



Footrest arc support with 3 levers



Göpel leg holders with 3 levers



Göpel leg holders with 2 levers and vertical adjustment



Infusion armrest



Instrument tray



Magnetic headrest cushion



Memory Foam upholstery



Movable patient handles



Paper roll



Pushbutton control



Quick Release



Stainless steel 2L basin



Tool holder lateral bars



Urology basin

SERENITY NEXT



BLOOD TEST CHAIR 9SP0047P

Standard features

- Independent backrest and legrest adjustment
- Movement served by selflocking gas springs, with lateral levers
- Anatomic paddings covered with washable fire-proof fabric
- Adjustable armrest with friction clutch
- Collapsible clamps for IV pole at the 2 sides
- Frame in steel tube painted with epoxypolyester powders
- 4 wheels with brake – 75 mm diameter



Options



9SP9115 (OPTION)
Headrest pillow



Stainless or plastic IV pole
with 2 hooks
(OPTION)



9SP9110 (OPTION)
Paper roll holder



Stainless or plastic IV pole
with 4 hooks
(OPTION)