

Ameda® Mya™ Joy Double Breast Pump

with Large Tote

131A41



Trusted by Hospitals, Loved by Moms.

The Mya Joy has been designed with versatile features, safety and comfort, offering the best value for moms who are looking for a convenient pumping experience. Compact, lightweight and quiet, the Mya Joy also features the world's only FDA-cleared pump kit that helps protect breastmilk from contaminants. This pump delivers on the ultimate promise of comfort by offering a wide range of settings, mom-friendly features, and flange sizes. Whether pumping on-the-go, at home or at work, Mya Joy offers pumping without compromise.



Safe

Ameda's HygieniKit true closed system not only helps protect breastmilk from backing up into the tubing but also prevents the system from potential contaminants like bacteria, viruses and mold growth



Comfort

wide range of suction and speed settings, functions, and flange sizes



Ultra-Quiet

operates at less than 45dB and allows for discreet pumping anywhere



Powerful

Hospital Strength Performance and two modes of pumping



Portable

battery operated, lightweight and clips right onto a waist band, purse pocket or favorite yoga pants

MyaJoy Features:

- Hospital strength performance
- Up to 250 mmHg suction strength
- Backlit digital LCD display
- Timer display easily tracks pumping time
- Smart auto power off shuts down after 45 minutes
- 2-phase adjustable suction modes:
- Stimulation (6 Levels) & Expression Mode (12 Levels)
- Double or single electric pumping options
- 2-year limited warranty with Ameda ParentCare Support

Contents Include:

- Ameda Mya Joy Breast Pump
- AC Power Adapter
- Dual HygieniKit® Milk Collection System
- CustomFit™ Flanges (25mm, 28.5mm)
- 2 Extra Valves (for HygieniKit)
- Tote Bag
- 2 Extra Diaphragms (for HygieniKit)
- Quick Start Guide



SKU 131A41 UPC 810725025079