

# PROMOTAL SWIFI

Swifi® is a transfer chair with an ultra modern and welcoming design that contributes to reducing patient anxiety on arrival at a hospital. The patient's privacy is respected when manipulating the footrest due to its bilateral control. The chair allows for easy access to all patients, notably reduced mobility or corpulent patients.



**MADE IN FRANCE**

## OVERVIEW

- A** Epoxy finished push bar
- B** Upholstery widths of 47cm or 60cm
- C** Patient weight capacity of up to 200kg
- D** Integrated clothes basket
- E** Retractable footrest with bilateral control
- F** Rear caster diameter of 300mm with centralised braking system
- G** Oxygen bottle holder
- H** Adjustable arm rests with stoppers



Swifi® was designed exclusively for a hospital environment: greater rigidity and a longer life span, easy manoeuvrability and a large range of accessories including Swing IV®, an exclusive system of rotating IV pole.

## FEATURES



### Accessibility

Liftable armrests, Open-sided seat, Retractable footrest with bilateral control



### Hygiene

Seamless seat and back rests and a mono block metallic frame contribute to making cleaning easier.



### Safety

Up to 200kg patient weight capacity, Rear casters (300mm) with centralised braking system, Front shock absorbers



### Comfort

Built in clothes basket, Patient file holder, Full width push bar, 300mm rear casters, various accessories. Available in 2 widths to adapt to different patient body shapes.