

WWW.MIDMED.COM.AU



# SEPSIS TROLLEY

## Product Images



## Description

- Trolleys designed for use in the treatment of Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- Two locking options
- High Security Bolt Lock - Supplied with two keys
- Electronic Push Button - Programmable four-digit code
- Two adjustable door shelves & an A4 document holder
- Six labelled drawers (Oxygen/Blood Cultures/Intravenous Antibiotics/Intravenous Fluid/Haemoglobin & Lactates/Urine Measurement)
- Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting Plate' (fits to the l/h side of the trolley)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- 125mm polymer castors
- When not in use, trolleys storing medicine should be secured to the wall. See option MT/1/6F Security Clamp



## Overall Dimensions (W x D x H) & Colours

- 511 x 542 x 1072mm
- Drawers - 335 x 295 x 113mm (Internal)
- Paint - Grey White
- Work Surface - Grey White

Code	Description	Overall Dimensions (W x D x H)	Weight
BM.SEPS/6/HSB	Sepsis Trolley - Six Drawer, High Security Bolt Lock	511 x 542 x 1072mm	53.4kg
BM.SEPS/6/EPB	Sepsis Trolley - Six Drawer, Electronic Push Button Lock	511 x 542 x 1072mm	53.4kg

