

## Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column

### Features

- Standard level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
- Top Panel - Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
- Side Panels - Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions
- Single Column - 550 x 550 x 993mm
- Double Column - 1005 x 550 x 993mm
- Plastic Top
- Single Column - 455 x 455 (internal) at height 968mm
- Double Column - 915 x 455 (internal) at height 968mm
- Trays
- Shallow Tray - 355 x 405 x 100mm (internal)
- Deep Tray - 355 x 405 x 150mm (internal)

Item	Description	Weight
BM.CTM1/10/1S3D	1 x 100mm, 3 x 150mm	27.5kg
BM.CTM1/10/3S2D	Trays 3 x 100mm, 2 x 150mm	28.6kg
BM.CTM1/10/4D	150mm Trays 4 x 150mm	27.9kg
BM.CTM1/10/4S1D	Trays 4 x 100mm, 1 x 150mm	28.3kg
BM.CTM1/10/6S	Trays 6 x 100mm	29.4kg



Easy disassembly & cleaning



BM.CTM1/10/3S2D