

# ENTERAL FEEDING



[www.midmed.com.au](http://www.midmed.com.au)



1300 643 633

## WHY WAS ENFit™ CREATED?

In 2006, following several incidences of tubing misconnections within the clinical setting, the International Organisation, GEDSA was formed. After years of research and development, a unique locking system was developed called ENFit™. The recommendations for the system was ISO:80369 and this has now evolved to the Standard ISO:80369-3.


GEDSA is a non-profit trade association, formed to establish a voice for addressing issues that faced enteral device manufacturers, suppliers and distributors, and to help international standards for healthcare tubing connectors comprised of manufacturers, clinical and medical staff and suppliers worldwide, GEDSA facilitates information flow about the initiative, which is designed to increase patient safety and optimal delivery of enteral feeding by reducing the risk of tubing misconnections. Our staff are GEDSA trained to provide you with information in accordance with ISO:80369-3



## SUMMARY

Together with GEDSA, and our trusted suppliers, Midmed is pleased to offer our range of Enteral Feeding solutions to ensure the risk of tubing misconnections are minimised.

# ENFit SINGLE AND HOME USE SYRINGES



## STERILE ENFit™ HOME-USE SYRINGES - Up to 14 Day Use

- MS.SE01H - 1ml | 100 Per Box - Low Dose Syringe
- MS.SE03H - 3ml | 100 Per Box - Low Dose Syringe
- MS.SE05H - 5ml | 100 Per Box - Low Dose Syringe
- MS.SE10H - 10ml | 100 Per Box
- MS.SE20H - 20ml | 50 Per Box
- MS.SE60H - 60ml | 30 Per Box
- ISO 80369-3 standard tip for enteral use connection
- Silicone O-ring gasket design for clean and multiple use
- Available for Low Dose Tip design on 1ml, 3ml, 5ml
- Purple plunger for ENFit™ enteral identification
- Permanent and heat-etched dual graduations with more safety
- Safe plunger backstop to prevent withdrawal
- Sterilised by EO or Gamma, non-toxic, non-pyrogenic



## STERILE ENFit™ SINGLE-USE SYRINGES

- MS.SE01 - 1ml | 100 Per Box - Low Dose Syringe
- MS.SE03 - 3ml | 100 Per Box - Low Dose Syringe
- MS.SE05 - 5ml | 100 Per Box - Low Dose Syringe
- MS.SE10 - 10ml | 100 Per Box
- MS.SE20 - 20ml | 50 Per Box
- MS.SE60 - 60ml | 30 Per Box
- ISO 80369-3 standard tip for enteral use connection
- Available for Low Dose Tip design on 1ml, 3ml and 5ml
- Purple plunger for ENFit™ enteral identification
- Permanent and heat-etched dual graduations with more safety
- Safe plunger backstop to prevent withdrawal
- Sterilised by EO or Gamma, non-toxic, non-pyrogenic



# SHORT TERM FEEDING TUBES



## Product Features:

- Range of Fr sizes (4Fr - 14 Fr)
- Range of lengths (50cm - 92cm)
- 2 non-opposing outlet holes
- Fully radiopaque
- Made of PVC
- Use for up to 7-10 days and 7 days for neonates
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA Free
- 25 units per box

4Fr



5Fr



6Fr



8Fr



10Fr



12Fr



14Fr




### STOCK CODE:

ENT.CF01.04.050ISOSAF	Nasogastric Feeding Tube 4Fr x 50cm
ENT.CF01.04.075ISOSAF	Nasogastric Feeding Tube 4Fr x 75cm
ENT.CF01.05.050ISOSAF	Nasogastric Feeding Tube 5Fr x 50cm
ENT.CF01.05.075ISOSAF	Nasogastric Feeding Tube 5Fr x 75cm
ENT.CF01.06.050ISOSAF	Nasogastric Feeding Tube 6Fr x 50cm
ENT.CF01.06.075ISOSAF	Nasogastric Feeding Tube 6Fr x 75cm
ENT.CF01.08.050ISOSAF	Nasogastric Feeding Tube 8Fr x 50cm
ENT.CF01.08.075ISOSAF	Nasogastric Feeding Tube 8Fr x 75cm
ENT.CF01.08.092ISOSAF	Nasogastric Feeding Tube 8Fr x 92cm
ENT.CF01.10.075ISOSAF	Nasogastric Feeding Tube 10Fr x 75cm
ENT.CF01.10.092ISOSAF	Nasogastric Feeding Tube 10Fr x 92cm
ENT.CF01.12.075ISOSAF	Nasogastric Feeding Tube 12Fr x 75cm
ENT.CF01.12.092ISOSAF	Nasogastric Feeding Tube 12Fr x 92cm
ENT.CF01.14.092ISOSAF	Nasogastric Feeding Tube 14Fr x 92cm



# FEEDING TUBES: PAEDIATRIC



## Product Features

- Transparent tube with 2 radiopaque stripes
- 29-day use
- Made of polyurethane
- Clearly numbered in centimetres
- Closed rounded tip
- 2 lateral distal openings
- ENFit™ male connector with tethered cap
- 40 units per box




STOCK CODE	CH/FR	LENGTH (CM)
CA.NCE04PU	4	40
CA.NCE104PU	4	125
CA.NCE05PU	5	50
CA.NCE105PU	5	125
CA.NCE06PU	6	50
CA.NCE106PU	6	120
CA.NCE306PU	6	110
CA.NCE08PU	8	50
CA.NCE108PU	8	125
CA.NCE308PU	8	110
CA.NCE10PU	10	50
CA.NCE110PU	10	125
CA.NCE310PU	10	110
CA.NCE312PU	12	110

# LONG TERM NUTRICARE INFANT FEEDING TUBES



## Product Features

- Range of Fr sizes (4Fr - 10 Fr)
- Range of lengths (50cm - 125cm)
- Full radiopaque
- Made of polyurethane
- Use up to 90 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- 25 units per box



### STOCK CODE:

ENT.NCI04.050ISOSAF	Nasogastric Feeding Tube 4Fr x 50cm
ENT.NCI04.075ISOSAF	Nasogastric Feeding Tube 4Fr x 75cm
ENT.NCI04.120ISOSAF	Nasogastric Feeding Tube 4Fr x 120cm
ENT.NCI05.050ISOSAF	Nasogastric Feeding Tube 5Fr x 50cm
ENT.NCI05.075ISOSAF	Nasogastric Feeding Tube 5Fr x 75cm
ENT.NCI05.125ISOSAF	Nasogastric Feeding Tube 5Fr x 125cm
ENT.NCI06.050ISOSAF	Nasogastric Feeding Tube 6Fr x 50cm
ENT.NCI06.075ISOSAF	Nasogastric Feeding Tube 6Fr x 75cm
ENT.NCI06.125ISOSAF	Nasogastric Feeding Tube 6Fr x 125cm
ENT.NCI08.050ISOSAF	Nasogastric Feeding Tube 8Fr x 50cm
ENT.NCI08.075ISOSAF	Nasogastric Feeding Tube 8Fr x 75cm
ENT.NCI08.120ISOSAF	Nasogastric Feeding Tube 8Fr x 120cm
ENT.NCI10.075ISOSAF	Nasogastric Feeding Tube 10Fr x 75cm
ENT.NCI10.120ISOSAF	Nasogastric Feeding Tube 10Fr x 120cm



# LONG TERM FEEDING TUBES


## Product Features

- Range of Fr sizes (6Fr - 14Fr)
- Range of lengths (50cm - 120cm) Fully radiopaque
- Made of polyurethane
- Use up to 90 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- 5 units per box





### STOCK CODE:

ENT.NC06.050ISOAF-ROW	Nasogastric Feeding Tube 6Fr x 50cm
ENT.NC06.080ISOAF-ROW	Nasogastric Feeding Tube 6Fr x 80cm
ENT.NC06.110ISOAF-ROW	Nasogastric Feeding Tube 6Fr x 110cm
ENT.NC08.050ISOAF-ROW	Nasogastric Feeding Tube 8Fr x 50cm
ENT.NC08.080ISOAF-ROW	Nasogastric Feeding Tube 8Fr x 80cm
ENT.NC08.092ISOAF-ROW	Nasogastric Feeding Tube 8Fr x 92cm
ENT.NC08.110ISOAF-ROW	Nasogastric Feeding Tube 8Fr x 110cm
ENT.NC08.120ISOAF-ROW	Nasogastric Feeding Tube 8Fr x 120cm
ENT.NC10.092ISOAF-ROW	Nasogastric Feeding Tube 10Fr x 92cm
ENT.NC10.110ISOAF-ROW	Nasogastric Feeding Tube 10Fr x 110cm
ENT.NC12.092ISOAF-ROW	Nasogastric Feeding Tube 12Fr x 92cm
ENT.NC12.110ISOAF-ROW	Nasogastric Feeding Tube 12Fr x 110cm
ENT.NC12.120ISOAF-ROW	Nasogastric Feeding Tube 12Fr x 120cm
ENT.NC14.092ISOAF-ROW	Nasogastric Feeding Tube 14Fr x 92cm
ENT.NC14.092ISOAF-ROW-GW	Nasogastric Feeding Tube 14Fr x 92cm
ENT.NC14.110ISOAF-ROW	Nasogastric Feeding Tube 14Fr x 110cm
ENT.NC14.120ISOAF-ROW	Nasogastric Feeding Tube 14Fr x 120cm
ENT.NC16.092ISOAF-ROW	Nasogastric Feeding Tube 16Fr x 92cm
ENT.NC16.092ISOAF-ROW-GW	Nasogastric Feeding Tube 16Fr x 92cm
ENT.NC16.120ISOAF-RW	Nasogastric Feeding Tube 16Fr x 120cm
ENT.NC08.140ISOAF	Nasogastric Feeding Tube 8Fr x 140cm
ENT.NC10.140ISOAF	Nasogastric Feeding Tube 10Fr x 140cm



## NUTRICARE FLO FEEDING TUBES - INFANT/PAEDIATRIC POLYURETHANE OPEN END USE



### Product Features

- Range of Fr sizes (4 - 12Fr)
- Range of lengths (40cm - 125cm)
- Radiopaque stripe
- Made of polyurethane
- Use up to 30 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- Polyurethane feeding tube with rounded open distal end
- The open end eliminates the collection of feed in the tip of the tube
- Less chance of bacterial colonisation
- 25 units per box

#### STOCK CODE:

ENT.NCF04.40	Nasogastric Feeding Tube 4Fr x 40cm
ENT.NCF04.50	Nasogastric Feeding Tube 4Fr x 50cm
ENT.NCF04.75	Nasogastric Feeding Tube 4Fr x 75cm
ENT.NCF05.40	Nasogastric Feeding Tube 5Fr x 40cm
ENT.NCF05.50	Nasogastric Feeding Tube 5Fr x 50cm
ENT.NC05.60	Nasogastric Feeding Tube 5Fr x 60cm
ENT.NC05.75	Nasogastric Feeding Tube 5Fr x 75cm
ENT.NC05.90	Nasogastric Feeding Tube 5Fr x 90cm
ENT.NC05.125	Nasogastric Feeding Tube 5Fr x 125cm
ENT.NC06.40	Nasogastric Feeding Tube 6Fr x 40cm
ENT.NC06.50	Nasogastric Feeding Tube 6Fr x 50cm
ENT.NC06.60	Nasogastric Feeding Tube 6Fr x 60cm

#### STOCK CODE:

ENT.NC06.75	Nasogastric Feeding Tube 6Fr x 75cm
ENT.NC06.90	Nasogastric Feeding Tube 6Fr x 90cm
ENT.NC06.125	Nasogastric Feeding Tube 6Fr x 125cm
ENT.NC08.40	Nasogastric Feeding Tube 8Fr x 40cm
ENT.NC08.50	Nasogastric Feeding Tube 8Fr x 50cm
ENT.NC08.60	Nasogastric Feeding Tube 8Fr x 60cm
ENT.NC08.75	Nasogastric Feeding Tube 8Fr x 75cm
ENT.NC08.90	Nasogastric Feeding Tube 8Fr x 90cm
ENT.NC08.125	Nasogastric Feeding Tube 8Fr x 125cm
ENT.NC10.75	Nasogastric Feeding Tube 10Fr x 75cm
ENT.NCF10.110	Nasogastric Feeding Tube 10Fr x 110cm
ENT.NCF12.110	Nasogastric Feeding Tube 12Fr x 110cm



# FEEDING TUBES: ADULT LEVIN TUBES - POLYURETHANE



Enteral Naso-Gastric Tubes, Levin Type Duration of use:  
Up to 29 days



## Product Features

- Transparent tube with 2 radiopaque stripes
- Clearly numbered in centimetres
- Closed rounded tip
- ENFit™ male connector with tethered cap
- 20 units per box
- 4 lateral distal openings, allows for drainage and feeding

Enteral Naso-Gastric Tubes, Levin Type with Guide-Wire  
Duration of Use: Up to 29 days



STOCK CODE	CH/FR	LENGTH (CM)
CA.NCE106AU	6	110
CA.NCE506AU	6	130
CA.NCE108AU	8	110
CA.NCE508AU	8	130
CA.NCE110AU	10	110
CA.NCE510AU	10	130
CA.NCE112AU	12	110
CA.NCE512AU	12	130
CA.NCE114AU	14	110
CA.NCE514AU	14	130
CA.NCE116AU	16	110
CA.NCE516AU	16	130
CA.NCE118AU	18	110
CA.NCE518AU	18	130
CA.NCE120AU	20	110
CA.NCE520AU	20	130

STOCK CODE	CH/FR	LENGTH (CM)
CA.NCE106G	6	110
CA.NCE506G	6	130
CA.NCE108G	8	110
CA.NCE508G	8	130
CA.NCE110G	10	110
CA.NCE510G	10	130
CA.NCE112G	12	110
CA.NCE512G	12	130
CA.NCE114G	14	110
CA.NCE514G	14	130
CA.NCE116G	16	110
CA.NCE516G	16	130

## FEEDING TUBES: ADULTS DUODENAL/JEJUNAL - POLYURETHANE



**Enteral Naso-Duodenal Tubes  
with Weight and Guide-Wire  
Duration of use: Up to 29 Days  
20 Units Per Box**

STOCK CODE	CH/FR	LENGTH (CM)	WEIGHT
CA.NCE206G	6	110	1.1g
CA.NCE208G	8	110	1.1g
CA.NCE210G	10	110	1.7g
CA.NCE212G	12	110	2.8g
CA.NCE214G	14	110	4.0g

**Enteral Naso-Jejunal Tubes with  
Weight and Guide-Wire  
Duration of use: Up to 29 Days  
20 Units Per Box**

STOCK CODE	CH/FR	LENGTH (CM)	WEIGHT
CA.NCE306G	6	130	2.8g
CA.NCE308G	8	130	2.8g
CA.NCE310G	10	130	4.0g
CA.NCE312G	12	130	4.0g
CA.NCE314G	14	130	4.0g

**Enteral Naso-Duodenal Tubes  
with Weight  
Duration of use: Up to 29 Days  
20 Units Per Box**


STOCK CODE	CH/FR	LENGTH (CM)	WEIGHT
CA.NCE206AU	6	110	1.1g
CA.NCE208AU	8	110	1.1g
CA.NCE210AU	10	110	1.7g
CA.NCE212AU	12	110	2.8g
CA.NCE214AU	14	110	4.0g





**ENFit™ SAMPLING STRAW (5cm & 15cm)**

Stock Code: CA.NCE15S - Wide Bore  
 Stock Code: CA.NCE05S - Fine Bore  
 Single Use  
 150 Per Box



**ENFIT™ FILTER NEEDLE CAP**

Stock Code: ENT.FNISO-50  
 Single Use 5µm 18G x 50mm  
 30 Per Box



**SOFT WIDE FUNNEL WITH CAP**

Stock Code: CA.NCE101A  
 100 Per Box



**ENFIT™ STEPPED CAP**

Stock Code: ENT.SBAISO  
 Stepped Cap for Oral/ENFit Syringes  
 50 Per Box



**ENFIT™ STEPPED CAP**

Stock Code: CA.NCE100C  
 10mm to 24mm  
 100 Per Box



**ENFit™ SCREW CAP**

Stock Code: CA.NCE302C  
 28mm  
 50 Per Box



**SOFT FUNNEL WITH CAP**

Stock Code: CA.NCE100A  
 150 Per Box




**4 - STEPPED CONNECTOR**

Stock Code: CA.NCE104A  
 100 Per Box



**5 - STEPPED CONNECTOR**

Stock Code: CA.NCE105A  
 100 Per Box



**Internal bottle caps - Allows for tamper proof cap to be replaced**

**Stock Code: ENT.BAISO1-S**

Single Use Size 1-12-15mm  
100 Per Box

**Stock Code: ENT.BAISO2-S**

Single Use Size 2-14-16.5mm  
100 Per Box

**Stock Code: ENT.BAISO3-S**

Single Use Size 3-17-20.5mm  
100 Per Box

**Stock Code: ENT.BAISO35-S**

Single Use Size 4-19-23mm  
100 Per Box

**Stock Code: ENT.BAISO4-S**

Single Use Size 5-24-26mm  
100 Per Box



**Y-PORT CONNECTOR  
(24 HOUR USE)**

**Stock Code: CA.NCE100Y**

Y-Port - 2 Male, 1 Female  
100 Per Box



**ENFit™ 3-WAY STOPCOCK  
(24 HOUR USE)**

**Stock Code: CA.NCE300R**

3-Way 2 - Male, 1 Female  
100 Per Box



**ACE CONNECTOR  
(30 DAY USE)**

**Stock Code: DA.485 - Box of 50**

**Stock Code: DA.485-10 Bag of 10**

Ace Connector - Closed Connection System  
30 day use - Dual Access Port



**DRAINAGE BAG (50ml)**

**Stock Code: ML.DBISO-50**

50ml Drainage Bag  
EVA Plastic - Non PVC  
20 Per Box



**DRAINAGE BAG (500ml)  
& (1000ml)**

**Stock Code: MIDMED-500**

500ml Drainage Bag  
EVA Plastic - 100cm tubing - Non PVC  
20 Per Box

**Stock Code: MIDMED-1000**

1000ml Drainage Bag  
EVA Plastic - 100cm tubing - Non PVC  
20 Per Box



**ENTERAL EXTENSION**

**Stock Code: ENT.ESCF.075-ISOSAF**

75cm  
25 Per Box  
Priming Volume 1.15ml

**Stock Code:**

**CA.NCE1315E**  
2.8mm Ø | 150cm  
25 Per Box

**Stock Code:**

**CA.NCE1215E**  
1.6mm Ø | 150cm  
50 Per Box



**PH STRIPS**

Stock Code: MN.92118



**SINGLE USE CLEANING BRUSH**

Stock Code: ENCS2017100



**MULTI USE CLEANING BRUSH**

Stock Code: ENC125212

## Enteral Syringe Delivery Driver

### Product Features

- Can be calibrated to meet your specific ENFit™ syringe technical requirements
- Bolus function
- Various visual and audible alarms
- Only compatible with Enteral syringe
- Rate or time delivery modes
- Battery running time - more than 6h at a rate of 5ml/h



Pre-calibrated to manage the following ENFit™ Syringes:

- MS.SE05
- MS.SE10
- MS.SE20
- MS.SE60

# TUMMY BUTTONS




**A button can be a comfortable solution, because of its small size and the fact that you don't have to carry a tube. We are committed to listen to the end-users and specialists and we show this by innovating our products year in, year out.**

**The MediLime Tummy Button has one sole purpose, and that is to make living with a feeding system more pleasant and efficient.**

## Product Features

- The MediLime Tummy Button consists of pure silicone, is latex- and DEHP free, resulting in a healthy skin around the gastro stoma.
- An apple shaped balloon improves an ideal shape for closing the stoma channel.
- Best compatibility with all MediLime Extension Tubes.
- Various balloon diameters (per size 12FR up to 24FR, from 3 up to 8ml volume) and is available in all lengths (cm) (0.8, 1.0, 1.2, 1.5, 1.7, 2.0, 2.3, 2.5, 2.7, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5 and 6.0).
- The Tummy Button comes with an input guide. This can be a useful tool when placing the button if there is scar tissue at the fistula opening or with a larger size of the Tummy Button.



## Package Contents

The Tummy Button is packaged in a box along with the following accessories for easy placement:

- Input Guide
- Syringe
- Lubricating gel 5 grams
- Gauze



# TUMMY BUTTONS



## STOCK CODE:

ML.0320-12-0.8-02	12Fr	0.8cm
ML.0320-12-1.0-02	12Fr	1.0cm
ML.0320-12-1.2-02	12Fr	1.2cm
ML.0320-12-1.5-02	12Fr	1.5cm
ML.0320-12-1.7-02	12Fr	1.7cm
ML.0320-12-2.0-02	12Fr	2.0cm
ML.0320-12-2.3-02	12Fr	2.3cm
ML.0320-12-2.5-02	12Fr	2.5cm
ML.0320-12-2.7-02	12Fr	2.7cm
ML.0320-12-3.0-02	12Fr	3.0cm
ML.0320-12-3.5-02	12Fr	3.5cm
ML.0320-12-4.0-02	12Fr	4.0cm

## STOCK CODE:

ML.0320-18-1.0-02	18Fr	1.0cm
ML.0320-18-1.2-02	18Fr	1.2cm
ML.0320-18-1.5-02	18Fr	1.5cm
ML.0320-18-1.7-02	18Fr	1.7cm
ML.0320-18-2.0-02	18Fr	2.0cm
ML.0320-18-2.3-02	18Fr	2.3cm
ML.0320-18-2.5-02	18Fr	2.5cm
ML.0320-18-2.7-02	18Fr	2.7cm
ML.0320-18-3.0-02	18Fr	3.0cm
ML.0320-18-3.5-02	18Fr	3.5cm
ML.0320-18-4.0-02	18Fr	4.0cm
ML.0320-18-4.5-02	18Fr	4.5cm
ML.0320-18-5.0-02	18Fr	5.0cm

## STOCK CODE:

ML.0320-14-0.8-02	14Fr	0.8cm
ML.0320-14-1.0-02	14Fr	1.0cm
ML.0320-14-1.2-02	14Fr	1.2cm
ML.0320-14-1.5-02	14Fr	1.5cm
ML.0320-14-1.7-02	14Fr	1.7cm
ML.0320-14-2.0-02	14Fr	2.0cm
ML.0320-14-2.3-02	14Fr	2.3cm
ML.0320-14-2.5-02	14Fr	2.5cm
ML.0320-14-2.7-02	14Fr	2.7cm
ML.0320-14-3.0-02	14Fr	3.0cm
ML.0320-14-3.5-02	14Fr	3.5cm
ML.0320-14-4.0-02	14Fr	4.0cm
ML.0320-14-4.5-02	14Fr	4.5cm
ML.0320-14-5.0-02	14Fr	5.0cm
ML.0320-14-5.5-02	14Fr	5.5cm
ML.0320-14-6.0-02	14Fr	6.0cm

## STOCK CODE:

ML.0320-20-0.8-02	20Fr	0.8cm
ML.0320-20-1.0-02	20Fr	1.0cm
ML.0320-20-1.2-02	20Fr	1.2cm
ML.0320-20-1.5-02	20Fr	1.5cm
ML.0320-20-1.7-02	20Fr	1.7cm
ML.0320-20-2.0-02	20Fr	2.0cm
ML.0320-20-2.3-02	20Fr	2.3cm
ML.0320-20-2.5-02	20Fr	2.5cm
ML.0320-20-2.7-02	20Fr	2.7cm
ML.0320-20-3.0-02	20Fr	3.0cm
ML.0320-20-3.5-02	20Fr	3.5cm
ML.0320-20-4.0-02	20Fr	4.0cm
ML.0320-20-4.5-02	20Fr	4.5cm
ML.0320-20-5.0-02	20Fr	5.0cm

## STOCK CODE:

ML.0320-16-0.8-02	16Fr	0.8cm
ML.0320-16-1.0-02	16Fr	1.0cm
ML.0320-16-1.2-02	16Fr	1.2cm
ML.0320-16-1.5-02	16Fr	1.5cm
ML.0320-16-1.7-02	16Fr	1.7cm
ML.0320-16-2.0-02	16Fr	2.0cm
ML.0320-16-2.3-02	16Fr	2.3cm
ML.0320-16-2.5-02	16Fr	2.5cm
ML.0320-16-2.7-02	16Fr	2.7cm
ML.0320-16-3.0-02	16Fr	3.0cm
ML.0320-16-3.5-02	16Fr	3.5cm
ML.0320-16-4.0-02	16Fr	4.0cm
ML.0320-16-4.5-02	16Fr	4.5cm
ML.0320-16-5.0-02	16Fr	5.0cm

## STOCK CODE:

ML.0320-24-1.5-02	24Fr	1.5cm
ML.0320-24-1.7-02	24Fr	1.7cm
ML.0320-24-2.0-02	24Fr	2.0cm
ML.0320-24-2.3-02	24Fr	2.3cm
ML.0320-24-2.5-02	24Fr	2.5cm
ML.0320-24-2.7-02	24Fr	2.7cm
ML.0320-24-3.0-02	24Fr	3.0cm
ML.0320-24-3.5-02	24Fr	3.5cm
ML.0320-24-4.0-02	24Fr	4.0cm
ML.0320-24-4.5-02	24Fr	4.5cm
ML.0320-24-5.0-02	24Fr	5.0cm

## STOMA MEASURING DEVICE

The MediLime Stoma Measuring Device is used to obtain the correct size of the MediLime Tummy Button. The stoma canal is measured with this catheter. The use of the Stoma Measuring Device will result in a more precise Tummy Button placement. We strongly recommend using one for each (new) placement of a button in order to optimize the perfect fit.



## Product Features

- Silicone Material
- PHT and BPA free
- Latex Free
- Disposable
- Sterile

### STOCK CODE

ML.0325-00

MediLime Measuring Stoma Device with Syringe





# GASTRIC FEEDING TUBES



The MediLime Gastric Measuring Device with ENFit™ connection is made for administering enteral nutrition and medication. The Gastric Feeding Tube is made of pure silicone and is therefore skin-friendly and pliable. The size of the tube and the balloon capacity are clearly stated on the shaft.

STOCK CODE:				
ML.0503-12	Gastric Feeding Tube	1 - Way	12 Fr	2 Pieces
ML.0503-14	Gastric Feeding Tube	1 - Way	14 Fr	2 Pieces
ML.0503-16	Gastric Feeding Tube	1 - Way	16 Fr	2 Pieces
ML.0503-18	Gastric Feeding Tube	1 - Way	18 Fr	2 Pieces
ML.0503-20	Gastric Feeding Tube	1 - Way	20 Fr	2 Pieces
ML.0503-22	Gastric Feeding Tube	1 - Way	22 Fr	2 Pieces
ML.0503-24	Gastric Feeding Tube	1 - Way	24 Fr	2 Pieces



# JEJUNAL FEEDING TUBES & STOMA PLUGS



The MediLime Jejunal Feeding Tube is a 100% silicone device for delivery of enteral nutrition into the distal duodenum or jejunum. The device has a double lumens tube, an air-vented skin disk and a fixation balloon.

<b>STOCK CODE:</b>			
ML.103080031280A	1 - Way	12 Fr	5ml
ML.103080031480A	1 - Way	14 Fr	5ml
ML.103080031680A	1 - Way	16 Fr	10ml




## Product Features

- Stops stoma sites from closing
- Saves the stoma site until a new gastrostomy tube can be placed
- Includes a range of French gauge sizes in each pack which provides options even when stoma has started to close
- Effective, safe and inexpensive
- Saves staff in emergency situations
- Promptly dealing with gastrostomy emergencies prevents use of other incorrect or inappropriate feeding devices
- Made from high grade medical silicone for flexibility
- Tapered tip for ease of insertion



ENPLUGs are designed for emergency use to keep stomas from closing when a gastrostomy device has fallen out. The ENPLUG packs are available in 2 different lengths and various French sizes in case the stoma has already started to close and the currently used device no longer fits the stoma.

ENT.EN-PLUGS-S	10Fr / 12Fr / 14Fr / 16Fr x 4cm (Small)
ENT.EN-PLUGS-L	10Fr / 12Fr / 14Fr / 16Fr x 7cm (Large)
ENT.EN-PLUGS-XL	18 Fr / 20Fr x 7cm (Xtra Large)



# ORAL SYRINGES



## Oral Syringes

- Special tip design will not accept Luer connectors, preventing medication errors and accidental intravenous injections of oral liquids
- Permanent and heat-etched dual graduations with more safety
- Safe plunger backstop to prevent plunger withdrawal
- Not manufactured with PVC, BPA, Latex or DEHP
- Sterilised by EO or Gamma, non-toxic, non-pyrogenic
- Certificate: FDA, CE, MDSAP, MHRA, TGA, ISO13485, ISO9001



## Sterile Syringes

- MS.SO01.LDT.S - 1ml | 100 Per Box - Low Dose
- MS.SO.03 - 3ml | 100 Per Box
- MS.SO.05 - 5ml | 100 Per Box
- MS.SO.10 - 10ml | 100 Per Box
- MS.SO.20 - 20ml | 50 Per Box

## Non-Sterile Syringes

- MS.SO.03.ORAL - 3ml | 100 Per Box - Non Sterile
- MS.SO.05.ORAL - 5ml | 100 Per Box - Non-Sterile
- MS.SO.10.ORAL - 10ml | 100 Per Box - Non-Sterile
- MS.SO.20.ORAL - 20ml | 50 Per Box - Non-Sterile
- MS.SO.60.ORAL - 60ml | 40 Per Box - Non-Sterile





[www.midmed.com.au](http://www.midmed.com.au)



1300 643 633