

# Ameda® Pearl™ Hospital Grade Breast Pump



## Powerful Performance in the Hospital and Home.

Combines powerful performance inside with a protective shell on the outside. Pearl is designed to build milk supply amongst pump dependent, late pre-term and term moms. Pearl's modern design combines the durability needed for the multi-user/hospital environment with lifestyle features moms expect. Pearl will serve the needs of hospitals and is the ideal multi-user product for milk banks and rental stations. Pearl is perfectly positioned to support continuity of care in the Hospital and Home.


## Ameda Pearl – Proven Hospital Grade Pumping

*Clinical study in process.*

- **Proven:** Over 90% of moms surveyed were highly satisfied with Pearl
- **Powerful:** Hospital Grade / Multi-User Performance with two modes of pumping
- **Portable:** Lightweight with rechargeable battery and integrated handle
- **Smart:** Vacusense™ Technology ensures consistent vacuum performance in varying altitudes or changing atmospheric conditions
- **Comfortable:** Fully customizable expression mode with 72 possible combinations of speed and suction settings
- **Durable:** Protective housing designed for easy cleaning between users; 3-year limited warranty
- **Whisper Quiet:** Discreet pumping anywhere
- **Safe:** Compatible with Ameda's Universal HygieniKit™, the only FDA-cleared closed system pump kit proven to protect against contaminants like bacteria, viruses and mold growth

## Ameda Pearl – Modern Design

### Hospital Grade/Multi-User Proven Durable Performance



## Ameda Pearl Technical Specifications

Suction Range	Stimulation: 50 - 150mmHg; 6 user-selectable settings Expression: 30 - 250mmHg; 12 user-selectable settings
Cycles per Minute (CPM)	Stimulation: 80 - 120 CPM; 3 preprogrammed settings Expression: 30 - 48 CPM; 6 user-selectable settings
Power Supply	Medical AC Adapter – Output: 15VDC, 2A, 30W Input: 100-240VAC, 50/60Hz, 1A 2600mAh Lithium-ion rechargeable battery; 3 hours per charge
Unit Weight	1.6kg
Unit Dimensions	20 x 18 x 15cm (L x W x H)
Preprogrammed pumping	Mode auto switches from stimulation to expression at 2 mins Auto Power off at 60 mins Vacusense Technology
Warranty	3-year limited



For more information  
about Ameda Pearl



## Pearl Pump Options and Accessories

Product	SKU	UPC
Pearl Breast Pump, AC power adapter, North America prong adapter, EN/FR/SP IFU, Quick Start Card	<b>AM.201W01</b>	<b>810725025864</b>
Trolley	<b>AM.201H01</b>	<b>810725025987</b>
Hard Shell Case	<b>AM.501M03</b>	<b>810725025932</b>
Multi-User Breast Pump Soft Carrying Bag – Detachable shoulder strap and insulated cooling compartment	<b>AM.501M02</b>	<b>810725025611</b>
Universal Diaphragm 2-Pack – 2 Universal Diaphragms	<b>AM.000X02</b>	<b>810725024867</b>
AC Power Adapter North America/AU/Philippines	<b>AM.201X05</b>	<b>810725026045</b>
Ameda Pearl Tubing Port Cover	<b>AM.201X08</b>	<b>810725024324</b>
Ameda Pearl Replacement Rubber Feet (set of 4)	<b>AM.201X03</b>	<b>810725026021</b>
Ameda Pearl Replacement Screen Overlay Panel	<b>AM.201X04</b>	<b>810725026038</b>
Ameda Pearl AU Prong Adapter	<b>AM.201X01</b>	<b>810725025895</b>



Ameda Pearl Trolley  
**AM.201H01**




Ameda Pearl Replacement Screen  
Overlay Panel  
**AM.201X04**



Ameda Pearl Hard Shell  
Transport & Storage Case  
**AM.501M03**



Ameda Pearl AC Power Adapter  
EU/UK/Asia  
**AM.201X06**



Ameda Pearl Replacement  
Rubber Feet  
**AM.201X03**



Ameda Pearl Tubing  
Port Cover  
**AM.201X08**



Ameda Pearl AC Power Adapter  
North America/AU/Philippines  
**AM.201X05**



For more info about  
Ameda's patented  
HygieniKits