



 **Amēda**<sup>®</sup>  
Australia



Trusted by Hospitals, Loved by Moms.



Ameda is the first and original name in breast pumps. Founded in 1942 by Swedish engineer Einar Egnell, the Ameda mission is to support all mums throughout their breastfeeding journeys with **safe, protective, trusted and accessible** products. Ameda is Who Code Compliant and consistently supports the International Code of Marketing Breast Milk Substitutes.

With a full portfolio of clinically-based solutions that includes the Platinum Multi-User Breast Pump, HygieniKit collection kit and Penguin Nutritional Warmer, Ameda is in the **hospital** ensuring mum get started right with products that are trusted by clinicians and lactation experts.

As mums transition **home** to begin life with their babies, Ameda is there with a full range of favorites that include Mya Joy Personal Pumps, Breast Care items and Breast Milk Collection & Storage options.

At Ameda, we recognize that each journey is unique, and we dedicate our efforts to providing safe, protective, and proven solutions and support for **all** mums.

**That's the Ameda Difference.**

# TABLE OF CONTENTS

Multi User Electric Breast Pumps.....	4/5
HygienKit Milk Collection System.....	6/7
CustomFit Flange System.....	8/9/10/11/12
Pump Replacement & Spare Parts.....	13
Personal Electric Breast Pumps.....	14/15
Milk Storage and Cleaning Products.....	16/17
Breast Care Products.....	18
Penguin Warmers & Accessories.....	19



# Ameda® Platinum™ Hospital Grade Breast Pump

Clinically Proven Performance for the Most Critical Patients. The **Ameda Platinum** is the only multi-user hospital breast pump to record milk volumes more than 1,000 mL/day during a 14-day clinical study. The Platinum combines clinically proven performance with a highly durable design that is perfect for a multi-use, hospital environment. The Ameda Platinum breast pump is clinically proven to be the most effective pump for building milk production and maintaining milk supply for all babies, including preterm and full-term infants.

## Ameda Platinum – Proven Design for Hospitals

CustomFit™ Flange System offers seven size options

Best flange fit for increased efficiency and comfort

3 easy to ready LCD screens that include level indicators and a timer

Independent Speed and Vacuum controls allow for complete customization

Patented HygieniKit™ Double Pump Kit featuring the only FDA-cleared protective viral barrier to protect breast milk

Built in bottle Holders secure collection bottles and milk bags to prevent spills

Timer makes following pumping protocols easier. Automatic shut-off after 60 minutes

Soft touch control panel for easy cleaning and intuitive operation

Product	Item Number	Sale Unit	UPC
Platinum Breast Pump with AU Power Cord	<b>17805T</b>	1 each	810725020876
Trolley Stand with locking casters and Mounting Kit for Ameda Platinum & Elite Breast Pumps	<b>17811</b>	1 each	810725021248
Breast Pump Soft Bag – for use with Platinum, Pearl & Elite	<b>501M02</b>	1 each	810725025994

Items sold separately. Contact your Ameda Sales Representative for additional information



**Platinum Trolley Stand and Mounting Kit 17811**  
Trolley stand with locking casters for Platinum & Elite Breast Pump not included.



**Platinum Power Cord Australia 24502149AU**  
(1) Power Cord



**Soft Bag for Platinum, Pearl & Elite 501M02**  
Lightweight, Microfiber Carrying Bag with Handle; Detachable Shoulder Strap; Insulated Cooling Compartment

# Ameda® Pearl® Hospital Grade Breast Pump



Powerful Performance in the Hospital and Home. Combines powerful performance inside with a protective shell on the outside. Pearl is designed to build milk supply amongst pump dependent, late preterm and term mums. Pearl's modern design combines the durability needed for the multi-user/hospital environment with lifestyle features mums expect. Pearl will serve the needs of hospitals and is the ideal multi-user product for WIC and rental stations. Pearl is perfectly positioned to support continuity of care in the Hospital and Home.

## Ameda Pearl – Modern Design

**Hospital Grade/Multi-User  
Proven Durable Performance**

Adjustable/Dimmable nightlight for late night sessions

Bottle Holders secure collection bottles and milk bags of all sizes up to 8oz. to prevent spills

Large, backlit LCD display with timer

Integrated carry handle

Adaptable for single or double pumping

6 suction levels in stimulation mode and 12 in expression mode. 6 independently controlled speeds in expression mode provide 72 separately configurable settings

CustomFit® Flange System offers seven size options

Built in memory and push button controls

Patented HygieniKit™ Double Pump Kit – universal solution featuring the only FDA-cleared protective viral barrier to protect breastmilk

Protective tubing port cover

2-minute auto switch to Expression Mode

Product	Item Number	Sale Unit	UPC
Pearl Electric Breast Pump	<b>201W01</b>	1 each	810725025864
Pearl Trolley	<b>201H01</b>	1 each	810725025987
Pearl Hard Case	<b>501M03</b>	1 each	810725024843
Pearl Power Adapter Australia	<b>201X01</b>	1 each	810725025895

Items sold separately. Contact your Ameda Sales Representative for additional information



Ameda Pearl Hard Shell Transport & Storage Case  
**501M03**



Ameda Pearl Power Adapter Australia  
**201X01**



Ameda Pearl Trolley  
**201H01**  
Breast Pump not included.

## Ameda HygieniKit® the only FDA-cleared Milk Collection System

### Ameda's Closed System = Safe Milk

Features a patented diaphragm in a closed system design – a proven barrier to bacterial and viral penetration.

- Ameda's HygieniKit System has the only Proven Airlock Protection™ protective barrier (laboratory-tested) reducing risk and simplifying pump use and care.
- Without the barrier, mold can easily grow in breast pump tubing when moisture becomes trapped inside it. Its solid


barrier keeps narrow pump tubing dry so no need to clean or replace it.

- A true “closed-system” has an effective barrier that also protects collected breast milk from potential contaminants, such as viruses and bacteria that may be present in the pump and/or kit tubing.

### An Industry First – Fewer Parts to Clean

All Ameda pumps are designed as a closed system with an **FDA-cleared** barrier to help keep milk safe.

The Ameda HygieniKit Milk Collection System comes pre-assembled and can be used immediately if purchased as a sterile kit, which keeps training simple and saves time for hospital staff, resulting in a reduction of nursing staff time.



### Clinically proven to protect against:

VIRUSES	BACTERIA	MOISTURE
Hepatitis A Hepatitis B HIV CMV Phi-X174 (Smallest known virus)	Staphylococcus aureus Streptococcus pneumoniae Pseudomonas aeruginosa Escherichia coli	Prevention of moisture eliminates the risk of dangerous mold growth Condensation is very common in the lead competitive kit and requires time/attention to remove

The Ameda HygieniKit has a patented silicone diaphragm that creates a barrier protecting the pump and tubing from penetration of bacterial, 2 and viruses (Hepatitis Band C and Human Immunodeficiency Virus HIV) 2, 3 from contaminated liquids inside the collection bottle system. The diaphragm in this product creates a barrier protecting liquids collected in the collection bottle system from penetration of bacteria, 1, 2 and viruses (Hepatitis Band C and Human Immunodeficiency Virus HIV) 2, 3 from the pump and tubing in the event of contamination of these components during a prior use or by a prior user. Note that the Ameda HygieniKit neither removes bacteria and viruses (Hepatitis Band C and Human Immunodeficiency Virus HIV) from the breast milk being collected into the collection bottle nor does the Ameda HygieniKit make it safe for infected mothers to feed their infants with milk collected using this device. 1,3 Testing performed in situ using actual HygieniKits and 5MB Piston Breast Pump with suspension of Phi-X174 bacteriophage and Escherichia coli, Staphylococcus pneumoniae and Pseudomonas aeruginosa, Procedure was adapted from ASTM 1671-97b "Standard Test method Resistance of materials Used in Protective Clothing 10 Penetration by Blood-borne Pathogens Using Phi-X174 Bacteriophage as a test System • 2 The use of this product neither guarantees or warranties against the transmission of Hepatitis Band C and Human Immunodeficiency Virus HIV



*Sterile*  
**Single HygieniKit Milk Collection System**

- ~ (2) Valves
- ~ (1) Diaphragm
- ~ (1) Tubing
- ~ (1) 4oz/118ml Bottles
- ~ (1) Tubing Adapter & Cap
- ~ (1) 25mm Flange
- ~ (1) Locking disc & Ring

Item Number	Sale Unit	UPC
AM.200H.02	1	810725024102



*Sterile*  
**Dual HygieniKit Milk Collection System**

- ~ (4) Valves
- ~ (2) Diaphragm
- ~ (2) Tubing
- ~ (2) 4oz/118ml Bottles
- ~ (1) Tubing Adapter
- ~ (2) Tubing Adapter Caps
- ~ (2) 25mm Flange
- ~ (2) Locking disc & Ring

Item Number	Sale Unit	UPC
AM.200H.05	1	810725024119



*Sterile*  
**Single HygieniKit Milk Collection system with One-Hand Breast Pump Adapter**

- ~ (2) Valves
- ~ (1) Diaphragm
- ~ (1) Tubing Adapter
- ~ (1) Silicone Tubes
- ~ (1) 4oz/118ml Bottles
- ~ (1) Tubing Adapter
- ~ (1) Tubing Adapter Caps
- ~ (1) 25mm Flange
- ~ (1) Locking disc & Ring

Item Number	Sale Unit	UPC
17150	1	810725025932



*Sterile*  
**Dual HygieniKit Milk Collection system with One-Hand Breast Pump Adapter**

- ~ (4) Valves
- ~ (2) Diaphragm
- ~ (1) Tubing Adapter
- ~ (2) Silicone Tubes
- ~ (2) 4oz/118ml Bottles
- ~ (2) Tubing Adapter
- ~ (2) Tubing Adapter Caps
- ~ (2) 25mm Flange
- ~ (2) Locking disc & Ring

Item Number	Sale Unit	UPC
17152	1	810725020593

11mm

13mm

15mm

17mm

19mm

21mm

22.5mm


25mm

# Ameda Flange Sizing Guide

Print this sheet at 100% scale for accuracy. Cut along the dotted lines to measure your nipple diameter. Be sure to measure BEFORE pumping.

## Step 1: Measure Your Nipple

Use the ruler on the tool to measure your nipple in millimeters as shown below. Do not measure the areola.



## Step 2: Find the Closest Hole

Gently place your nipple into the holes, starting with the smallest one that fits without squeezing.

## Step 3: Choose the Right Flange Size

Pick a flange that is 2-3mm larger than your nipple measurement.

If you're between sizes, start with the larger flange for comfort.

27mm

28.5mm


30.5mm

32.5mm

36mm


## Did you know?

Incorrect flange fit can affect breast pump suction! Using the right size breast flange can help make pumping more comfortable and efficient.




### Good Fit

During pumping, your nipple moves freely in the breast flange tunnel. You see space around the nipple. Not much areola is drawn into the tunnel with the nipple.



### Too Small

During pumping, some or all of your nipple rubs against the sides of the breast flange tunnel.



### Too Large

During pumping, more areola is drawn into the breast flange with your nipple. Your areola may rub against the side of the breast flange tunnel.



## CustomFit Flange System™

The Ameda CustomFit Flange System makes it easy to find a good flange fit by offering seven sizes that adapt to all Ameda Breast Pumps and Kits. Research has found that nearly one third of mothers who pump report greater comfort and more milk with a larger-than-standard breast flange size<sup>1,2</sup>. Ameda CustomFit™ Flanges have been designed to work with any Ameda HygieniKit® Milk Collection Kit. Simply transfer the valves and diaphragms to the new, better-fitting breast flange size and you are ready to pump.



Flexshield Areola Stimulator



Reducing Insert



XS

**Flexishield™ Areola Stimulator | 21.0mm**  
(Fits Inside Standard 25.0 mm Breast Flange)  
AM.10010

S

**Reducing Insert | 22.5mm**  
(Fits Inside Standard 25.0 mm Breast Flange)  
AM.10011.A

XXL • Custom Breast Flange



XL • Custom Breast Flange Insert



XL

**Custom Breast Flange Insert | 32.5mm**  
(Fits Inside XX-Large 36.0mm Breast Flange)

XXL

**Custom Breast Flange | 36mm**  
AM.17358M

1 Jones, E. et al. A randomised controlled trial to compare methods of milk expression after preterm delivery. Arch Dis Child Fetal Neonatal Ed 2001; R91-F95.

2 Wilson-Clay, B. and K. Hoover. The Breastfeeding Atlas, 4th ed. Austin, TX: LactNews Press, 2008.

L • Custom Breast Flange



M • Custom Breast Flange Insert



M

**Custom Breast Flange Insert | 28.5mm**  
(Fits Inside Large 30.5mm Breast Flange)

L

**Custom Breast Flange | 30.5mm**  
AM.17148PM

Standard Breast Flange



STANDARD

**Standard Breast Flange | 25.0mm**  
(Included with all Ameda HygieniKit Systems and Breast Pumps)

AM.10012.A

The Ameda Universal Soft Breast Pump Flange Inserts provide assured quality, comfort, and optimal breast pumping performance. Designed to provide a reliable fit for every mum, these authentic silicone inserts come in a variety of sizes to meet your unique needs. Whether you're just beginning your breastfeeding journey or are an experienced mum, Ameda's flange inserts help create a secure and comfortable pumping experience. Avoid the risks associated with knock-off products—trust Ameda for the comfort and efficiency you deserve in every session.

## Key Features:

The Ameda Soft Flange Inserts, available in sizes 11mm to 21mm, fit into your 25mm flange for a better, more comfortable fit. A proper fit can enhance milk expression effectiveness.

Made from soft, flexible silicone, these inserts provide enhanced comfort throughout your pumping session, reducing any discomfort often associated with pumping.

The customizable sizing helps create a secure seal during use, reducing the chances of leaks and promoting more effective pumping.

Designed for simple and quick insertion into existing breast flanges, the Ameda Soft Flange Inserts can be installed effortlessly for convenient adjustments whenever needed.

Made without BPA and DEHP, ensuring a safe and healthy pumping experience for both mum & baby.

Designed by Ameda, a leader in breast pump innovation, to deliver the performance and reliability you deserve.



SIZE	CODE
11mm	AM.120G09
13mm	AM.120G10
15mm	AM.120G11
17mm	AM.120G12
19mm	AM.120G13
21mm	AM.120G14

## CustomFit Flange

~ (1) Breast flange

Item Number	Description	Sale Unit	UPC
<b>10010</b>	21.0mm Flexishield	1 each	810725023112
<b>10011</b>	22.5mm Reducing Insert	1 each	810725023129
<b>10012</b>	25.0mm Standard Breast Flange	1 each	810725023143



## Dual Pack Flexishield™ Areola Stimulator

Gentle massaging action, softens the feel of the breast flange, reduces standard 25.0mm flange size to 21.0mm for maximum effectiveness for smaller nipples.

~ Without BPA or DEHP

Item Number	Sale Unit	UPC
<b>17212</b>	1 each (1 pair)	810725020128



## CustomFit Breast Flanges & Inserts

- ~ (2) Breast flanges
- ~ (2) Reducing inserts
- ~ Without BPA or DEHP

Item Number	Description	Sale Unit	UPC
<b>17148PM</b>	Large/Medium – 30.5mm flange/28.5mm insert	1 each	810725020104
<b>17358M</b>	XX-Large/X-Large – 36.0mm flange/32.5mm insert	1 each	810725020111





### Breast Pump Valves

~ (2) Valves

Item Number	Sale Unit	UPC
17146M	1 each (1 pair)	810725020012



### Adapter Cap

~ (2) Adapter caps

Item Number	Sale Unit	UPC
10005	1 each (1 pair)	810725023099

### Bottle – Polypropylene

~ (1) 4 oz./120 ml. bottle  
~ Without BPA or DEHP



Item Number	Sale Unit
24200017	1 each

### Silicone Tubing

~ (2) Tubes



Item Number	Sale Unit	UPC
10003	2 per Pack	810725023027

Full Case Quantity: 50/case



### Tubing Adapter

~ (2) Adapters

Item Number	Sale Unit	UPC
10004	1 each (1 pair)	810725023082



### Locking Ring and Disc

~ (2) Rings and discs

Item Number	Sale Unit	UPC
80001	1 each (2 sets)	810725023174

### Silicone Diaphragms

~ (2) Diaphragms



Item Number	Type	Compatible with	Sale Unit	UPC	Case Count
17147M	Standard	Platinum, Elite	1 each (1 pair)	810725020029	24 each/case



## Small AND Mighty!

Meet the Ameda Mya Joy PLUS – the newest member of the Mya Joy family of pumping solutions. Mya Joy PLUS combines the hospital strength performance mums demand with the incredibly small pump design they need. The thoughtfully designed pump comfortably fits in the palm of a hand and has whisper quiet operation. We” PLUSED” up the new Mya Joy PLUS with additional lifestyle features that include a long-lasting rechargeable battery, USB charging cord and convenient wrist lanyard. All of this from Ameda – the brand Trusted by Hospitals and Loved by Mums for over 80 years.



## Ameda® Mya™ Joy PLUS Double Electric Breast Pump with Tote

- ~ (1) Ameda Mya Joy PLUS Breast Pump
- ~ (1) AC Power Adapter with USB-C cord
- ~ (1) Dual HygieniKit Milk Collection System
- ~ (2) CustomFit™ flanges (25mm)
- ~ (2) CustomFit™ flanges (28.5mm)
- ~ (2) Extra valves and diaphragms (forHygieniKit)
- ~ (1) Tote bag

Item Number	Sale Unit	UPC
131A41	1 each	810725025581



## Pump'N Protect™ 6oz Breast Milk Storage Bags

Save time by pumping directly into the bag using integrated adapter holes to attach flange for direct pumping.

- ~ Tamper-proof seal for safety
- ~ Integrated adapter holes allow mum to pump directly into the bags with use of a pump adapter
- ~ Easy-close, double zip seal with leak-proof technology
- ~ Reinforced side seams for ultimate protection
- ~ Gusset bottom allows bag to stand or lay flat for storing
- ~ Easy to handle bag, BPA free
- ~ Easy to read markings, hold up to 6 fluid ounces

Item Number	Description	Sale Unit	UPC
<b>800M01*</b>	50-Count Bags (No Adapters)	50 ct/each	810725024997
<b>800M02*</b>	100-Count Bags (No Adapters)	100 ct/each	810725025000
<b>800A01**</b>	150-Count Bags with 2 Adapters (Bulk)	1 case	810725025017
<b>800X02***</b>	Standard Milk Bag Adapters	1 each (1 pair)	810725025321

## CleanEase Microwaveable Steam Sanitizing Bags

Ameda CleanEase bags are a convenient and efficient way to ensure breast pump parts are thoroughly sanitized within minutes.

- ~ Sanitize breast pump kit parts such as flanges, flange adapters, diaphragm, valves (no need to clean Ameda tubing, tubing adapter and cap)
- ~ Can also be used to store kit parts while on-the-go

Item Number	Sale Unit	UPC
<b>800M03*</b>	1 each (7 ct)	810725022085
<b>800H03**</b>	1 each (100 ct)	810725022030





## Breast Milk Storage Bottles

Bottles are not compatible with Ameda Mya. Universal threads allow for compatibility with all Ameda breast pumps and most other brands.

- ~ (4) 4 oz. polypropylene bottles with lock-tight sealing lids
- ~ Without BPA or DEHP
- ~ Dishwasher and freezer safe

Item Number	Sale Unit	UPC
<b>17244M</b>	1 each (4)	810725020036



## Cool 'N Carry Milk Storage Tote

Ideal for mothers who need another tote for transporting milk. Designed for use with standard and wide-mouth containers.

- ~ (1) Insulated carry bag
- ~ (3) Cooling elements
- ~ Without BPA or DEHP

Item Number	Sale Unit	UPC
<b>501M01</b>	1 each	810725025147



## Cool 'N Carry Milk Storage Kit

Bottles are not compatible with Ameda Mya. Ideal for mothers who need another tote for transporting milk to day care provider.

- ~ (1) Insulated carry bag
- ~ (6) bottles and lock-tight sealing lids
- ~ (3) Cooling elements
- ~ Milk storage guidelines
- ~ Bottles are without BPA or DEHP

Item Number	Sale Unit	UPC
<b>17076PM</b>	1 each	810725020043

Full Case Quantity: 6 each/case

Non-Sterile

## Cut-Out Nipple Shields

Ameda nipple shields feature a unique cut-out design which allows for an increased skin-to-skin contact between mother and baby.

Item Number	Description	Sale Unit	UPC
<b>60003</b>	Non-Sterile w/cut out - 16mm	1 each	810725025338
<b>60009</b>	Non-Sterile w/cut out - 20mm	1 each	810725025352
<b>60010</b>	Non-Sterile w/cut out - 24mm	1 each	810725025382



## DuoShell Breast Shells

- ~ (2) Aerated shells
- ~ (2) sets of backs (1 pair for sore nipples and 1 pair for inverted nipples)

Item Number	Sale Unit	UPC
<b>17231M</b>	1 set	810725020210



## Ameda MoistureGuard™ Disposable Nursing Pads

- ~ Extremely thin to fit discreetly into your bra to help prevent embarrassing leaks
- ~ The super absorbent core keeps your clothing dry, day or night by absorbing up to 35x its weight
- ~ Pads are contoured for a better fit and have an adhesive strip to keep them securely in place without showing

Item Number	Sale Unit	UPC
<b>602M01</b>	1 each/50 per box	810725024270

## Ameda Contoured Washable Pads

Ameda Contoured Washable Nursing Pads help prevent embarrassing wet stains on clothing with multi-layered construction for leak prevention. Pads feature a contoured shape for a discreet appearance and are machine washable for a repeated use.

Item Number	Sale Unit	UPC
<b>601M08</b>	1 each (8 pads)	810725024966



### Penguin Therma-Liner™

Updated design features:

- ~ Dual pouch liner separates nutrition from contact with water in a dry pouch surrounded by a sealed wet pouch
- ~ Improved, easy to seal with simplified instructions for use
- ~ Label space for patient information reducing chances of errors
- ~ Disposable and 100% recyclable
- ~ Note: For single patient/12 hour use

Item Number	Sale Unit	UPC
PNWZIP2DCSB	1 case (6 boxes of 50)	810725022115

### Penguin Deluxe Warmer Single Well

Delivering breast milk at ideal temperature is critical to infant health and nutrition. Ameda Penguin is the ONLY nutritional warmer that uses:

- ~ Smart Warming Technology™
  - 360° surround warming achieves ideal temperature
  - Gentle vibration mixes feeds throughout the warming cycle
  - Dual pouch liner separates nutrition from contact with water in a dry pouch surrounded by a sealed wet pouch
- ~ This eliminates environmental pathogens that can be introduced in a traditional water bath warming method
- ~ Blocks airborne environmental pathogens that can be introduced by fan blowing contaminants from the air

Item Number	Sale Unit	UPC
PNW00115AU	1 each	810725025031



### Penguin Perche

- ~ The Ameda Penguin Perche safely secures the entire Single Well Penguin warming system decreasing chances of accidental tipping over/ falls

Item Number	Sale Unit	UPC
PNWTKSHLV15	1 each	810725024072



CONNECT WITH US

Ameda Australia

Unit 4/62  
Borthwick Ave,  
Murarrie QLD 4172

1300 643 633

[www.ameda.com.au](http://www.ameda.com.au)