

Baby Crib - Variable Height, Gas Strut Operation

Features

- Height adjustable using a gas damper (710 - 1025mm)
- Push handle either end
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- Supplied with a clear plastic crib (PETG) & flame-retardant mattress (max capacity 6.0kg)
- 75mm twin-wheel braking castors

Dimensions (w x d x h)

- Overall Dimensions - 930 x 476 x 775 - 1025mm
- Mattress - 635 x 330 x 75mm

Standard Colours

- Paint - Grey White / Spearmint
- Work Flap - Spearmint

Accessories / Spares


- BA902 - Slide Out Debris Bin
- BA903 - Wire Basket - Internal Dimensions 420 x 170 x 220mm
- BA904 - Lift Up Flap - Dimensions 400 x 300mm
- 2PLCO014 - Clear Plastic Crib (Spare / Replacement)
- 2FACO052 - Mattress, Flame Retardant, Washable (Spare / Replacement)



BA902 - Debris Bin



BA903 - Wire Basket



BA904 - Lift Up Flap



BM.BCV/G

Item	Description	Weight
BM.BCV/G	Baby Crib - Variable Height, Gas Strut Operation	20.2kg