


Manufacturer



RCR MEDICAL
WHEN SECONDS COUNT



X8TTM



t o u r n i q u e t

EFFECTIVE. SAFE. EASY. LIGHT WEIGHT. LOW PROFILE.

1 MEDTM Technology

The patent-pending Maximum Efficiency DialTM provides you the most rapid and simplified windlass dial on the market. Self-locking, low profile, extreme adjustability for use of only gross motor skills.



2 D-Ring

X8T allows you to quickly tighten the strap and is developed for simple one-handed operation of the tourniquet.

3 Quick Release Buckle

X8T allows you to quickly wrap the tourniquet around a trapped limb or for quick release by medical professionals.



TourniqJuet

XS.T

CA'J



Width $\geq 1.5''$	1.47''	L.56',
*Effective Width $\geq 1.5''$	1.41''	1...00'
*Length $\geq 37.5''$	42..00''	37.80t'
Weight IPackaged	4.50 oz	2.63. oz
*Cubic Size ≤ 25.6 cubic inches	24...75''	28.18'J'
Time Required to Tighten	10.6 seconds	26..4 seconds
Occulsion, of Vessel (Visual)	✓	✓
Occulsion, of Vessel (Pressure)	Good	Good
Application using Gross Motor Skills, Only	✓	✗
Tactile Feedback	✓	✗
Incremental Adjustability	✓	✗
Audible Precision Tension Adjustment	✓	✗
Self-Engaging	✓	✗
Easy to Use in Low Light	✓	✗
Easy to Use with Gloves	✓	✗
Slip Resistant	✓	✗
Meets all Military Specs	✓	✗
Fall Safe	✓	✗

*military specifications

