

# CT-7 LEG TRACTION SPLINT

## NEXT GENERATION LEG KIT

The wait is over. After two years of development, the next generation CT-7 is here. Tested under the harshest conditions around the world, the CT-7 traction splint is the best choice for any medic.

REPLACING THE CT-6



## FEATURES

- A** Acme Screw Traction Technology  
Precise, Strong and Accurate
- B** Carbon Tubing
- C** Compact, Lightweight Design
- D** Super Clean & Easy to Use
- E** Applied In 30 Seconds



## IMPROVEMENTS & SPECIFICATIONS



### Improvements:

- No longer extends beyond the patients foot
- Simplified layout
- Less velcro
- One piece ankle hitch
- Compact & lighter than any other traction splint on the market
- Increased diameter of carbon tubing
- Applied in under 30 seconds
- Easier to repack, uses new Tactical Bag



### Specifications:

Size: 30.48cm x 10.16cm x 10.16cm  
Weight: 0.45kg  
Materials: Carbon Tubing, Nylon & Elastic Webbing  
Colors: Black, Orange (Straps)