

Sofia automatic recliner



the *Sofia* Automatic recliner chair combines maximum comfort and versatility with optimal functionality suitable for the most demanding care giving environments. Sofia is more than a patient recliner; it's a transportation chair, a sleep recliner and a treatment chair.

○ ● ● ● ● positions

reclined position

thanks to the 2 independent motors, Sofia, features multiple reclining positions controlled with an easy to use, hand held, wire controller



upright position

Sofia offers an ergonomic design with a reinforced internal structure made of steel that provides outstanding durability



sleeper

when reclined flat, Sofia ensures the highest degree for the user. Patients and other users will feel that they are resting in a bed rather than a chair



Lifting position (optional)

Patients often suffer from reduced mobility. Helping patients get out of the chair will not only make life easier for the patients but will also reduce caregiver injuries



○ ● ● ● ● standard features

push bar

an innovative design that facilitates the movement and rotation of the chair

urethane arm caps

provide users a comfortable grip and protection for the arms

moisture barrier

protection from fluids such as blood, sweat and urine, with minimum maintenance

bar caps

allows for fixing a multiple accessories bar which enables attaching features such as an IV pole or blood donor arm support

moveable armrests

the armrests can be lifted up to transfer the patient to or from the bed, without any effort

side pocket

useful for holding extras such as the chair controller or a cell phone or television remote

controller

easy to use, hand held, wire controller

push bar and arm cap palette
personalize your chair with wengue, orange or anthracite colours



pull out food rest

for additional support for legs when seated

zero turning radius

eases mobility in small spaces

casters

silent and resistant double casters with individual locks

○ ● ● ● ● optional features

battery

a rechargeable battery will provide user with a much greater level of autonomy

multiple accessories bar

allows for the attachment of many different accessories such as hanging or drainage bags, an articulating table, a blood donor armrest or an IV pole holder

swing away table

mechanism included

Decam multiple accessories bar needs to be previously installed in the chair

oxygen tank holder

attaches to the back of the chair to carry an oxygen cylinder

IV pole

stainless steel pole that fits into the adjustable IV Pole holder

IV Pole holder

attaches to the multiple accessories bar. It is height, tilt and direction adjustable

blood donor arm holder

the adjustable donor arm holder provides a comfortable arm position for a donor or patient
multiple accessories bar must be installed

Head rest

the adjustable head rest is removable, and provides greater comfort to the user

○ ● ● ● ● technical specifications

Sofia automatic

dimensions	overall		seat		arm height		full recline	
width	62,5 cm	24,4"	49,5 cm	19,5"	70 cm	27,6"		
depth	75 cm	29,5"	49,5 cm	19,5"			170 cm	66,9"
height	127 cm	50"	52 cm	20,5"				
weight	70 Kg	154 lb						
max load	150 Kg	330 lb						

Sofia automatic plus (wider seat)

dimensions	overall		seat		arm height		full recline	
width	67 cm	26,4"	55,8 cm	22"	72 cm	28,3"		
depth	75 cm	29,5"	47,5 cm	18,7"			170 cm	66,9"
height	127 cm	50"	52,5 cm	20,7"				
weight	89 Kg	196 lb						
max load	150 Kg	330 lb						

○ ● ● ● ● fabrics

The Decam textile collection, based on Spradling® Diamante and Silvertex patterns is designed specifically to meet the strictest requirements of the healthcare environment by creating **effectives barriers against germs, abrasion and stains**



* visit our website for a full list of available colors

Our standard collection meets the following features:

- | | | | |
|--------------------------|------------------|------------------|--------------------------|
| Abrasion resistance | Urine protection | Sweat Protection | 100% Polyester |
| Antibacterial Protection | Blood Protection | Mildew Resistant | Sulphide Stain Resistant |
| Antifungal Protection | Stain Protection | | |

This brochure provides an overview of the main and benefits of the DECAM product. For a full list of products available, including specifications, fabrics, and finishes, please visit www.decamlux.com. For pricing enquiries contact your local DECAM sales representative.