



“Your HotBag saved her life.

Our crew began warming a hypothermia victim with a blanket from your HotBag after she fell through the ice. By the time they reached the hospital her core temperature had been raised from the low ~26.6°C to 35°C.”

Brett Murphy
Regional EMS

BLANKET WARMERS

DESIGNED FOR PRE-HOSPITAL TRAUMA AND MEDICAL CARE



3 Sizes To Choose From

2-3
blankets

SMALL

Order #
MED-BLANKET-1LS

4-5
blankets

MEDIUM

Order #
MED-BLANKET-XWLS

6-7
blankets

LARGE

Order #
MED-BLANKET-2LS

All of our Blanket Warmers are

IEC 60601-1

Medical Accessory Compliant

- Interior Air Temperature of 54-60°C
- **Low Voltage** (12vDC) <50W
- Removable Heating Elements
- Flat Velcro Closure
- Durable Nylon Construction (Bag)
- **Washable** (Bag Only)
- Lightweight/Portable
- 12 Month Warranty

Always ready

When plugged into an available power source The Electric HotBag™ keeps blankets at a safe, constant temperature, always ready for use. (Power supply and adaptor sold separately.)

AC Power Supply

Order #
HBPS-1



12v Cord Adaptor

Order #
HBPC-AMP-72



Orders, Quotes & Questions

A real person is waiting to speak with you. Contact us today!
8am-4:30pm AEST, Mon-Fri
www.midmed.com.au
sales@midmed.com.au

1300 643 633

Office Number



 **midmed**
providing positive solutions