

BREWER ACCESS

EXAM TABLE



Large Step

- Most secure step
- 45cm x 52cm surface
- Stable four leg design
- Fully adjustable
- Allows patients to safely and easily get on and off table



Comfortable, Ergonomic Upholstery

- Seamless, easy to clean surface
- Resistant to medicinal type stains
- Available in 15 standard colors
- Cal-133 fire code compliant upholstery available
- Simple pneumatic lift backrest



Pass-Through Drawers

- Patented pass-through side drawers allow right or left access
- Eliminates right and left hand setup problems



Ergonomic Stirrups

- Extendable, adjustable
- Conveniently stored in the table
- Manually extended through infinite positions and four lateral positions



Paper Roll Storage

- Convenient table paper storage
- Accepts 53.3cm x 8.9cm and 45.7cm x 8.9cm paper rolls
- Storage for up to 6 additional paper rolls under backrest



Standard Outlets

- Duplex
- Hospital-grade
- 115VAC, 60Hz
- 5.8 amps max

Advanced Safety

- Large patient step
- Stable four leg step design
- Largest step surface available
- Step surface area: 45 cm x 51.8 cm

Patient Comfort

- Simple pneumatic backrest for easy patient positioning
- Multi-density, bolstered support upholstery
- Leg extension/plastic debris tray
- Adjustable ergonomic stirrup

Abundant Storage

- Patented pass-through side drawers accessible from right or left side of the table
- Drawer stops (provided) allow drawer to open from only one side if desired
- Front drawer: 50.3 cm x 48.8 cm x 9.7 cm
- Side drawers: 57.4 cm x 51.3 cm x 17.8 cm
- Covered storage for up to 6 additional paper rolls
- Paper roll storage capacity: 58.4 cm x 58.4 cm x 8.3 cm
- Slide through holder accepts 53.3 cm x 7.6 cm and 45.7 cm x 7.6 cm paper rolls

Table Specifications

- 500 lb (226 kg) weight capacity
- Standard 5-year warranty
- Height: 81.3 cm
- Width: 71.1 cm
- Length: 145.5 cm - 180.3 cm
- Leg extension length: 38.1 cm
- Step Height: 20.3 cm
- Step Length: 45 cm
- Step Width: 51.8 cm
- Hospital-grade duplex electrical outlets
- 115VAC, 60Hz, 5.8 amps max
- Outlets can be put on left or right side
- Optional 5° pelvic tilt / 60-Watt front drawer heater
- Easy to clean, seamless upholstery
- 15 standard upholstery color options
- Cal-133 fire code compliant upholstery available in 4 colors
- Custom colors also available

Abundant Storage Capacity

- Largest storage capacity available
- Patented pass-through side drawers
- Two additional spacious front drawers
- Convenient table paper roll storage

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

UL 60601-1



CONTACT INFORMATION

Unit 4, 62 Borthwick Avenue, Murarrie, 4172, Queensland, Australia

 1300 643 633

sales@midmed.com.au



www.facebook.com/midmedau/



www.midmed.com.au



www.linkedin.com/company/midmed