

54 Series

RECLINER/TRANSPORTER

RATED UP TO 136.08 KG

Part #547 MANUAL RECLINER
Part #54P POWER RECLINER*



Shown with optional second side table with cup holder.

Patient comfort, ease of use, dependability and long life make our classic 54 Series the “workhorse” recliner/transporter of the industry. The 54 Series is engineered for your most demanding applications, including around-the-clock use. Durable five-inch casters roll easily during transport.



Z-spring construction for greater comfort.

All steel, powder-coated frame for durability.

Patients recline easily by themselves.

STANDARD FEATURES

- Stationary arms
- Fold-away side table (one)
- Push handle
- Retractable foot tray
- 12.7cm thermoplastic casters
- Wall-hugger mechanism (infinite recline positions)
- Quick-release seat




Wall-hugger recline mechanism



Quick-release seat

STANDARD VINYL SELECTION



Natural T09

Aquamarine T19

Bonnie Blue T21

Iced Mint T26

Tearose T43

Colonial Blue T45


Gray T56

Cranberry T96

Printed color swatches may not be an exact match to actual fabrics and are enlarged to show detail.


*The 54P (without the heat & massage option) is UL 60601-1 certified.

SPACE REQUIREMENTS



SPECIFICATIONS

Weight capacity	136.08 KG
Back height (from floor/seat)	121.92cm/73.66cm
Overall width (2 tables down/up)	81.28cm/132.08cm
Seat width	57.15cm
Arm height (from floor, front/back)	76.2cm/68.58cm
Seat height/depth	57.15cm/45.72cm
Upright overall depth (w/handle)	119.38cm
Thermoplastic casters	12.7cm



OPTIONS

Trendelenburg (T)

Break-away arm – Note patient left (L) or right (R) side

Additional fold-away side table (S7) –
 not available with break-away arm option

Cupholder – on all tables (C), on one of two tables,
 specify CUP-HOLDER (right or left)

Premium casters – Set of 4 (K)

Heat (J) – (547 manual model only)

Heat and massage (JM)

CAL 133 fire barrier (3)**

Vinyl upgrades

Champion, graded vinyl selections***

Customer's own material (COM)

PVC-free material

Dual-hand control (DUALACT-PWR) – (54P model only)

Battery backup (BATTERY) – (54P model only)

**Contact Midmed for a list of fire-tested vinyls.

***For available patterns and colors for graded vinyl selections, please contact your Midmed sales representative.



BREAK-AWAY ARM



TRENDLENBURG



VINYL UPGRADES

ACCESSORIES

Provide your patient the best care and comfort with an over-the-lap tray, cup holders, head pillow/lumbar support, cushion wedge/arm support, footrest extender, TV mount, nocturnal sleep pad, side guards, IV poles

or an oxygen tank holder. For a complete list of our accessories, please see our *Accessories* brochure or ask your Midmed sales representative for more information.