

# MARSDEN M-999

Patient Transfer Scale

UK PATENT  
Patent number  
2543835

EU PATENT PENDING


US PATENT PENDING

CANADA PATENT PENDING

HONG KONG PATENT PENDING



- One device for lateral transfer and weighing of patients
- Weigh immobile patients in seconds: less stress, more dignified
- Administer drug dosages with confidence
- Class III Approved and MDD Approved





## **WEIGHING CRITICAL PATIENTS IS A CRUCIAL PROCESS. BUT IT'S ALSO A CHALLENGING ONE.**

Weighing patients correctly is essential to ensuring accurate drug dosages and administering the correct treatment. But for patients requiring urgent, immediate medical attention this is not always easily achievable. How do you weigh critical and immobile patients accurately, quickly and easily?

## **THE PATIENT TRANSFER SCALE: WELCOME TO THE FUTURE OF PATIENT WEIGHING**

The idea of a transfer board with integrated scale was born of frustration: An acute nurse working for NHS Lanarkshire was determined to find an easier way to accurately weigh immobile patients. The result, the M-999, revolutionises the weighing of patients in hospitals. It doesn't just make getting an accurate weight reading easy: it actually makes weighing part of the already-existing process of transferring a patient from stretcher to bed, or bed to bed. This means it's quick and easy, and saves you time.

## **GET AN INSTANT, ACCURATE WEIGHT – WHEN YOU NEED IT MOST**

Healthcare professionals have a duty of care to ensure a safe environment for patients. When knowing the weight of a patient is absolutely crucial to delivering the right care and treatment, the M-999 embeds confidence and accuracy in your processes. Use the M-999 throughout the hospital pathway to enhance outcomes – particularly for trauma patients in Emergency Departments where time is crucial, Stroke patients to ensure a quick 'door to needle' time, plus Cardiac, ICU, HTU and X-ray patients.

## EASY TO USE MEANS MINIMAL TRAINING AND EFFECTIVE RESULTS

We've made the M-999 as lightweight as possible so it's easy to take it to exactly where you need to and get an instant, accurate weight. Grab handles conveniently positioned down each side mean it's easier to move. And because it's all about a simple, precise weight reading, there are no complicated functions to get your head around.

## TESTED BY 30 HOSPITALS - AND WE'VE SHAPED THE DESIGN TO ANSWER EVERYONE'S NEED

As part of the initial research and development of the M-999, over 30 hospitals saw, used and fed back their thoughts on the scale. Trialling the scale as it would be used in practice, and getting feedback from relevant healthcare professionals and the wider multi-disciplinary team, contributed massively to the design and functionality of the production version. Subsequent intensive testing has meant the M-999 delivers precision weighing to Class III medical standards.


## KEY FEATURES

- MDD Approved/Class III Approved
- Powered by rechargeable battery
- Large, bright weight display
- Grab handles for portability
- Easy to clean - wipe clean surface
- Hold function

## SPECIFICATION

- Capacity/graduations: 250kg x 500g
- Power: Rechargeable battery
- Dimensions: 1800mm x 700mm x 30mm
- Weighing area dimensions: 1800mm x 510mm





## CONTACT INFORMATION

Unit 4, 62 Borthwick Avenue, Murarrie, 4172, Queensland, Australia

 **1300 643 633**

[sales@midmed.com.au](mailto:sales@midmed.com.au)



[www.facebook.com/midmedau/](https://www.facebook.com/midmedau/)



[www.midmed.com.au](http://www.midmed.com.au)



[www.linkedin.com/company/midmed](https://www.linkedin.com/company/midmed)