

Bottles

Safe & easy

A close-up photograph of a woman with blonde hair, wearing a white shirt, feeding a baby. She is holding a clear plastic Materni bottle with a white cap and nipple. The baby is lying down, wearing a light blue shirt, and is drinking from the bottle. The woman has a ring on her finger. The background is a soft, out-of-focus white.

Materni's ready to use
bottle range, designed for
the collection and storage of
expressed or formula milk



Collection and storage bottles for the hospital environment

- Wide range: 50, 80, 120 and 240 ml
- Indelible and accurate scaling
- Universal bottle, hospital pump compatible
- Ready to use (RTU)
- Clean by manufacturing (CBM)
- Freezer safe
- Zero-chemicals: ETO, BPA, phthalate and latex free
- Hygienic: hermetically closed till point of use
- Conveniently packed
- Recyclable
- Food and medical approved polymers



Materni safe & easy standard product range RTU bottle and closure

Article	Article code	Screw thread	Microbiology	Bundles/Box	Items/Box
Materni RTU* bottle with closure - 50 ml	MEPA0106	EU	CBM**	3 x 70	210
Materni RTU* bottle with closure - 80 ml	MEPA0107	EU	CBM**	6 x 35	210
Materni RTU* bottle with closure - 120 ml	MEPA0108	EU	CBM**	6 x 35	210
Materni RTU* bottle with closure - 240 ml	MEPA0109	EU	CBM**	4 x 35	140

* RTU: Ready to use.

**CBM: Clean by manufacturing, clean room ISO 6 produced and packed. Microbiologically clean.