

FlowSpoon™ FlowCup™

Breastfeeding aids

Materni's innovative safe and easy
breastfeeding aids to support the
mother and child from the very
first minute



Innovation in breastfeeding support

Materni's proprietary simple but efficient feeding aids support breastfeeding for hospital and home use from the very first moment.

Colostrum administration

The soft and flexible FlowSpoon™ allows safe collection and efficient administration of colostrum milk. The spoon is individually packed and protected until point of use. It can hold up to 6 ml of milk.

Safe and controlled cupping

Materni FlowCup™ is a unique, simple and affordable cupping device for the administration of expressed breast milk. FlowCup™ helps to avoid nipple-teat confusion and allows feeding at the natural pace of the baby, on top of being cheiloschisis friendly. No leakage or spillage. No milk-saliva mixing. Its easy design and use offers feeding flexibility for the professionals and both parents.

We offer 2 designs: M1 was originally designed with a bridge to help control the milk flow in order to follow the baby's pace. The M2 cup was developed in close collaboration with healthcare professionals. Its design permits administration till the last drop while still offering a slow and controlled release.



Materni FlowSpoon™ and Materni FlowCup™ standard product range

Article	Article code	Screw thread	Microbiology	Items/Box
Materni FlowCup™ - M1 (bridge) (assembled)	MEPE0111	EU	CBM*	140
Materni FlowCup™ - M2 (assembled)	MEPE0113	EU	CBM*	140
Materni FlowSpoon™	MEPA0012	n/a	CBM*	125

* CBM: Clean by manufacturing, clean room ISO 6 produced and packed. Microbiologically clean.