

INTRODUCING
The FIRST
Non-absorbent
Rapid Drain
Decon/Hazmat
Sled with
100% Sealed
Decon Straps



+ MED SLED

NO ONE LEFT BEHIND

DECON / HAZMAT

BE PREPARED. Med Sled® ensures easy movement of a non-ambulatory person through the Decon Shower process, and safely transports to a cleared zone.
(Not used for vertical evacuation.)

DECON / HAZMAT MED SLED®

- + Decon Showers
- + Eliminates hard backboards
- + Eliminates pinch zones
- + Easy on roller beds
- + Chemicals
- + Pathogens Spills
- + Caustic
- + Dirty Bomb

#MS36DC Standard Decon Med Sled® 36"W x 87"L
Bariatric and Youth Sizes also available.

Patent Pending



KEY BENEFITS

- + **NEW** 100% Fully Decon capable
- + Easy to lift onto Decon Shower rollers
- + Can be used as a two or four-man stretcher
- + Expels Decon water rapidly
- + Eliminates hard backboards and pinch zones



#MSSD4DC

Decon / Hazmat 4-Sled
Storage Unit
(Sold Separately)
16"W x 74"H x 8 3/4"D

- + **Accommodates:**
(4) MS36DC Sleds
stacked vertically 2
on top of 2