

LifeBox50

“BRIEFCASE”

The next generation of EMS transport containers is here. The key lies within the Carbon panels that produce insulation values of R50; 2x more powerful than anything else on the market. This allows for more payload space and greater efficiency in both hot and cold environments.

The LifeBox is designed for the First Responder to transport anything he needs that is temperature sensitive; without having to worry about a heavy, tethered powered unit. The Responder is not restrained by temperature changes over time. The LifeBox is able to maintain its core temperature for unprecedented durations in the most extreme ambient conditions. The standard box can hold 8 1 Liter IV solutions at temp for over 50 hours. Superior insulation, combined with rugged cases for protection, make it the number one choice for any medic.



Solid phase change packs to maintain temperature.



	LifeBox, Standard
Payload Capacity	6,709 (cm cubed)
Payload Dimensions	35.5 x 21 x 9cm
Exterior Dimensions	46 x 34 x 17cm
Tare weight (lbs)	5 kgs.
Outer Container	Rugged Pelican/Storm Case
Insulator	Carbon Aerogel Vacuum Panels
Before Use	Ready to Go
Use for	Blood, Cold Therapy, Organs, Tissue, IV's, Meds

Available in several sizes and colors, to best fit your needs. Please call us with any questions or concerns. Patent Applied For.



Equipped with digital temperature display; battery and solar powered!



PH-07 3348 9155 FAX-07 3348 9950
 WEB-www.midmed.com.au
 EMAIL-steve@midmed.com.au

