

# QinFL<sup>W</sup>

**NO MORE COLD BLOOD & FLUIDS**

Standardising Patient warming throughout your facility with a system that delivers exactly what is required - the Qinflow Blood Fluid Warmer the *ONLY* Warmer of Choice. Now validated to effectively warm intense intermittent flows generated via various intermittent transfusion techniques.

## KEY BENEFITS OF WARMING PATIENTS:

- A) Reduce infection rates
  - B) Reduce the length of stay
  - C) Lower medication costs
  - D) Will help save lives
  - E) Intermittent flow capability
  - F) Visible temperature display & alarm
- Savings thousands of dollars and improving overall patient outcomes.

Warrior  
w/- battery



The Warrior Battery solution for critical care transport. Suitable for Emergency, Theatres & ICU

Warrior AC



The Warrior AC solution for Emergency Operating Rooms & ICU

Compact DU



A disposable warming unit that fits all IV giving sets in the market, no additional extension sets required.

Extension Cable



Ability to have the warming units away from the patient for increase patient accessibility, whilst maintaining warming capacity.

## CONFIDENCE IN WARMING TECHNOLOGY:

- ▶ How are you guaranteed that the warming device you are currently using is doing what it claims.
- ▶ The Qinflow Warmer is the Blood Fluid Warming Unit that can deliver what it says.
- ▶ Under all situations Pre Hospital, Emergency, Theatres and ICU, this warmer can perform and has been tried and tested in the most extreme of events.
- ▶ Delivering low volumes to the very young and elderly to the highest of volumes required for trauma and surgical conditions the Qinflow Warmer will perform.
- ▶ The disposable cartridge is capable of withstanding pressures of up to 700mmHg and the fluid pathway is uninterrupted thus offering a safe delivery of fluids and blood.
- ▶ It has the options of AC Power and Battery so the warming journey is uninterrupted for transport.
- ▶ The Qinflow warming cartridge is made from Stainless Steel and can be disposed of easily in general contaminated waste.
- ▶ ANZSBT Australian College of Nurses guidelines for the administration of blood products.
  - Visible temperature display
  - Audible warning alarm for recording temperature during transfusion.

For more information contact us on **1300 643 633**  
or **sales@midmed.com.au**

