



**TECHLEM**  
STRETCHERS

 1300 643 633

 [sales@midmed.com.au](mailto:sales@midmed.com.au)  [www.midmed.com.au](http://www.midmed.com.au)

# TABLE OF CONTENTS

<b>1</b>	<b>ABOUT US</b> OUR HISTORY & OUR MISSION
<b>2-3</b>	<b>FEATURES &amp; ACCESSORIES</b> ANATOMY OF THE TECHLEM STRETCHER . OPTIONAL ACCESSORIES
<b>4-8</b>	<b>SERIES 5000 STRETCHERS</b> TRANSPORT . BARIATRIC . EXTENDED STAY . FIXED HEIGHT . X-RAY & TRAUMA
<b>9</b>	<b>SERIES 7600 STRETCHER</b> PATIENT LIFT & TRANSFER
<b>10-12</b>	<b>SERIES 4500 STRETCHERS</b> FOLD-RAIL . COMPACT . GYNAECOLOGY
<b>13</b>	<b>CADAVER STRETCHERS</b> CADAVER . BARIATRIC CADAVER
<b>14</b>	<b>COLOUR GUIDE</b> BUMPER ACCENTS . I.D. ACCENTS
<b>15</b>	<b>I.V. STAND</b> MODEL 802A
<b>16-17</b>	<b>ACCESSORIES GUIDE</b> FOR ALL TECHLEM STRETCHERS



1300 643 633

PO BOX 508 Morningside Q 4170

# ANATOMY OF THE TECHLEM STRETCHER



- 1 "Drop Under" Side-Rails
- 2 Protective Shroud
- 3 Four-Sided Access, Hands-Free Lift/Lower/Trendelenburg/Reverse Trendelenburg
- 4 8" All-Surface Casters
- 5 Personal Belongings & Horizontal O2 Tank Holder
- 6 "Accu-Steer" Dual 5th Wheel



- 7 Pneumatic or Crank Fowler
- 8 "Fold Down" Side-Rails
- 9 Bumper and Colour Accents
- 10 "Posi-Lock" Braking

# OPTIONAL ACCESSORIES



- 1 Removable Food Tray (5007)
- 2 Permanent (5010) or Removable (3113) I.V. Pole
- 3 Push & Pull Handles (5009)

- 7 Vertical Oxygen Tank Holder (4012)
- 8 Knee-Lift (5000KG)
- 9 4" Viscolastic Mattress (5174P)
- 10 "Posi-Lock" Braking

# TRANSPORT & GENERAL USE 5000

THE PREMIER TECHLEM STRETCHER FOR TRANSPORT AND GENERAL USE FEATURES AN INTEGRATED OXYGEN TANK HOLDER AND PATIENT STORAGE, 5TH WHEEL TECHNOLOGY FOR STEERING ASSIST, AND MORE THAN 20 CUSTOMIZABLE OPTIONS FOR DESIGN AND FUNCTIONALITY.



## MOVEMENT

**HEIGHT ADJUSTMENT:** 56CM TO 90CM  
**HANDS-FREE TRENDLENBERG:** +18° TO -18°  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
20CM ALL-SURFACE, LOCKING CASTERS  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
DROP UNDER SIDE-RAILS

## DIMENSIONS

**TOTAL DIMENSIONS:** 206CM X 83CM  
**LITTER DIMENSIONS:** 206CM X 75CM  
**FULL LENGTH SIDE-RAILS:** 191CM X 32CM  
**HIGH DENSITY FOAM MATTRESS:** 193CM X 72CM X 8CM

**WEIGHT CAPACITY:** 700LBS / 318KG

**STANDARD:**  
HIGH DENSITY FOAM  
MATTRESS WITH CONDUCTIVE COVER

**OPTIONAL:**  
VISCOLASTIC, PRESSURE  
SENSITIVE MATTRESS, WITH FOUR-WAY,  
ELASTATEX STRETCH COVER

# BARIATRIC 5000BR

SAFE, RELIABLE, AND NIMBLE. EQUIPPED WITH TECHLEM'S 5TH WHEEL TECHNOLOGY AND HIGH TENSILE COMPONENTS, OUR BARIATRIC STRETCHER BOASTS A 750 POUND WEIGHT CAPACITY.



## MOVEMENT

**HEIGHT ADJUSTMENT:** 56CM TO 81CM  
**HANDS-FREE TRENDELENBERG:** +12° TO -12°  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
20CM ALL-SURFACE, LOCKING CASTERS  
DROP UNDER SIDE-RAILS  
OPTIONAL FOUR-SIDED ACCESS

## DIMENSIONS

**TOTAL DIMENSIONS:** 206CM X 90CM  
**LITTER DIMENSIONS:** 206CM X 81CM  
**FULL-LENGTH SIDE RAILS:** 191CM X 41CM  
**HIGH-DENSITY FOAM MATTRESS:** 193CM X 81CM

**WEIGHT CAPACITY:** 750LBS / 341KG

**STANDARD:**  
HIGH DENSITY FOAM  
MATTRESS WITH CONDUCTIVE COVER

**OPTIONAL:**  
VISCOLASTIC, PRESSURE  
SENSITIVE MATTRESS, WITH FOUR-WAY,  
ELASTATEX STRETCH COVER

# EXTENDED STAY 5000W

VERSATILE, FUNCTIONAL, AND EXCEPTIONALLY COMFORTABLE. OUR EXTENDED STAY STRETCHER IS EQUIPPED WITH AN EXTRA-WIDE FRAME, DROP UNDER SIDE-RAILS, AND TECHLEM'S DISTINGUISHED VISCOLASTIC MATTRESS TECHNOLOGY.



PRESSURE SENSITIVE FOAM  
MATTRESS WITH VISCOLASTIC,  
FOUR-WAY STRETCH COVER

## MOVEMENT

**HANDS-FREE TRENDELENBERG:** +15° TO -15°  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
20CM ALL-SURFACE, LOCKING CASTERS  
DROP UNDER SIDE-RAILS  
OPTIONAL FOUR-SIDED ACCESS

## DIMENSIONS

**TOTAL DIMENSIONS:** 206CM x 90CM  
**LITTER DIMENSIONS:** 206CM x 81CM  
**MATTRESS DIMENSIONS:** 193CM x 81CM

**WEIGHT CAPACITY:** 700LBS / 318KG

# FIXED HEIGHT 1100 / 1000N

TECHLEM'S VALUE STRETCHER OPTION, OUR FIXED HEIGHT STRETCHER IS EQUIPPED WITH MANY OF THE SAME FEATURES AS OUR TRANSPORT STRETCHER, WITHOUT THE ADJUSTABLE HEIGHT.



FULL WIDTH, HIGH DENSITY  
FOAM MATTRESS WITH  
A CONDUCTIVE COVER

## MOVEMENT

**HANDS-FREE TRENDELENBERG:** +18° TO -18°  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
20CM ALL-SURFACE, LOCKING CASTERS  
DROP UNDER SIDE-RAILS

## DIMENSIONS

**REGULAR**  
**TOTAL DIMENSIONS:** 206CM X 83CM  
**MATTRESS DIMENSIONS:** 191CM X 72CM  
**NARROW**  
**TOTAL DIMENSIONS:** 206CM X 74CM  
**MATTRESS DIMENSIONS:** 191CM X 64CM

**WEIGHT CAPACITY:** 700LBS / 318KG

# PATIENT LIFT 7600

SIMPLE AND EFFICIENT IN DESIGN, OUR PATIENT LIFT AND TRANSFER STRETCHER IS ADAPTABLE TO YOUR CLINICAL ENVIRONMENT. HANDCRAFTED FROM HIGH GRADE MATERIALS, THE 7600 IS HIGHLY DURABLE, MAINTAINING ITS FORM AND FUNCTION OVER DECADES OF USE.



HIGH DENSITY PATIENT TRANSFER SHEET, HIGH GRADE STEEL SHEET BUCKLES, AND TELESCOPING I.V. POLE INCLUDED

## MOVEMENT

**HEIGHT ADJUSTMENT:** 30" TO 42"  
13CM ALL-SURFACE, LOCKING CASTERS  
5TH WHEEL FOR STEER ASSIST  
CRANK HANDLE FOR EASE OF USE

## DIMENSIONS

**OUTSIDE DIMENSIONS:** 196CM X 79CM  
**INSIDE DIMENSIONS:** 188CM" X 62CM

**WEIGHT CAPACITY:** 300LBS / 136KG

# TRAUMA & X-RAY 5000XR

FOR USE IN TRAUMA AND X-RAY ENVIRONMENTS. EQUIPPED WITH A FIBRE RESIN FILM SURFACE AND HIGH-DENSITY RADIOLUCENT MATTRESS. A TECHLEM TOTAL-VIEW DESIGN, WITH FULL SURFACE ACCESSIBILITY FOR CASSETTE PLACEMENT.



**STANDARD:**  
HIGH DENSITY FOAM  
MATTRESS WITH CONDUCTIVE

## MOVEMENT

**HEIGHT ADJUSTMENT:** 56CM TO 90CM  
**HANDS-FREE TRENDLENBERG:** +18° TO -18°  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
FOLD DOWN SIDE-RAILS  
20CM ALL-SURFACE, LOCKING CASTERS  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM

## DIMENSIONS

**TOTAL DIMENSIONS:** 206CM X 83CM  
**LITTER DIMENSIONS:** 206CM X 72CM  
**MATTRESS DIMENSIONS:** 191CM X 74CM X 8CM  
**DROP-RAIL SIDE-RAILS:** 154CM X 32CM

**WEIGHT CAPACITY:** 700LBS / 318KG

# FOLD RAIL 4500

RELIABLE, MOBILE, AND ACCESSIBLE. OUR FOLD RAIL STRETCHER IS BUILT ON THE SAME, DURABLE CHASSIS AS THE 5000, AND IS EQUIPPED WITH COLLAPSING SIDE-RAILS FOR UNOBSTRUCTED, FOUR SIDED ACCESS FOR MEDICAL PROFESSIONALS.



#### STANDARD:

HIGH DENSITY FOAM  
MATTRESS WITH CONDUCTIVE COVER

#### OPTIONAL:

VISCOLASTIC, PRESSURE  
SENSITIVE MATTRESS, WITH FOUR-WAY,  
ELASTATEX STRETCH COVER

## MOVEMENT

**HEIGHT ADJUSTMENT:** 56CM TO 90CM  
**HANDS-FREE TRENDELENBERG:** +18° TO -18°  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
FOLD DOWN SIDE-RAILS  
20CM ALL-SURFACE, LOCKING CASTERS  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM

## DIMENSIONS

**TOTAL DIMENSIONS:** 206CM X 83CM  
**LITTER DIMENSIONS:** 206CM X 72CM  
**MATTRESS DIMENSIONS:** 191CM X 74CM X 8CM  
**FOLD DOWN SIDE-RAILS:** 154CM X 32CM

**WEIGHT CAPACITY:** 700LBS / 318KG

# COMPACT 4500N / 4500NS

HIGHLY ADAPTABLE, OUR COMPACT STRETCHER IS SUITED TO YOUR SMALL-SPACE, CLINICAL ENVIRONMENT. PROVIDES MANY OF THE SAME FEATURES AS THE 4500, INCLUDING COLLAPSING SIDE-RAILS. AVAILABLE IN NARROW AND SHORT & NARROW.



**STANDARD:**

HIGH DENSITY FOAM  
MATTRESS WITH CONDUCTIVE COVER

**OPTIONAL:**

VISCOLASTIC, PRESSURE  
SENSITIVE MATTRESS, WITH FOUR-WAY,  
ELASTATEX STRETCH COVER



## MOVEMENT

**HEIGHT:** 56CM TO 90CM  
**HANDS-FREE TRENDELENBERG:** +18° TO -18°  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
20CM ALL-SURFACE, LOCKING CASTERS  
FOLD DOWN SIDE-RAILS

## DIMENSIONS

**NARROW**  
**TOTAL DIMENSIONS:** 206CM X 74CM  
**LITTER DIMENSIONS:** 206CM X 64CM  
**MATTRESS DIMENSIONS:** 198CM X 64CM X 8CM

**SHORT & NARROW**  
**TOTAL DIMENSIONS:** 188CM X 74CM  
**LITTER DIMENSIONS:** 188CM X 64CM  
**MATTRESS DIMENSIONS:** 188CM X 64CM X 8CM

**WEIGHT CAPACITY: 500LBS / 227KG**

# GYNAECOLOGY 4500GYN

FOR USE IN OBSTETRICS AND GYNAECOLOGY, THE 4500GYN STRETCHER IS EQUIPPED WITH A TWO-PIECE SECTIONAL, VISCOLASTIC MATTRESS. VERSATILE, IT BOASTS SUPINE TO FULL-CHAIR POSITIONING WITH PNEUMATIC ASSIST FOR EASE OF USE.



SECTIONAL, PRESSURE SENSITIVE  
FOAM MATTRESS WITH VISCOLASTIC,  
FOUR-WAY STRETCH COVER

## MOVEMENT

**SLIDING TOP DECK:** STRETCHER TO CHAIR  
SUPINE TO FULL CHAIR POSITIONING  
INTEGRATED, ADJUSTABLE HEEL STIRRUPS  
**HEIGHT ADJUSTMENT:** 61CM TO 93CM  
**HANDS-FREE TRENDELENBERG:** +18° TO -18°

**PNEUMATIC ASSIST FOWLER:** 0° TO 90°  
**PNEUMATIC ASSIST LEG SECTION:** 0° TO 105°  
20CM ALL-SURFACE, LOCKING CASTERS  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM

## DIMENSIONS

**TOTAL DIMENSIONS:** 198CM X 81CM  
**FOLD DOWN SIDE RAILS:** 122CM X 30CM SECTIONAL  
**MATTRESS DIMENSIONS:** 188CM X 69CM X 10CM

**WEIGHT CAPACITY:** 700LBS / 318KG

# CADAVER & BARIATRIC CADAVER 5000C & 5000BRC

FEATURING FLAP ACCESS FOR PERSONAL BELONGINGS,  
AND AN EASY-CLEAN BASE WITH INTEGRATED STORAGE,  
OUR CADAVER STRETCHER IS A SIMPLE, FUNCTIONAL DESIGN  
FOR USE IN CADAVER HANDLING UP TO 600LBS (C)  
OR 750LBS (BRC).



## MOVEMENT

**HEIGHT:** 56CM TO 90CM (C)  
OR 56CM TO 81CM (BRC)  
ACCU-STEER 5TH WHEEL  
POSI-LOCK BRAKING SYSTEM  
20CM ALL-SURFACE, LOCKING CASTERS  
OPTIONAL FOUR-SIDED ACCESS

## DIMENSIONS

**5000C**  
TOTAL DIMENSIONS: 198CM X 74CM  
RECESSED LITTER: 196CM X 74CM  
CANOPY: 193CM X 71CM X 36CM  
**5000BRC**  
TOTAL DIMENSIONS: 206CM X 81CM  
RECESSED LITTER: 203CM X 81CM  
CANOPY: 203CM X 81CM X 61CM

## FEATURES

- VINYL CANOPY COVER WITH FLAP ACCESS
- LIGHTWEIGHT, REMOVABLE CANOPY FRAME
- STAINLESS STEEL RECESSED LITTER SHEET

**WEIGHT CAPACITY (C):** 700LBS / 318KG  
**WEIGHT CAPACITY (BRC):** 750LBS / 341KG

# I.V. STAND MODEL 802A

MANUFACTURED WITH HEAVY GAUGE  
COMPONENTS FOR STABILITY AND DURABILITY



## FEATURES

**STAINLESS STEEL  
CONSTRUCTION**

**LOWEST COST  
OF OWNERSHIP  
BASED ON YEARS  
IN SERVICE &  
DURABILITY**

**INTEGRATED  
STAINLESS STEEL  
PATIENT PULL**

**TWO-PIECE  
CONSTRUCTION  
FOR EASY  
ASSEMBLY AND  
MAINTENANCE**

**FOUR INTEGRATED  
HOOKS FOR  
MULTIPLE  
ATTACHMENTS**

**OPTIONAL HOOKS  
AND MONITOR  
STANDS**

## SPECIFICATIONS

58CM 5-LEG BASE  
HIGH DURABILITY 8CM RUBBER CASTERS  
20CM HOOKS  
FULL EXTENSION TO 254CM  
FULL COMPRESSION TO 152CM

# COLOUR GUIDE BUMPERS & ACCENTS



1 Wraparound Bumper

2 Colour Accents



Antique Rose



Dark Rose



Classic Burgundy



Tropical Fuschia




Cranberry Red



Platinum Grey



Mint Green



Sea Mist



Turquoise Blue



Graphite



Teal Blue



Nautical Blue

\* Colour swatches are approximate. Techlem Medical Corporation reserves the right to change specifications without notice

# ACCESSORIES GUIDE

Transforming  
the standard of stretcher  
support surfaces.

Techlem patient surfaces and OFM replacement mattresses are handcrafted in North America with patients and medical professionals in mind. We use high grade materials sourced domestically to produce top quality products that are the envy of the industry

Techlem is proud to offer a catalogue of customizable sizing and material options to suit your clinical requirements.

- **WARRANTY PROTECTED**
- **INDUSTRY BEST VALUE OVER TIME**
- **HIGHLY DURABLE**
- **CHEMICAL-RESISTANT**
- **CUSTOM+STANDARD DESIGNS**
- **EASY TO CLEAN**
- **TEAR-RESISTANT**
- **AVAILABLE IN 8, 10, AND 13CM**



# MATTRESSES



## **LATEX-FREE**

Safe for use in hospital and clinical settings.



## **ANTI-BACTERIAL**

Foam and cover material.



## **DOUBLE-INSULATED**

Keeping the core material intact and fully functional.



## **FLAME-RETARDENT**

Meet or exceeds CTB 117.

Users typically experience  
8+ years of useful life  
from Techlem patient  
support surfaces.  
References are available  
upon request.



# ACCESSORIES GUIDE

S - STANDARD    O - OPTIONAL    N/A - NOT AVAILABLE

ACCESSORIES	5000	5000BR	5000W	1100/N	5000XR	7600	4500-N/NS	4500GYN	5000C/BRC
Aereus Shield Side Rail Upgrade (B500A)	O	O	O	O	O	N/A	N/A	N/A	N/A
Knee Lift (5000KG)	O	O	O	O	O	N/A	O-N/A	N/A	N/A
Vertical O2 Tank Holder (4012)	O	O	O	O	O	S	O	O	N/A
Removable, Adjustable I.V. Pole (3113)	O	O	O	O	O	S	O	O	N/A
Push-Pull Handles Set (5009)	O	O	O	O	O	N/A	O	O	O
Push-Pull Handles, Aereus Set (5009A)	O	O	O	O	O	N/A	O	O	O
Permanent I.V. Pole (5010)	O	O	O	O	O	N/A	O	O	N/A
Permanent I.V. Pole, Aereus Shield (5010A)	O	O	O	O	O	N/A	O	O	N/A
Head/Footboard/Chart Holder (3101)	O	O	O	O	O	N/A	O	O	N/A
Litter Extension (3102)	O	O	O	O	O	N/A	O	N/A	N/A
I.V. Armboard (5003)	O	O	O	O	O	N/A	O	O	N/A
Articulating Armboard (5005)	O	O	O	O	O	N/A	O	O	N/A
I.V. Tow Clamp (5004)	O	O	O	O	O	N/A	O	O	N/A
Restraint Straps (3104)	O	O	O	O	O	N/A	O	O	N/A
Safety Side Padding (5110)	O	O	O	O	O	N/A	N/A	N/A	N/A
Paper Holder (2109)	O	O	O	O	O	N/A	O	N/A	N/A
Food/Instrument Tray (5007)	O	O	O	O	O	N/A	O	O	N/A
Department I.D. Labels (5008)	S	S	S	S	S	O	S	S	N/A
Monitor Stand (5011)	O	O	O	O	O	N/A	O	O	N/A
Stirrups With Adapter (4108)	O	O	O	O	O	N/A	O	S	N/A
Integrated Utility Tray (5019)	S	N/A	N/A	S	S	N/A	S	S	N/A
Double Articulating Headrest (4015)	O	O	O	O	N/A	N/A	O	O	N/A
Utility Hitch For I.V. Pole (5004)	O	O	O	O	N/A	N/A	O	O	N/A
Lateral Wrist Rest (Optional for 4015) (401502)	O	O	O	O	N/A	N/A	O	O	N/A
Superior Wrist Rest (Optional For 4015) (401501)	O	O	O	O	N/A	N/A	O	O	N/A
Crank Fowler To Replace Pneumatic (5000C)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
X-Ray Cassette Holder - Back Section (5014)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
X-Ray Cassette Holder - Lateral C-Spine (5015)	O	N/A	N/A	O	O	N/A	O	N/A	N/A
Crank Fowler To Replace Pneumatic (5000XLC)	N/A	O	O	N/A	N/A	N/A	N/A	N/A	N/A
Extra Height Side-Rails Set (5111)	O	O	O	O	O	N/A	N/A	N/A	N/A
Head/Footboard/Chart Holder- Extra Wide (3101XL)	N/A	O	O	N/A	O	N/A	N/A	N/A	N/A
Litter Extension - Extra Wide (3102XL)	N/A	O	O	N/A	O	N/A	N/A	N/A	N/A
Integrated Utility Tray - Extra Wide (5019XL)	N/A	S	S	N/A	O	N/A	N/A	N/A	N/A
Powered Actuator - Fowler (5020)	O	O	O	O	O	N/A	O	O	N/A
Powered Actuator - Knee Lift (5021)	O	O	O	O	O	N/A	O	N/A	N/A

# ACCESSORIES GUIDE

S - STANDARD    O - OPTIONAL    N/A - NOT AVAILABLE

ACCESSORIES	5000	5000BR	5000W	1100/N	5000XR	7600	4500-N/NS	4500GYN	5000C/BRC
X-Ray Fowler With Cassette Holder (5018)	O	N/A	N/A	O	O	N/A	O	O	N/A
X-Ray Fowler With Cassette Holder (5018XL)	N/A	O	O	N/A	N/A	N/A	N/A	N/A	N/A
Lift Sheet - Standard (7601C)	N/A	N/A	N/A	N/A	N/A	S	N/A	N/A	N/A
Lift Sheet - Cysto-Obstetrical (7602C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Lift Sheet - Cataract Sheet (7603C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Lift Sheet - Bath (7608C)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Rod For Cataract Sheet (7621)	N/A	N/A	N/A	N/A	N/A	O	N/A	N/A	N/A
Integrated Scale (5022)	O	O	O	O	O	N/A	O	O	O



1300 643 633

PO BOX 508 Morningside Q 4170