



Materni CBM Containers For Hospital Use

Screw top milk container of transparent plastic for clear visual inspection and for accurate volume collection and deliveries.

Suitable for all styles of milk container warmers up to 38° Dishwasher and freezer safe.

Available in four sizes: 50 ml – 80 ml – 120 ml – 240 ml

CBA- Clean by manufacturing. TGA Approved .
TGA Number - 452305

- Made from Polypropylene, fully recyclable
- Safe for pasteuriser
- Bisphenol A and Phthalates FREE
- The milk container for in hospital use is compatible with all breast pumps in the Australian Market.

