

[www.midmed.com.au](http://www.midmed.com.au)



SEPSIS  
TROLLEYS

## Features

- Trolleys designed for use in the treatment of Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- Two locking options
- High Security Bolt Lock - Supplied with two keys
- Electronic Push Button - Programmable four-digit code
- Two adjustable door shelves & an A4 document holder
- Six labelled drawers (Oxygen/Blood Cultures/Intravenous Antibiotics/Intravenous Fluid/Haemoglobin & Lactates/Urine Measurement)
- Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting Plate (fits to the l/h side of the trolley)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- 125mm polymer castors
- When not in use, trolleys storing medicines should be secured to the wall. See option MT/1/6F Security Clamp Sepsis Trolley Features

### Dimensions (w x d x h)

• Overall Dimensions - 511 x 542 x 1072mm • Drawers - 335 x 295 x 113mm (internal)

### Standard Colours

- Paint - Grey White
- Work Surface - Grey White



Electronic Push Button Lock



High Security Bolt Lock



BM.SEPS/6/HSB

High Security Bolt Lock	Electronic Push Button Lock	Description	Weight
BM.SEPS/6/HSB	BM.SEPS/6/EPB	Sepsis Trolley - Six Drawer	53.4kg

## Sharps Box Holder

### Features

- Accepts up to seven litre capacity sharps boxes
- Sharps box secured using a rubber band fitted with a toggle
- Attaches to the left-hand side of the unit (customer to fit - fixings supplied)



Sharps box not supplied

Item	Description	Weight
BM.SEP/SBH	Sharps Box Holder	0.4kg

## Sharps Box Mounting Plate - Universal

### Features

- Compatible with the following...
  - Daniels Healthcare
    - POUDS® Wall-Bracket (up to 5L) - SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
    - POUDS® Wall-Bracket (7L & above) - SHARPSGUARD® 7, 11.5 & 22
  - Frontier (Sharpsafe) - 2, 3, 4 & 7
  - Sharpak Healthcare - Sharpak 18, 35, 60, 120 & 220 PLUS
- Customer to supply the plastic brackets & fit (fixings supplied)



Trolleys are compatible with most of the popular UK sharps box support brackets.



Frontier (Sharpsafe)\*\*



Sharpak Healthcare\*\*



Daniels Healthcare  
Up to 5L\*\*



Daniels Healthcare  
7L & Above\*\*


\*\*Sharps box & bracket not supplied

Item	Description	Weight
BM.SEP/USBMP	Sharps Box Bracket Mounting Plate - Universal	0.2kg

## Security Clamp

### Features

- Allows the trolley to be secured to the 'fabric' of the building
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed and locked
- Recommended wall fitting height 1000mm



#### Important Note - Security Clamps

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home".


"Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB-March 2005"

Item	Description	Cable Length	Weight
BM.MT/1/6F	Security Clamp	950mm	0.8kg

## Drawer Inserts

### Features

- Plastic drawer inserts, height 100mm
- Assists segregation of drawer contents
- Insert supplied with removable & repositionable dividers providing 4 separate areas
- Supplied for self-assembly
- Clip on label holders & plastic labels available to assist identification




BM.SEP/DIS

Item	Description	Weight
BM.SEP/DIS	Drawer Insert Set - Single Door Trolleys	0.5kg

## Label Holders & Labels

### Features

- Plastic label holders & labels
- Assists identification of drawer contents
- Label holders snap onto drawer inserts/dividers
- Accepts ISO labels 40 x 60mm (w x h)
- Sold in packs of ten



Item	Description
BM.SEP/LH	Label Holder
BM.SEP/L	Label

## Label - Sepsis Tiger

### Features

- Sammy the sepsis tiger label developed for use on Paediatric trolleys
- Compatible with both single & double door trolleys



Item	Description
BM.SEP/TIG01	Label - Sammy Sepsis Tiger

## Features

- Trolleys designed for use in the treatment of Sepsis
- Supports delivery of 'The Sepsis Six' which clearly defined locations for equipment used in each of the stages
- Two locking options
- High Security Bolt Lock - Supplied with two keys
- Electronic Push Button - Programmable four-digit code
- Two adjustable door shelves & an A4 document holder
- Six labelled drawers (Oxygen/Blood Cultures/Intravenous Antibiotics/Intravenous Fluid/Haemoglobin & Lactates/Urine Measurement)
- Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting Plate' (fits to the l/h side of the trolley)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- 125mm polymer castors
- When not in use, trolley storing medicines should be secured to the wall.  
See option MT/1/6F Security Clamp Sepsis Trolleys

## Dimensions

- Overall Dimensions - 550 x 550 x 993mm
- Plastic Top - 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray - 355 x 405 x 100mm (internal) •
  - Deep Tray - 355 x 405 x 150mm (internal)

## Standard Colours

- Paint - Orange

Description	Item	Weight
6 x 100mm Trays, High Security Bolt Lock	BM.CTSEPSM1/10/6S/HSB	36.8kg
6 x 100mm Trays, Electronic Push Button Lock	BM.CTSEPSM1/10/6S/EPB	36.8kg



BM.CTSEPSM1/10/6S/HSB

Other products that you may be interested in...

### Resuscitation Trolleys & Carts

