



Premier Series 5000BR Bariatric Stretcher



Features & Specifications:

- Techlem "Posi-Lock" Braking System
- Four Corner Access Brake / Steer Pedals
- 8" (20cm) All Surface, Four Wheel Locking Casters
- Techlem "Accu-Steer" Fifth Wheel Steering
- Full Length Side Rails
- Side Rail Release at Head & Foot End
- Hands Free Trendelenberg / Reverse Trendelenberg:
+12 Degrees to -12 Degrees
- Pneumatic Head Section: 0 Degrees to 90 Degrees
- Raise / Lower Pedals: Side Access
- Weight Capacity - 650 lbs
- Wraparound Bumper / Accents Available in 14 Colors**

Dimensions:

- Height Adjustment: 22" to 32" (56cm to 81cm)
- Total Dimensions: 81" X 35.5" (206cm X 90cm)
- Litter Dimensions: 81" X 32" (206cm X 81cm)
- Full Length Side Rails: 75" X 16" (190cm X 41cm)
- Mattress Dimensions: 76" X 32" (193cm X 81cm)
Optional: 78" X 32" (198cm X 81cm)

Optional Accessories:

- Raise / Lower Pedals: 4 Side Access
- Mattress: 4" (10cm) or 5" (13cm) Thickness
- Viscolastic Pressure Reduction Mattress:
4" (10cm) or 5" (13cm) Thickness
- Push / Pull Steering Handles
- I.V. Pole - Stainless Steel Permanent Folding
Head / Foot End
- I.V. Pole - Removable, Adjustable Telescoping
- I.V. Stand Utility Hitch
- I.V. Armboard
- Knee Gatch
- Head / Foot Rail
- Head / Footboard with Chart Holder
- Restraint Straps
- Vertical Oxygen Tank Holder
- Department Identification Labels: Head / Foot End

** See the "Techlem Wraparound Bumper and Accents Color Guide"

Techlem Medical Corporation reserves the right to change material and specifications without notice.
All sizes are approximate.

Techlem
MEDICAL CORPORATION

6890 Pacific Circle, Mississauga
Ontario, Canada, L5T 1N8
Phone: 1-800-387-4848, 905-670-5000
Fax: 905-670-3925
E-mail: techlem@istar.ca

5000BR - 0905